

Product datasheet for **SC301403**

LST 3TM12 (SLCO1B7) (NM_001009562) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	LST 3TM12 (SLCO1B7) (NM_001009562) Human Untagged Clone
Tag:	Tag Free
Symbol:	LST 3TM12
Synonyms:	LST-3; LST-3TM12; LST3; SLC21A21
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC301403 representing NM_001009562.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTA
CCGAGGAGATCTGCCGCCGCGATCGCCGGCGCGCC
ATGAAAATATCAACCACTCAAATAGAAAGGAGATTTGAGATATCCTCTTCTCTTGTGTTAATTGAT
GGAAGCTTCGAAATAGGAAATTTGTTTGTGATTGATTTGTAAGTTACTTTGGATCTAACTACACAGA
CCGAAGTTAATTGGAATTGGTTGTTTCTTATGGGAACCTGGAAGTATTTGATGGCTTTACCACATTTC
TTCATGGGATATTACAGGTATTCTAAAGAAACCAATATTGATCCATCAGAAAATTCAACATCAAACCTTA
CCAAACTGTTTAATTAATCAAATGTTATCACTCAATAGAACACCGTCTGAGATAATAGAAAGAGTTGT
GTGAAGGAATCTGGGTCACACATGTGGATCTATGTCTTCATGGGAATATGCTTCGTGGCATAGGGGAA
ACCCCATAGTACCATTGGGATTTCTTACATTGATGATTTTGCAAAGAAGGACATTCTCCTTGAT
TTAGGTACTGTGAATGAATGGGAATGACTGGTCTAGTTTTGCCTTTATGCTGGGATCTCTGTTTGCT
AAAATGTATGTGGATATCGGATATGTGGATCTGAGCACTATCAGGATAACTCCTAAGGACTCTCGTTGG
GTTGGAGCTTGGTGGCTTGGTTTCTTGTGTCTGGAATAGTATCCATTATTTCTTCTATACCATTCTT
TTCTTGCCTCTAAATCCAAATAAACACAGAAAGAAAGGAAAGTTTCACTATTTTGCATGTGCTAAAA
ACTAATGATAAAAGGAATCAAATAGCTAATTTGACCAACCGAAGAAAATATATTACCAAAAATGTGACT
GGCTTTTTCCAGTCTTTGAAAAGCATCCTTACCAATCCCCTGTATGTTATATTTGTAATTTTACATTG
TTACACATGAGCAGCTACATTGCTTCTTACTTATATCATTAAATGGTGGAGCAACAGTATGGTTGG
TCTGCATCTAAGCTAATTTTTGTTGGGAGTCTCGCCCTACCTGCTGTTGCAATTGGCATGTTTTCA
GGAGGATATATCATTAAAAAATCAAATTGTCTTTAGTTGGACTTGCCAAATTGGCATTGTTCTGCA
ACAGTGCATCTCTTATCTCAAGTTTTATTTCTTTCTAATCTGTGAAAGCAAATCAGTTGCCGGCCTA
ACCTTGACCTATGATGGAATAGTCCAGTAAGATCTCATGTAGATGTACCACTTTCTTATTGCAACTCA
GAGTGCAATTGTGATGAAAGTCAATGGGAACCGCTCTGTGAAACAATGGAATAACTTACCTGTACCT
TGTCTAGCAGGATGCAAATCTTCAAGTGGTAATAAAGAGCCCATAGTGTGTTTATACTGTAGCTGTGTG
GAAGTAATTGGTCTCCAGAACAAAAATTACTCAGCGCACTTGGGTGAATGCCAAGAGATGATGCTTGT
ACAAGGAAATCTTACGTTATTTTGTAAATCAAGTCTTAGATGCTTCTCTGTGCAGTTGGACTTACC
TCATATTCCGTGCTGGTGATTAGGATTGTTCAACCTGAATTGAAAGCACTTGCAATCGGCTTCCATTCA
ATGATTATGCGATCGCTAGGAGGTATTCTAGTTCCAATATATTTGGGGCTGATTGATACAACGTGT
ATGAAGTGGTCCACCAACAGCTGTGGAGCACGAGGGGCTGTAGGATATATAATTCCACATATTTGGGA
AGAGCCTTCTTCGGCTTGAAGGTAGCCCTAATATTTCCAGTACTGTTTTACTTACTGTATTTATTTT
GTTGTAAGGAAAAATCCCATGGAAGGATACCAAGTATTAGAAAATGAAAGACAAGTAATGGATGAA
GCAAACCTAGAATTCTTAAACGACAGTGAACATTTGTACCTTCTGCTGAAGAACAGTAA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGCGC
  
```

Restriction Sites: AscI-MluI

ACCN: NM_001009562

Insert Size: 1923 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001009562.4](#)

RefSeq Size: 1923 bp

RefSeq ORF: 1923 bp

Locus ID: 338821

UniProt ID: [G3V0H7](#)

Cytogenetics: 12p12.2

Protein Families: Transmembrane

MW: 71.2 kDa