

Product datasheet for **SC301435**

C12orf29 (NM_001009894) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	C12orf29 (NM_001009894) Human Untagged Clone
Tag:	Tag Free
Symbol:	C12orf29
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC301435 representing NM_001009894. Blue=Insert sequence Red=Cloning site Green=Tag(s)

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCC**CGATCGCC**
ATGAAGCGCTTGGGCTCCGTGCAGCGGAAAATGCCGTGTGTGTTGTGACGGAGGTGAAAGAGGAGCCT
TCCTCCAAAGGGAGCATCAGCCATTTAAAGTTTTGGCAACTGAACTGTAAGTCACAAGGCATTAGAT
GCAGATATATACAGTGAATTCACAGAAAAAGTGGATGGAACATGTTGTTATGTTACTACCTATAAA
GATCAGCCATACCTTTGGGCTCGACTAGATAGAAAACCAACAACAAGCTGAAAAAGATTTAAAAAT
TTTCTACATTCAAAGAAAACCCAAAAGAATTTTTTGGAACGTTGAGGAGGACTTCAACCAGCTCCG
GAGTGCTGGATACCAGCAAAGGAGACAGAACAATAAATGGGAACCCAGTGCCTGATGAAAATGGACAC
ATTCTCTGGTTGGGTACCGTAGAGAAAAACAACAACAGTATTGCTGGCATTCTCTGTAGTTAATTAT
GAATTTGAAATTGCCCTGGTACTAAAACATCATCCTGATGATTCTGGACTTTTGGAAATTAGTGCAGTG
CCACTTTCAGATCTCTTAGAACAAACACTGGAACATAGGAACAAATATCAATGGAAACCCATATGGG
TTAGGGAGCAAGAAACATCCATTACATCTTCTTATACCACATGGAGCATTTCAAATAAGAAATCTACCT
TCATTGAAGCACAATGATCTCGTGTCTGGTTTGAAGATTGCAAAGAGGGTAAATTTGAAGGAATAGTG
TGGCATTGCAGTGATGGCTGTTAATAAAGGTCCATCGCCATCATCTGGCTTATGCTGGCAATTCCA
GATACTTACATGAATTCAGACAGTTATTATCAACATGAACCTGAACAAATGTGACTCTGCCTTTGAT
ATTAAGTGTGTTTAAATCATTTTTTAAAAATAGATAATCAGAAATTTGTTAGACTCAAAGATATAATA
TTTGATGTAA**TA**
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC

Restriction Sites:	SgfI-MluI
ACCN:	NM_001009894
Insert Size:	978 bp



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001009894.2</u>
RefSeq Size:	2895 bp
RefSeq ORF:	978 bp
Locus ID:	91298
UniProt ID:	<u>Q8N999</u>
Cytogenetics:	12q21.32
MW:	37.5 kDa