

Product datasheet for **SC301465**

ZNF648 (NM_001009992) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ZNF648 (NM_001009992) Human Untagged Clone
Tag:	Tag Free
Symbol:	ZNF648
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC301465 representing NM_001009992.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGCACAAAGTGGACTCCAGGACAGGTGGGAGAGGCGTCTCCTCTCAGCAGCTTGACTGAGGAAGCT
CATGACACCCAGATGCTGAGCATGAACCTTAGAGAGTGACGATGAAGATGGTGGGAGGCCGAAAAAGAG
GGCACCCTGACCCGGTGGCCTGTCCAAGGGGACGCTCCCAAGTAACACACGAAAATCCTGACTTGCCA
TGGCCCATCCACTGGGCAAAGAGGAAGAGAAATTCTCTGACTCCTCCAGTGCTGGGGGCATGGGGCAG
AAACAGTGGAATGTCTGGGAAAGCCAGTTGGAGCAGAGATGTGACAAAGATCAACGAGACCCAGGGT
TCCCCGGGAGCAAGCAGAGCTCTGGGTTCCCTTCCAGTGGTCTCGCACACAAATTGTTAGGTGAGATG
CAACCTCTTGGGACCGACTACCTGCGGGTGATGATGGATACTCGGGGGCAAACAGGACGAGCTTTG
GATGTCCACCCAGCTTCCCAAGCAATGGAAAGTATCTCTGTGCGCACAAAGTGTAGACACGTCCGCA
GGGAACCTTCTCTGTTGTGTTTCCCAAGGCCGGGAGCAACTGGGACCTTCCACGCAAGAGACACAT
ACACCAGCCAGGCGTCGCCACCCAGCCAGCCTGGTGCCGCGGTCTGGCAAAGCGCGGAACAGC
AGGAAAGTACAGAACCAGGCGGGCCGGCGCAGGGCGGAGAGGCTGAGGCGCGTCCCTACAGGTGCCTG
CGGGGCGGGCGGCCCTTTCAGAAGCCAGCAAGCCGCTGAGCCCCGCGGAGACGCGCGGCCGCCGCC
AAGCGCTACGCGTGCGAGCTATGCGGAAGGCCTACTCCACCGCGGCACACTCCAGCAGCACAGGCGC
CTGCACACGGGCGAGCGGCCCTACCAAGTCTCTTCTGCGACAAGGCCTACACCTGGTCTCCGACCAC
CGGAAGCACATCCGACCCACACAGGCGAGAAACCCTACCCGTGTCCAGACTGCGGAAGGCCTTCGTG
CGCTCTTCGGACCTGCGCAAACACAGCGCAACATGCACAGCAACAATAAGCCCTTCCCGTGTCCGAG
TGCGGCGTGACCTTCAACAAGCCGCTGTGCTGCTGCGCCACAGCGCACGACCTGGGCGCCAAGCCC
TTCCGCTGCCCGCCTGCGACCGGAGTTCGCTGTGGCCAGCCGATGGTGGAGCACCAGCGCGTGAC
TCGGGCGAGCGGCCCTTCCCTGCCCCACCTGCGGCAAGTGCTTACCAAGTCTCCAATCTGTCCGAG
CACCAGACGCTGCACACCGGCCAGAGGCCCTTCAAGTGCCTGACTGCGCGTGCCCTTCGCGCAGCCC
TCGCGCCTCGTGCGCCACAGCGCATCCCACTGGCGAGAGGCCCTTTCCTTGACGCGAGTGCGCCAG
GCCTTTGCCCGCTCTTCGACCCTGAAGCGGCACCAACAGATCCACTCCGGGGAGAAGGGATTCTCTGT
GCCGAGTGCGGCAGGGCCTTCCGATTGCCTCTGAGTTGGCCAGCACATACGAATGCACAACGGAGAG
AGGCCCTACCAAGTGAGGACTGCGGCCAGGCCCTTACCAGGTCCAATCACCTCCAACGACACCGAGCC
AAGCACGGCACCTGCAAGAAGGAGCCCATCCCTTCTCTCTGACGAGTGA
AGCGGACCGACGCTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGAT
ATCCTGGATTACAAGGATGACGACGATAAGGTTTAA
  
```

Restriction Sites: SgfI-RsrII

ACCN: NM_001009992

Insert Size: 1707 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001009992.1</u>
RefSeq Size:	3649 bp
RefSeq ORF:	1707 bp
Locus ID:	127665
UniProt ID:	<u>Q5T619</u>
Cytogenetics:	1q25.3
MW:	62.3 kDa
Gene Summary:	May be involved in transcriptional regulation.[UniProtKB/Swiss-Prot Function]