

## Product datasheet for **SC301472**

### **ARHGAP28 (NM\_001010000) Human Untagged Clone**

#### **Product data:**

Product Type:	Expression Plasmids
Product Name:	ARHGAP28 (NM_001010000) Human Untagged Clone
Tag:	Tag Free
Symbol:	ARHGAP28
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

**Fully Sequenced ORF:** >SC301472 representing NM\_001010000.  
 Blue=Insert sequence Red=Cloning site Green=Tag(s)


```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGCATCGCC
ATGAGGAAAAAGGATAAGCAATCTATCAGGGATGTCAGAGACATTTTTGGAGTCAGTGAATCTCCTCCT
CGTGATACCTGTGGCAACCACACTAATCAGCTGGATGGCACCAAGGAAGAAAGAGAGCTTCCAAGAGTT
ATCAAGACAAGTGGTCCATGCCAGATGATGCTTCTCTCAACAGTACTACCCTGTCTGACGCATCCCAG
GATAAAGAAGGGAGTTTTGCGGTTCCAGGAGTACTCTGTGGCTATACTTGAGACCATTCCAGTTCTA
CCAGTTCAATCCAAATGGATCACCGGAGCCTGGACAGCCAGTTCAGAATGCGATAAGTGAATGATGATTTT
CTGAAAAAGAACATCCACCAGAGGCTGAAGAGCTGTCATTTGAAAGTGTCTTATTCAGAAATGGTTACG
GAGGCTCTAAAAAGAAATAAACTTAAGAAATCAGAGATTAAGAAAGAAGACTATGTTTTAACTAAATTT
AATGTTCAGAAAACCAGATTTGGCTTAAGTGAAGCAGGAGATCTGTCTGCTGAAGACATGAAGAAAAATC
CGCCATCTCTCTGATTGAATTGACTGCCTTTTTGATGCCTTTGGAATCAACTGAAAAGGAACAAA
ACAGAGAAAGTAAAAGGACGAGACAATGGGATTTTTGGAGTTCCACTTACAGTCTCCTGGACGGTGAC
CGAAAGAAAGACCTGGAGTGAAAGTCCCTGGTATTACAAAAATTTTTGAGAAAGTTGAGGAATCA
GGTCTGGAATCTGAAGGAATTTTCGACTTTCAGGATGTACTGCTAAAGTCAAGCAATACCGTGAAGAA
CTTGATGCCAAGTTAATGCTGATAAAATTTAAATGGGACAAAAATGTGCCATAGAGAAGCTGCAGTAATG
TTGAAAGCGTTTTTCAGAGAACTACCCACCTCTCTTCCCTGTGGAATATATACCTGCCTTCATCAGT
CTAATGGAAGAGGGCCTCACGTCAAAGTACAGTTTCAAGCCTTACACCTCATGGTCATGGCGCTGCCT
GATGCCAACAGAGATGCAGCTCAGGCCCTCATGACATTTCAATAAAGTGAATGCCAATGAATCAAAA
AACCGAATGAGTCTGTGGAACATTTCTACAGTGATGGCACCAACCTTTTCTTTCAGTAGAAGCAAACAC
TCTGATTATGAAGAATTACTGTTAGCAAACACTGCGGCCACATCATCCGCCTAATGCTTAAGTACCAG
AAGATTTTGTGGAAGTTCCATCTTTCTAATCACTCAAGTAAGAAGAATGAATGAAGCCACGATGCTA
TTGAAAGAAGCAGCTCCCAAGTGTGAGGAAGCTGCTCAGGAGGAAGACCCTCGAGCGGGAGACTGCAAGC
CCCAAGACTTCAAAGGTACTGCAAAAAATCACCCCTCGGCAAGACGAATGTCTGACGTGCCGGAAGGAGTC
ATACGGGTCCATGCTCCACTTCTCTCCAAGGTGTCATGGCCATTCAACTCAACAATCAAACCAAAGCC
AAAGACATATTGGCAAAATTTCAATATGAAAACAGTCAATGGTTTCATCAGAATGTATTAAGATTGAGAC
CAAAGTTATATGAAATTGGAGGAAATATAGGAGAGCATTGCTTGGATCCAGATGCTTATATATTGGAT
GTATATCGTATAAATCCTCAAGCAGAATGGGTGATTAACCCCAACAAAGTTCTTAA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGCGC
  
```

**Restriction Sites:**

SgfI-MluI

**Plasmid Map:**



<b>ACCN:</b>	NM_001010000
<b>Insert Size:</b>	1713 bp
<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>OTI Annotation:</b>	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"> <li>1. Centrifuge at 5,000xg for 5min.</li> <li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li> <li>3. Close the tube and incubate for 10 minutes at room temperature.</li> <li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li> <li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li> </ol>
<b>RefSeq:</b>	<u><a href="#">NM_001010000.2</a></u>
<b>RefSeq Size:</b>	5542 bp
<b>RefSeq ORF:</b>	1713 bp
<b>Locus ID:</b>	79822
<b>Cytogenetics:</b>	18p11.31
<b>MW:</b>	64.6 kDa