

Product datasheet for **SC301547**

CYP4A22 (NM_001010969) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	CYP4A22 (NM_001010969) Human Untagged Clone
Tag:	Tag Free
Symbol:	CYP4A22
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC301547 representing NM_001010969.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)


```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGAGTGTCTCTGTCTCTGAGCCCCAGACAGCCCTGGGTGGTGTCTCCGGGATCCTCCAAGTGACCTCC
CTGCTCATTCTGCTTCTGCTGTGATCAAGGCAGCTCAGCTCTACCTGCATAGGCAGTGGCTGCTCAA
GCCCTCCAGCAGTTCCTCGTGCCTCCCTCCCACTGGCTCTTCGGGCACATCCAGGAGTTCCAACAGCAG
CAGGAGCTACAACGGATTCAAGAACGGGTGAAGACATTCCAAGTGCCTGTCTTATTGGATATGGGGA
GGCAAAGTTCGTGTCCAGCTCTATGACCCTGACTATATGAAGGTGATTCTGGGGAGATCAGACCCGAAA
TCCCATGGATCTACAAATTCCTGGCTCCACGGATTGGGTACGGCTTGCTCCTGTTGAATGGGCAGACA
TGGTTCAGCATCGACGGATGCTGACCCAGCCTTCCACAATGACATCCTGAAGCCATACGTGGGGCTC
ATGGCAGACTCTGTACGAGTGTGCTGGACAAATGGGAAGAGCTCCTGGCCAGGATTCCCTCTGGAG
GTCTTTCAGCACGTCTCCTTGATGACCCTGGACACCATCATGAAGAGTGCCTTCAGCCATCAGGGCAGC
ATCCAGGTGGACAGGAATTCTCAGTCTACATCCAGGCCATTAGTGACCTGAACAGCCTGGTTTTTTGC
TGTATAGGAATGCCTTTCATGAGAATGACACCATCTACAGCCTGACCTGTGCTGGCCGCTGGACACAC
CGCGCCTGCCAGCTGGCCCATCAGCACACAGACCAAGTATCCAAGTGAAGGCTCAACTACAGAAG
GAGGGGAGCTGGAGAAGATCAAGAGGAAGAGGCACTTGGATTTTCTGGACATCCTCCTTGGCCAAA
ATGGAGAATGGGAGCATCTTGTGAGACAAGGACCTCCGTGCTGAGGTGGACACGTTTATGTTGAGGGC
CACGACACCACAGCCAGTGGGATCTCCTGGATCCTTATGCTCTGGCCACACACCCCAAGCATCAGGAG
AGGTGCCGGGAGGAGATCCATGGCCTCCTGGGTGATGGAGCCTCCATCACCTGGAACCACCTGGACCAG
ATGCCCTACACCACCATGTGCATTAAGGAGGCACTGAGGCTTACCCACCGGTGCCAGGCATTGGAAGA
GAGCTCAGCACTCCCGTCACTTCCCTGATGGCGCTCCTTGCCCAAAGGTATCATGGTCTCCTCTCC
ATTTATGCCTTACCACAAACCCAAAAGTGGCCCAACCTAGAGGTGTTTGACCTTCCCGTTTTGCA
CCGGGTTCTGCTCAACACAGCCACGCTTTCCTGCCCTTCTCAGGAGGATCAAGGAAGTGCATCGGGAAA
CAATTTGCCATGAACCAGCTGAAGGTGGCCAGGGCCCTGACCCTGCTCCGCTTTGAGCTGCTGCCTGAT
CCCACCAGGATCCCCATCCCCATGGCACGACTTGTGTTGAAATCCAAAAATGGAATCCACCTGCGTCTC
AGGAGGCTCCCTAACCTTGTGAAGACAAGGACCAGCTTTGA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites:

SgfI-MluI

Plasmid Map:



ACCN:

NM_001010969

Insert Size:	1560 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001010969.3
RefSeq Size:	1678 bp
RefSeq ORF:	1560 bp
Locus ID:	284541
UniProt ID:	Q5TCH4
Cytogenetics:	1p33
Protein Families:	P450, Transmembrane
Protein Pathways:	Arachidonic acid metabolism, Fatty acid metabolism, Metabolic pathways, PPAR signaling pathway, Retinol metabolism, Vascular smooth muscle contraction
MW:	59.2 kDa
Gene Summary:	This gene encodes a member of the cytochrome P450 superfamily of enzymes. The cytochrome P450 proteins are monooxygenases which catalyze many reactions involved in drug metabolism and synthesis of cholesterol, steroids and other lipids. This gene is part of a cluster of cytochrome P450 genes on chromosome 1p33. Alternative splicing results in multiple transcript variants encoding different isoforms. [provided by RefSeq, Apr 2015] Transcript Variant: This variant (1) represents the longer transcript and encodes the longer isoform (1).