

Product datasheet for **SC301602**

PNLIPRP3 (NM_001011709) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	PNLIPRP3 (NM_001011709) Human Untagged Clone
Tag:	Tag Free
Symbol:	PNLIPRP3
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC301602 representing NM_001011709. Blue=Insert sequence Red=Cloning site Green=Tag(s)


```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGCTTGGAAATTTGGATTGTTGCATTCTTGTCTTTGGCACATCAAGAGGAAAAGAAGTTTGCTATGAA
AGGTTAGGGTGTTCAAAGATGGTTTACCATGGACCAGGACTTTCTCAACAGAGTTGGTAGGTTTACCC
TGGTCTCCAGAGAAGATAAACACTCGTTTCTGCTCTACACTATACAAATCCCAATGCCTATCAGGAG
ATCAGTGGGTTAATTCTTCAACTATCCAAGCCTCATATTTTGGAAACAGACAAGATCACCCGTATCAAC
ATAGCTGGATGGAAAACAGATGGCAAATGGCAGAGAGACATGTGCAATGTGTTGCTACAGCTGGAAGAT
ATAAATGGCATTAAATTTAGATTGGATCAACGGTTCACGGGAATACATCCATGCTGTAACAATCTCCGT
GTTGTTGGTGCTGAGGTGGCTTATTTTATTGATGTTCTCATGAAAAAATTTGAATATTCCTTCTAAA
GTGCACTTGATTGGCCACAGCTTGGGAGCACACCTGGCTGGGGAAGCTGGGTCAAGGATACCAGGCCTT
GGAAAGAATAACTGGGTTGGACCCAGCTGGGCCATTTTCCACAACACTCCAAAGGAAGTCAGGCTAGAC
CCCTCGGATGCCAACTTTGTTGACGTTATTACATAAAATGCAGCTCGCATCCTCTTTGAGCTTGGTGTT
GGAACCATTGATGCTTGTGGTCACTTGACTTTTACCAAATGGAGGGAAGCACATGCCAGGATGTGAA
GACTTAATTACACCTTTACTGAAATTTAACTTCAATGCTTACAAAAAGAAATGGCTTCTTCTTTGAC
TGTAACCATGCCGAAGTTATCAATTTTATGCTGAAAGCATTCTTAATCCTGATGCATTTATTGCTTAT
CCTTGATAGATCCTACACATCTTTTAAAGCAGGAAATGCTTCTTTTGTCCAAAGAAGGTTGCCCAACA
ATGGGTCAATTTGCTGATAGATTTCACTTCAAAAATATGAAGACTAATGGATCACATATTTTTTAAAC
ACAGGGTCCCTTTCCCATTTGGCCGTTGGAGGCACAAATGTCTGTTAAACTCAGTGGAAGCGAAGTC
ACTCAAGGAAGTGTCTTTCTCGTGTAGCGGGGAGTTAGGAAAAGTGGGGAGTTTGCCATTGTCAGT
GGAAAAGTGGAGCCAGGCATGACTTACACAAAATTAATCGATGCAGATGTTAACGTTGGAAAACATTACA
AGTGTTCAGTTCATCTGGAAAAACATTTGTTTGAAGATTCTCAGAATAAGTTGGGAGCAGAAATGGTG
ATAAATACATCTGGAAATATGGATATAAATCTACCTTCTGTAGCCAAGACATTATGGGACCTAATATT
CTCCAGAACCTGAAACCATGCTAA
ACCGGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGCGC
```



[View online >](#)

Restriction Sites: Sgfl-MluI

Plasmid Map:



ACCN: NM_001011709

Insert Size: 1404 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001011709.2](#)

RefSeq Size: 2309 bp

RefSeq ORF: 1404 bp

Locus ID: 119548

UniProt ID: [Q17RR3](#)

Cytogenetics: 10q25.3
Protein Families: Druggable Genome
MW: 52.3 kDa