

Product datasheet for **SC301609**

Heterogeneous Nuclear Ribonucleoprotein (A1 like) (HNRNPA1L2) (NM_001011724) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Heterogeneous Nuclear Ribonucleoprotein (A1 like) (HNRNPA1L2) (NM_001011724) Human Untagged Clone
Tag:	Tag Free
Symbol:	Heterogeneous Nuclear Ribonucleoprotein (A1 like)
Mammalian Cell Selection:	None
Vector:	<u>pCMV6-XL5</u>
E. coli Selection:	Ampicillin (100 ug/mL)
Fully Sequenced ORF:	>OriGene sequence for NM_001011724 edited CACCTTGGCGTCATGTCTAAGTCAGCGTCTCCAAAAGAGCCCGAACAGCTGAGGAAGCTC TTCATTGGAGGGTTGAGCTTTGAAACAACACTGATGAGAGCCTGAGGAGCCATTTTGAAGCAA TGGGGAACGCTCACAGACTGTGTGGTAATGAGAGATCCAAACACCAAGCGCTCCAGGGGC TTTGGGTTTGTGCATATGCCACTGTGGAGGAGGTGGATGCAGCTATGAATACAACGCCA CACAAGGTGGATGGAAGAGTTGTGGAACCAAAGAGAGCTGTCTCCAGAGAAGATTCTCAA AGACCAGGTGCCCACTTAACTGTGAAAAAGATATTTGTTGGTGGCATTAAAGAAGACACT GAAGAACATCACCTAAGAGATTATTTTGAACAGTATGGAAAAATTGAAGTAATTGAAATC ATGACTGACCGAGGCAGTGGCAAGAAAAGGGGCTTTGCCTTTGTAACCTTTGACGACCAT GACTCCGTGGATAAGATTGTCATTTCAGAAATACCATACTGTGAAGGGCCACAACCTGTGAA GTTAGAAAAGCCCTGCCAAGCAAGAGATGGCTAGTGCTTCATCCAGCCAAAGAGGTCGA AGGGGTTCTGGAACTTTGGTGGTGGTGGTGGAGATGGTTTCGGTGGGAATGACAACCTTT GGTCGTGGAGGAACTTCAGTGGTGGTGGTGGCTTTGGTGGCAGCTGTGGTGGTGGTGGAA TATGGTGGCAGTGGGGATGGCTATAATGGATTTGGTAATGATGGAAGCAATTTTGGAGGT GGTGGAAAGCTACAATGATTTTGGCAATTACAACAATCAGTCTTCAAATTTTGGACCCATG AAGGGAGGAAATTTTGGAGGCAGAAGCTCTGGCCCTATGGCGGTGGAGGCCAATACTTT GCAAAACCACAAAACCAAGGTGGCTATGGCGTTTCCAGCAGCAGCAGTAGCTATGGCAGT GGCAGAAGATTTTAA
Restriction Sites:	Please inquire
ACCN:	NM_001011724
Insert Size:	1000 bp



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001011724.1</u> , <u>NP_001011724.1</u>
RefSeq Size:	2365 bp
RefSeq ORF:	963 bp
Locus ID:	144983
UniProt ID:	<u>Q32P51</u>
Cytogenetics:	13q14.3
Protein Pathways:	Spliceosome
Gene Summary:	<p>Involved in the packaging of pre-mRNA into hnRNP particles, transport of poly(A) mRNA from the nucleus to the cytoplasm and may modulate splice site selection.[UniProtKB/Swiss-Prot Function]</p> <p>Transcript Variant: This variant (1) represents the longer transcript. Variants 1 and 2 encode the same protein.</p>