

Product datasheet for SC301740

FAM209A (NM_001012971) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	FAM209A (NM_001012971) Human Untagged Clone
Tag:	Tag Free
Symbol:	FAM209A
Synonyms:	C20orf106; dj1153D9.3
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC301740 representing NM_001012971. Blue=Insert sequence Red=Cloning site Green=Tag(s)

GCTCGTTTGTAGTAACCGTCAGAATTTTGTAAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTG
 GATCCGGTACCGAGGAGATCTGCCGCC**CGATCGCC**
 ATGTGGACGCTGAAATCGTCCCTGGTCTGCTTCTGTGCCTCACCTGCAGCTATGCCTTTATGTTCTCT
 TCTCTGAGACAGAAAAGTAGCGAAGGCTGGCTTGGGAGCAATACGGAGAGCACTTCGGATTCGGCAG
 AATCTACCAGAGCACACCAAGGCTGGCTTGGGAGCAATGGCTCTGGCTTCTTTTGTGTTGTGCCG
 TTTGTGATACTGCAGTGTCAAAGAGACAGTGAGAAGAATAAGGAGCAGAGTCCTCCTGGCCTTCGAGGC
 GGCCAACTTCACTCTCCATTAAGAAAAAAGAAATGCTTCCCCAACAAAGACTGTGCATTCAATACC
 TTAATGGAACGAGGTGGAGCTTATGAAATTTGTGTCCAAAGTGCAGGAACTTAAACGTGCCATGGCA
 ACAGGTAGTGGCAGTAACCTCAGGCTTCGAAAGTCAGAGATGCCTGCAGATCCATACCATGTCACGATC
 TGTGAAATATGGGGAGAAGAAAGCTCTAGCTGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
 TACAAGGATGACGACGATAAGGTTAAACGGCCGGC

Restriction Sites:	SgfI-MluI
ACCN:	NM_001012971
Insert Size:	516 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).


[View online »](#)

OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001012971.3</u>
RefSeq Size:	689 bp
RefSeq ORF:	516 bp
Locus ID:	200232
UniProt ID:	<u>Q5JX71</u>
Cytogenetics:	20q13.31
Protein Families:	Transmembrane
MW:	19.6 kDa