

Product datasheet for SC301745

PABPC1L2A (NM_001012977) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	PABPC1L2A (NM_001012977) Human Untagged Clone
Tag:	Tag Free
Symbol:	PABPC1L2A
Synonyms:	RBM32A
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC301745 representing NM_001012977. Blue=Insert sequence Red=Cloning site Green=Tag(s)

GCTCGTTTAGTGAACCGTCAGAATTTTGAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTG
 GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
 ATGGCCTCCCTGTACGTGGGCGACCTGCACCCTGAGGTGACCGAGGCAATGCTGTACGAGAAGTTCAGT
 CCAGCTGGGCCCATCTCTCCATCCGCATCTGCAGGGACAAGATCACCCGCCGCTCATTGGGCTACGCG
 TATGTCAACTACCAGCAACCGGTGGACGCCAAGCGGGCCCTGGAGACCCTGAACCTTGTATGTCATAAAG
 GGCAGGCCAGTGCATCATGTGGTCCCAGAGGGACCCGTCGCTCCGCAAGAGCGGGGTGGCAACGTC
 TTCATCAAGAACCTGGCAAGACCATCGACAACAAGGCGCTGTACAACATCTTCTCGGCGTTCGGCAAC
 ATCCTCTCCTGCAAAGTGGCCTGTGACGAAAAGGGGCCAAGGGCTACGGGTTCTGTCACCTTCCAAAAG
 CAGGAATCTGCGGAGCGGCCATCGATGTGATGAATGGCATGTTCTGAACTACCGCAAAATTTTCGTC
 GGGAGATTCAAGTCGCATAAAGAACGAGAGGCCGAAAGGGGAGCCTGGGCCAGGCAGTCCACTAGTGT
 GACGTCAAGGATTTTCAGGAAGACACCGACGAGGAGGCCACCTTGCGATGA
 AGCGGACCGACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGAT
 ATCCTGGATTACAAGGATGACGACGATAAGGTTTAA

Restriction Sites:	Sgfl-RsrII
ACCN:	NM_001012977
Insert Size:	603 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).


[View online »](#)

OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001012977.2</u>
RefSeq Size:	2175 bp
RefSeq ORF:	603 bp
Locus ID:	340529
UniProt ID:	<u>Q5JQF8</u>
Cytogenetics:	Xq13.2
MW:	22.8 kDa