

Product datasheet for **SC301781**

PRAMEF5 (NM_001013407) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	PRAMEF5 (NM_001013407) Human Untagged Clone
Tag:	Tag Free
Symbol:	PRAMEF5
Synonyms:	PRAMEF5L; PRAMEF23
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC301781 representing NM_001013407.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGAGCATCCGGACTCCACCCAGACTCCTGGAGCTTGCAGGGCGGAGCCTGCTGAGGGACCAAGCCTTG
GCCATGTCCACCCTGGAGGAGCTGCCACAGAACTTTCCCCCACTGTTTCATGGAGGCCTTCAGCAGG
AGACGCTGTGAGGCCCTGAAGCTGATGGTGCAGGCCTGGCCCTTCGCGCCGCTCCCTCTGAGGCCTCTG
ATAAAGATGCCTTGTCTGGAGGCCTTCCAAGCTGTGCTCGATGGGCTGGATGCACTGCTTACCCAAGGG
GTTTCATCCAGGAGGTGGAACCTTCAAGTGTGGATTTACAGGATGTCTGTGAGAACTTCTGGATGGTT
TGGTCTGAAGCTATGGCCATGGGTGCTTCTCAATGCCAAGAGGAACAAAAACCAGTGCAGGACTGT
CCAAGGATGAGAGGACAGCAGCCCTTACTGTGTTCTGTAAGCTTTGGCTCAAGAACAGGACTCTGGAT
GAATACCTCACCTGCCTCCTTCTATGGGTCAAGCAGAGGAAAGATTTACTACACCTGTGCTGTAAGAAG
CTGAAAAATTTGGGAATGCCCTTCGCAATATCAGAAGCATCCTGAAAAATGGTGAACCTAGACTGTATC
CAGGAGGTGGAAGTGAATTGCAAGTGGTACTGCCCATCTGACACAGTTTACCCCATACCTGGGCCAC
ATGAGGAATCTTCAGAAGCTCGTTCTCTCCACATGGATGTCTCTCGCTACGTTTCCCAGAGCAGAAG
AAGGAGATTGTTACCCAGTTCACCACTCAGTTCCTCAAGCTGTGCTGCCTCCAAAAGCTTCTATGAAC
TCTGTTTCTTTCTCGAAGGCCACCTGGACCAGCTGCTCAGCTGTCTGAAGACCTCGTTAAAGGTCTC
ACAATAACTAAGTGTGTGCTTTTGAATCAGACTTGAAGCATCTATCCCAGTGCCCGAGTATCAGTCAA
CTAAAGACCCTGGACCTGAGTGGCATCAGACTGACCAATTACAGTCTTGTGCCTCTCAAATTCCTTA
GAAAAAGTTGCAGCCACCCTTGAGTACCTGGATTTAGATGACTGTGGCATCATAGACTCCCAAGTCAAC
GCCATCCTGCCTGCCCTGAGCCGCTGCTTTGAGCTCAACACCTTCAGTCTCTGTGGAATCCCATCTCC
ATGGCCACCCTGGAGAACCTGCTGAGCCACACAATCATACTCAAAACTTATGCGTGGAGCTGTATCCT
GCCCCCGGAGAGTTATGATGCTGATGGTACTCTCTGCTGGAGCAGATTCTCTCAAATTAGGGCTGAG
CTGATGAAGAGAGTGAGGACTTAAGGCACCCCAAGAGGATCTTGTCTGTACTGACTGCTGCCCTGAC
TGTGGCAACAGGTCATTTTATGACCTGGAGGCAGATCAATGCTGCTGTGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: SgfI-MluI

ACCN: NM_001013407

Insert Size: 1431 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001013407.2](#)

RefSeq Size: 1855 bp

RefSeq ORF: 1431 bp

Locus ID: 343068

UniProt ID: [Q5TYX0](#)

Cytogenetics: 1p36.21

MW: 54.9 kDa