

## Product datasheet for **SC301832**

### PTX4 (NM\_001013658) Human Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	PTX4 (NM_001013658) Human Untagged Clone
Tag:	Tag Free
Symbol:	PTX4
Synonyms:	C16orf38
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)


**Fully Sequenced ORF:** >SC301832 representing NM\_001013658.  
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGGGAGTGGAAACTGGGAGGTCACAGGTCACCTTGTGGCAGTCGGTGTGAACGGCCTCCCGGACTC
CCACGCGGGTCCACGGACAGTCCGGAGTCTGTGGCCTGTCGGCCCTGCGCTGGTTCGGAGATTTC
CAGGAGGTGACCTGGACACACCTGCAGAACATCGCCAGCAACTACAACGTCTCTACAACGTTGAGGTC
CGTTCGGAGCCTGGCGAAGAGAGCCAGGCTGTGGCTCAGGCAGTCAACCGGTCACAGGCCTCGGTG
CAGGGGGAGCTGGCGCAGCTCAAGGCTGGGTGAGGAAGCTGCAGCGCCGAGGCCGAAAGTAGACACG
CGGCTGCGGGCCTTGACCTCACCTGGGCGAGAGGAGCCAGCAGCGGGCCCGGAAAGGAAGGCACAC
AAGGCCAGAGGGACGCCCTGCAGGACTCACTGGCACGCCTGGAGGGCCTCGTCCACAGCCAGGGCGCC
AGGCTGGTGTCTGGAGGGGGCTGCCGTGGCCACCCTGGCACTGCAGCCCTGGGCCGGCCCTG
GTCCAAACCCCAACCCAGCCTGAGGAGCTGGGCCCGACCTCCCTGAAGCTTCAGAGGGACAGGCAGGAG
CTCCGAGCTGCCTCTGAGCACAGGGGCCCCACAGGACTCCTCGGCCCTCTCAAGGGAGGGCGGAG
CCTCCAGCCTCAGGCAGCCATCGGGTACTCAGTGGGACTGCCCAAAAGACCCTCGGCAGAGGCATGG
TCCCCCAAGTGCCAGGAGAGATTTGCGCGTGGGCCCCACCCTCGTTTTCCCAAACGCCTCCACCAGG
AACGTGGTCTTCTCAGCCCTGGTTTCGTCACTGCCCTGCGAGCCCTGTCTTCTGCAGCTGGGTCCGC
ACGGCCTCCGGCCGCTGGGCACCCTCTGTCTACGCCACCGAGGACAATGACAACAAGCTGGTGTG
CACGGCCGAGACTCCCTGCTGCCGGATCCATCCACTTCGTGATCGGGGACCCGGCCTCAGGGAGCTG
CCCTTGCAGCTGCTGTGGACGGCCAGTGGCACCACATCTGTGTATCTGGACGTCCACCAGGGCAGG
TACTGGCTCCACGTGGATCGCAGGCTGGTGGCCACCGGCTCCCGCTTCAGGGAGGGCTATGAGATCCCC
CCCGGAGGGTCCCTCGTGTGGCCAGGAACAAGACAGCGTGGGGGGCGGATTGCAGACTCCGAGGCC
TTCGTGGGGAGCATGTCTGGCTTGGCTATCTGGGATCGGGCGCTGGTTCCCGGGGAAGTTGCAAACCTT
GCCATCGGGAAAGAGTTCCCGACAGGTGCCATCCTGACGCTGGCCAATGTGCACTAGCAGGCGGATTT
GTGCAGGGGGCAACTGCACCTGCCTGGAACGCTGTCCTGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

**Restriction Sites:** SgfI-MluI

**Plasmid Map:**



**ACCN:** NM\_001013658

**Insert Size:** 1422 bp

<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>OTI Annotation:</b>	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"> <li>1. Centrifuge at 5,000xg for 5min.</li> <li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li> <li>3. Close the tube and incubate for 10 minutes at room temperature.</li> <li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li> <li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li> </ol>
<b>RefSeq:</b>	<u><a href="#">NM_001013658.1</a></u>
<b>RefSeq Size:</b>	1422 bp
<b>RefSeq ORF:</b>	1422 bp
<b>Locus ID:</b>	390667
<b>UniProt ID:</b>	<u><a href="#">Q96A99</a></u>
<b>Cytogenetics:</b>	16p13.3
<b>MW:</b>	51.3 kDa
<b>Gene Summary:</b>	<p>This gene belongs to the pentraxin superfamily, whose members encode highly conserved multifunctional proteins. The encoded protein, like other members of this family, contains a conserved pentraxin domain at the C-terminus. The highest levels of expression of the protein were observed in bone marrow, small intestine and testes. [provided by RefSeq, Jun 2016]</p> <p>Transcript Variant: This transcript (2) uses an alternate 5' terminal exon and an alternate start codon, compared to variant 1. It encodes isoform 2, which has a shorter and distinct N-terminus, compared to isoform 1.</p>