

Product datasheet for SC301949

DDX19B (NM_001014451) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	DDX19B (NM_001014451) Human Untagged Clone
Tag:	Tag Free
Symbol:	DDX19B
Synonyms:	DBP5; DDX19; RNAh
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC301949 representing NM_001014451. Blue=Insert sequence Red=Cloning site Green=Tag(s)


```
GCTCGTTTGTAGTAACCGTCAGAATTTTGTAAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGCCACTGACTCATGGCCCTGGCGGTGGACGAGCAGGAAGCTGCGGCTGAGTCGTTGAGCAACTTG
CATCTTAAGGAAGAGAAAATCAAACCAGATACCAATGGTCTGTTGTCAAGACCAATGCCAATGCAGAG
AAGACAGATGAAGAAGAGAAAGAGGACAGAGCTGCCAGTCTTACTCAACAAGCTGATCAGAAGCAAC
CTTGTTGATAACACAAACCAAGTGAAGTCTGCAGCGGGATCCAACTCCCCTCTGTACTCGGTGAAG
TCTTTTGAAGAGCTTCGGCTCCACAGAACTTAATTGCCAATCTCAGTCTGGTACTGGTAAAACAGCT
GCCTTCGTGCTGGCCATGCTTAGCCAAGTAGAACCTGCAAACAAATACCCCAAGTGTCTATGTCTCTCC
CCAACGTATGAGCTCGCCCTCAAACAGGAAAAGTGATTGAACAAATGGGCAAATTTACCTGAACTG
AAGCTAGCTTATGCTGTTGAGGCAATAAATTGAAAAGAGGCCAGAAGATCAGTGAGCAGATTGTCATT
GGCACCCCTGGGACTGTGCTGGACTGGTCTCCAAGCTCAAGTTCATTGATCCCAAGAAAATCAAGGTG
TTTGTCTGGATGAGGCTGATGTCATGATAGCCACTCAGGGCCACCAAGATCAGAGCATCCGCATCCAG
AGGATGCTGCCAGGAACTGCCAGATGCTGCTTTTCTCCGCCACCTTTGAAGACTCTGTGTGGAAGTTT
GCCAGAAAGTGGTCCCAGACCCAAACGTTATCAAACCTGAAGCGTGAGGAAGAGACCCTGGACACCATC
AAGCAGTACTATGTCCTGTGCAGCAGCAGAGACGAGAAGTCCAGGCCTGTGTAACTCTACGGGGCC
ATCACCATTGCTCAAGCCATGATCTTCTGCCATACTCGCAAACAGCTAGTTGGCTGGCAGCAGAGCTC
TCAAAGAAGGCCACCAGGTGGCTCTGCTGAGTGGGGAGATGATGGTGAACAGAGGGCTGCAGTGATT
GAGCGCTTCCGAGAGGGCAAAGAGAAGTTTTGGTGACCACCAACGTGTGTGCCCGCGGCATTGATGTT
GAACAAGTGTCTGTCATCAACTTTGATCTTCCCGTGGACAAGGACGGGAATCCTGACAATGAGACC
TACCTGCACCGGATCGGGCGCACGGCCGCTTTGGCAAGAGGGCCCTGGCAGTGAACATGGTGGACAGC
AAGCACAGCATGAACATCCTGAACAGAATCCAGGAGCATTTTAATAAGAAGATAGAAGATTGGACACA
GATGATTTGGACGAGATTGAGAAAATAGCCAAC TGA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGCGC
```



[View online »](#)

Restriction Sites: SgfI-MluI

Plasmid Map:



ACCN: NM_001014451

Insert Size: 1347 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001014451.2](#)

RefSeq Size: 1736 bp

RefSeq ORF: 1347 bp

Locus ID: 11269

UniProt ID: [Q9UMR2](#)

Cytogenetics: 16q22.1

MW: 50.5 kDa

Gene Summary: DEAD box proteins, characterized by the conserved motif Asp-Glu-Ala-Asp (DEAD), are putative RNA helicases. They are implicated in a number of cellular processes involving alteration of RNA secondary structure such as translation initiation, nuclear and mitochondrial splicing, and ribosome and spliceosome assembly. Based on their distribution patterns, some members of this family are believed to be involved in embryogenesis, spermatogenesis, and cellular growth and division. This gene encodes a DEAD box protein, which exhibits RNA-dependent ATPase and ATP-dependent RNA-unwinding activities. This protein is recruited to the cytoplasmic fibrils of the nuclear pore complex, where it participates in the export of mRNA from the nucleus. Multiple alternatively spliced transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jul 2008]

Transcript Variant: This variant (2) lacks an in-frame exon in the coding region, as compared to variant 1. The encoded isoform (2) lacks an internal segment, as compared to isoform 1.