

Product datasheet for SC301972

BANF2 (NM_001014977) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	BANF2 (NM_001014977) Human Untagged Clone
Tag:	Tag Free
Symbol:	BANF2
Synonyms:	BAF-L; BAF2; BAFL; C20orf179
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC301972 representing NM_001014977. Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTGTAAACGACTCACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGCATCGCC
ATGGACAACATGTCTCCAGGCTGAGAGCCTTCCTCTCCGAACCCATTGGAGAAAAGGATGTCTGCTGG
GTGGATGGCATCAGCCATGAGCTCGCGATCAATTTGGTACCAAAGGTATCAATAAGGCCTACATCCTG
CTGGGACAATTCCTTCTGATGCACAAGAATGAAGCCGAGTTTCAGAGGTGGCTCATTGTGTTTTGGT
GCCACTGAGTGTGAGGCCAGCAGACTTCTCACTGCCTCAAGGAGTGGTGTGCCTGCTTCTCTGTAG
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGCGC
```

Restriction Sites:	Sgfl-Mlul
ACCN:	NM_001014977
Insert Size:	273 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.



[View online »](#)

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001014977.3</u>
RefSeq Size:	466 bp
RefSeq ORF:	273 bp
Locus ID:	140836
UniProt ID:	<u>Q9H503</u>
Cytogenetics:	20p12.1
MW:	10.3 kDa
Gene Summary:	May play a role in BANF1 regulation and influence tissue-specific roles of BANF1. [UniProtKB/Swiss-Prot Function] Transcript Variant: This variant (2) differs in the 5' UTR compared to variant 1. Both variants 1 and 2 encode the same isoform (1).