

Product datasheet for SC302012

IGSF11 (NM_001015887) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	IGSF11 (NM_001015887) Human Untagged Clone
Tag:	Tag Free
Symbol:	IGSF11
Synonyms:	BT-IgSF; CT119; CXADRL1; Igsf13; VSIG3
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC302012 representing NM_001015887. Blue=Insert sequence Red=Cloning site Green=Tag(s)


```
GCTCGTTTGTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGACTTCTCAGCGTTCCCTCTGGCGCCTTTGCTGCTCCTCTCTGCACGGTGTTCAGCATCCCTG
GAAGTGTGAGAGGCCCTGGGAGTATCCAGGTGGCCCGGGTCCAGCCAGCAGTCTGCCCTGCACCTTTC
ACTACCAGCGTGCCTCATTAACTCAATGTCATTTGGATGGTCACTCCTCTCCAATGCCAACCAA
CCTGAACAGGTATCCTGTATCAGGGTGGACAGATGTTTGGATGGTCCCGGGTCCACGGTAGGGTA
GGATTTACAGGCACCATGCCAGTACCAATGTCTATCTTCAATAAACAACACTCAGTTATCAGACACT
GGCACCTACCAGTGCCTGGTCAACAACCTTCCAGACATAGGGGGCAGGAACATTGGGGTCCCGGTCTC
ACAGTGTAGTTCCCTTCTGCCCCACTGCCAAATCCAAGGATCCCAGGATATTGGCAGCGATGTC
ATCCTGCTCTGTAGTCAAGGAAGGCATTCTCGACCAACTTACCTTTGGGAGAAGTTAGACAATACC
CTCAAACACTCCAACAGCTACTCAGGACCAGGTCCAGGGAACAGTACCATCCGGAACATCAGTGCC
CTGTCTCAGGTTTGTACCAGTGCCTGCTTCTAATGCTATTGGAACAGCACCTGTCTTCTGGATCTC
CAGGTTATTTACCCAGCCAGGAACATTGGACTAATAGCTGGAGCCATTGGCACTGGTGCAGTTATT
ATCATTTTTGCATTGCACTAATTTAGGGCATTCTTTACTGGAGAAGCAAAAATAAAGAGGAGGAA
GAAGAAGAAATTCCTAATGAAATAAGAGAGGATGATCTCCACCAAGTGTCTTCTGCCAAAGCATT
CACACTGAGATTTCTCCTCGGACAACAACACTAACCTTTCCAATGCCTACAACAGTCGATACTGG
AGCAACAATCCAAAAGTTCATAGAAACACAGAGTCAGTCAGCCACTTCAGTGACTTGGGCAATCTTTC
TCTTCCACTCAGGCAATGCCAACATACCATCCATTTATGCTAATGGGACCCATCTGGTCCCGGGTCAA
CATAAGACTCTGGTAGTGACAGCCAACAGAGGGTATCACCACAGGTGATGTCCAGGAGCAATGGCTCA
GTCAGTAGGAAGCCTCGGCCTCCACACTCATTCTACACCATCAGCCACGCAACACTGGAACGAATT
GGTGCAGTACCTGTATGGTACCAGCCAGAGTCGGGCCGGTCTTGGTATAG
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```



[View online »](#)

Restriction Sites: SgfI-MluI

Plasmid Map:



ACCN: NM_001015887

Insert Size: 1296 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001015887.1](#)

RefSeq Size: 3446 bp

RefSeq ORF: 1296 bp

Locus ID: 152404

UniProt ID: [Q5DX21](#)

Cytogenetics: 3q13.32

Protein Families: Transmembrane

MW: 46.1 kDa

Gene Summary: IGSF11 is an immunoglobulin (Ig) superfamily member that is preferentially expressed in brain and testis. It shares significant homology with coxsackievirus and adenovirus receptor (CXADR; MIM 602621) and endothelial cell-selective adhesion molecule (ESAM).[supplied by OMIM, Apr 2005]
Transcript Variant: This variant (2) differs in the 5' UTR and coding region compared to variant 1. The resulting isoform (b) is one aa longer and has a distinct N-terminus compared to isoform a.