

Product datasheet for **SC302043**

USP1 (NM_001017415) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	USP1 (NM_001017415) Human Untagged Clone
Tag:	Tag Free
Symbol:	USP1
Synonyms:	UBP
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)




[View online »](#)

Fully Sequenced ORF: >SC302043 representing NM_001017415.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAATACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGGATCGCC
ATGCCTGGTGCATACCTAGTGAAGTAATGGACTTTCAAGAGGTAGCCCTTCAAAGAAAAACAGACTT
TCCTTAAAGTTTTTTTCAGAAAAAGGAACTAAGAGAGCTTTGGATTTACAGATTCTCAAGAAAAATGAA
GAAAAAGCTTCTGAATATAGAGCATCTGAAATTGATCAAGTTGTTCTGCAGCACAGTCTTCACCTATA
AACTGTGAGAAGAGAGAAAACTTGTACCATTTGTGGACTGAATAATCTCGGCAATACTTGCTATCTT
AATAGTATACTTCAGGTATTATATTTTTGTCCCGTTTTAAATCTGGAGTAAAGCACTTATTTAATATT
ATTTCAAGGAAGAAAGAAGCTCTAAAGGATGAAGCCAATCAAAAAGACAAGGGAAATTGCAAGAAGAT
TCTTTGGCAAGTTATGAATTGATATGCAGTTTACAGTCTTAATCATTTCGGTTGAACAGCTCCAGGCT
AGTTTTCTCTAAATCCAGAGAAATATACTGATGAACTTGCCACTCAGCCAAGGCGACTGCTTAACACA
CTGAGGGAACCAACCCTATGTATGAAGGATATCTACAGCATGATGCACAGGAAGTATTACAATGTATT
TTGGGAAACATTCAGAAACATGCCAATCCTAAAAAAGAAGAAGTAAAAATGTGGCAGAATTACCT
ACTAAGGTAGAAGAAATACCTCATCCGAAAGAGGAAATGAATGGTATTAACAGCATAGAGATGGACAGT
ATGAGGCATTCTGAAGACTTTAAAGAGAACTCCCAAAAGGAAATGGGAAAAAGAAAAAGTACACTGAA
TTTGGTAAACATGAAGAAAAAGTTAAATATCAAGGAACACCAGTCAATTGGAAGAGAACCAGAGACAA
ACTAGATCAAAAAGAAAAGCTACAAGTGATACATTAGAGAGTCCCTCTAAAAATAATTCCTCAAGTATATT
TCTGAAAAATGAGAGTCCAAGACCCTCACAAAAGAAATCAAGAGTTAAAATAAATTGGTTAAAGTCTGCA
ACTAAGCAACCCAGCATTCTTTCTAAATTTTGTAGTCTGGGAAAAATAACAACAAACCAAGGAGTCAAA
GGCAACTCTAAAGAAAAATGAATGTGATCCTGAAGAGGACTTGGGGAAGTGTGAAAGTATAACACAAT
AATGGTTGTGGACTTGAATCTCCAGGAACTGTTACACCTGTAATGTTAATGAAGTTAAACCCATA
AACAAAGGTGAAGAACAATTGGTTTTGAGCTAGTGGAGAAATTTTCAAGGTACAGCTGGTATTAAGG
ACGCGTTGCTTGAATGTGAAAGTTAACAGAAAGAGAAGATTTTCAAGACATCAGTGTGCCAGTA
CAAGAAGATGAGCTTTCCAAAGTAGAGGAGAGTCTGAAATTTCTCCAGAGCCAAAAACAGAAATGAAG
ACCCTGAGATGGGCAATTCACAATTTGCTTCAGTAGAAAGGATTGTAGGAGAAGATAAATTTTCTGT
GAAAACCTGCCATCATTATACTGAAGCTGAACGAAGTCTTTGTTTGACAAAATGCCTGAAGTTATAACT
ATTCATTTGAAGTGCTTTGCTGCTAGTGGTTTGGAGTTTGATTGTTATGGTGGTGGACTTCCAAGATC
AACACTCCTTTATTGACACCTCTAAATTTGCTACTAGAAGAATGGAGCACAAAGCCAACCTAACGACAGC
TATGGATTATTTGCGGTTGTGATGCATAGTGGCATTACAATTAGTAGTGGCATTACACTGCTTCTGTT
AAAGTCACTGACCTTAACAGTTTGAAGTATAGATAAAGGAAATTTTGTGGTTGACCAAATGTGTGAAATA
GGTAAGCCAGAACCATTGAATGAGGAGGAAGCAAGGGGTGGTTGAGAATTATAATGATGAAGAAGTG
TCAATTAGAGTTGGTGGAAATACACAGCCAAGTAAAGTTTTGAACAAAAAAATGTAGAAGCTATTGGA
CTTCTTGGAGGACAAAAGAGCAAAGCAGATTATGAGCTATACAACAAAGCCTCTAATCCTGATAAGGTT
GCTAGTACAGCGTTTGTGAAAAAGAAATCTGAGACTAGTACTACTGGGACCCATGAATCTGAT
AGAAACAAGGAATCCAGTGACCAACAGGCATTAATATTAGTGGATTTGAGAACAAAATTCATACGTA
GTGCAAAGCTTAAAGGAGTATGAGGGGAAGTGGTTGCTTTTTGATGATTCTGAAGTCAAAGTTACTGAA
GAGAAGGACTTTCTGAATCTCTTTCCCTTCTACATCTCCTACTCTACTCCTACTTCTACTTATTTAT
AAGAAATTAAG
AGCGGACCGACGCGTACGCGGCCCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGAT
ATCCTGGATTACAAGGATGACGACGATAAGGTTTAA
  
```

Restriction Sites: Sgfl-RsrII

Plasmid Map:


ACCN: NM_001017415

Insert Size: 2358 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001017415.1](#)

RefSeq Size: 3519 bp

RefSeq ORF: 2358 bp

Locus ID: 7398

UniProt ID: [O94782](#)

Cytogenetics: 1p31.3

Protein Families: Protease

MW: 88.2 kDa

Gene Summary: This gene encodes a member of the ubiquitin-specific processing (UBP) family of proteases that is a deubiquitinating enzyme (DUB) with His and Cys domains. This protein is located in the cytoplasm and cleaves the ubiquitin moiety from ubiquitin-fused precursors and ubiquitinated proteins. The protein specifically deubiquitinates a protein in the Fanconi anemia (FA) DNA repair pathway. Alternate transcriptional splice variants have been characterized. [provided by RefSeq, Jul 2008]
Transcript Variant: This variant (2) uses an alternate upstream transcriptional start site and contains an alternate exon in the 5' UTR, compared to variant 1. Variants 1, 2 and 3 encode the same protein.