

Product datasheet for **SC302062**

ARHGAP8 (NM_001017526) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ARHGAP8 (NM_001017526) Human Untagged Clone
Tag:	Tag Free
Symbol:	ARHGAP8
Synonyms:	BPGAP1; PP610
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)


Fully Sequenced ORF: >SC302062 representing NM_001017526.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAATACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGCTGGCCAGGATCCTGCGCTGAGCACGAGTCACCCGTTCTACGACGTGGCCAGACATGGCATTCTG
CAGGTGGCAGGGGATGACCGCTTTGGAAGACGTTGTACAGTTCAGCTGCTGCCGATGCCACCCTCC
CACGAGCTGGACCACCAGCGGCTGCTGGAGTATTTGAAGTACACACTGGACCAATACGTTGAGAAGCAT
TATACCATCGTCTATTTCCACTACGGGCTGAACAGCCGGAACAAGCCTTCCCTGGGCTGGCTCCAGAGC
GCATACAAGGAGTTCGATAGGAAAGACGGGATCTCACTATGTGGCCAGGCTGGTCTCGAACTCCAAG
CTCAAGCGATCTCCACCTCAGCCTCCCAAAGTACTGGGATTACAGGTACAAGAAGAACTTGAAGGCC
CTCTACGTGGTGACCCACCAGCTTCATCAAGGTCTGTGGAACATCTTGAAGCCCTCATCAGTCAC
AAGTTTGGGAAGAAAGTCATCTATTTCACTACCTGAGTGAGCTCCACGAACACCTTAAATACGACCAG
CTGGTCATCCCTCCGAAGTTTTCGGTACGATGAGAAGCTCCAGAGCCTGCACGAGGGCCGGACGCCG
CCTCCCAAGACACCACCGCCGCGGCCCGCTGCCACACAGCAGTTTGGCGTCAGTCTGCAATAC
CTCAAAGACAAAATCAAGGCAGACTATCCCCCTGTGCTGAGGTTACAGTGAGCTACCTGAGAGAG
AAAGGCCTGCGCACCGAGGGCCTGTTCCGGAGATCCGCCAGCGTGCAGACCGTCCGCGAGATCCAGAGG
CTCTACAACCAAGGGAAGCCCGTGAACCTTGTGACTACGGGACATTACATCCCTGCCGTGATCCTG
AAGACCTTCTGCGAGAGCTGCCCCAGCCGCTTCTGACCTTCCAGGCCTACGAGCAGATTCTCGGGATC
ACCTGTGTGGAGAGCAGCCTGCGTGTACTGGCTGCCGCCAGATCTTACGGAGCCTCCAGAGCACAAAC
TACGTCGTCTCCGCTACCTCATGGGCTTCTGCATGCGGTGTCCCGGAGAGCATCTTCAACAAAATG
AACAGCTAACCTGGCCTGTGTCTTCCGGCTGAATTTGATCTGGCCATCCAGGGGTCTCCTCCCTG
AGTGCCCTGTGCCCTGAACATGTTCACTGAAGTCTGATCGAGTACTATGAAAAGATCTTCAGCACC
CCGGAGGCACCTGGGAGCACGGCTGGCACCATGGGAACAGGGGAGCAGGGCAGCCCTTTCAGGAG
GCTGTGCCACGGACACAAGCCACGGCCTCACCAAGCCTACCCTACCTCCGAGTCCCTGATGGCAGCC
AGAAGACGTCTCTAG
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: SgfI-MluI

Plasmid Map:



ACCN: NM_001017526

Insert Size: 1395 bp

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001017526.1
RefSeq Size:	1725 bp
RefSeq ORF:	1395 bp
Locus ID:	23779
UniProt ID:	P85298
Cytogenetics:	22q13.31
MW:	53.5 kDa
Gene Summary:	<p>This gene encodes a member of the RHOGAP family. GAP (GTPase-activating) family proteins participate in signaling pathways that regulate cell processes involved in cytoskeletal changes. GAP proteins alternate between an active (GTP-bound) and inactive (GDP-bound) state based on the GTP:GDP ratio in the cell. This family member is a multidomain protein that functions to promote Erk activation and cell motility. Alternative splicing results in multiple transcript variants. Read-through transcripts from the upstream proline rich 5, renal (PRR5) gene into this gene also exist, which led to the original description of PRR5 and ARHGAP8 being a single gene. [provided by RefSeq, Nov 2010]</p> <p>Transcript Variant: This variant (1) represents the longest transcript and encodes the longest isoform (1).</p>