

## Product datasheet for **SC302077**

### LRIT2 (NM\_001017924) Human Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	LRIT2 (NM_001017924) Human Untagged Clone
Tag:	Tag Free
Symbol:	LRIT2
Synonyms:	LRRC22
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

**Fully Sequenced ORF:** >SC302077 representing NM\_001017924.  
 Blue=Insert sequence Red=Cloning site Green=Tag(s)


```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGCTTCAGTTTTTTCATTABlueCTCTCTGTTAGTTCTGGTCTTTCTGGATACACACGCGAGCTCAGCCTTTC
TGCTGCCAGGATGCACBlueTGCTCAGAGGAGATTTTGGCAGGACTCTGCAGTGCACATCTGTCTCCTTG
GAAAGATCCCTGGAACTTTCTGAAGAGTTCAAGCAAGTGAAGTGAATGAAAATTCACCCTATTTGAG
ATGCCCAAGGGTCTTTCATCAACATGAGCACCTTGAATACCTCTGGCTCAATTTAACAATATCAGT
GTGATCCACCTAGGAGCCCTGGAACACCTGCCAGAACTGAGGGAGCTGAGACTGGAGGGGAACAAGCTC
TGCTCAGTACCATGGACAGCGTTCCTGCCACCCCTCTCCTGAGGGTCTTGGATCTCAAACGCAACAAG
ATTGATGCACTCCCTGAGCTGGCTCTTCAATTCTTGGTCAGCCTGACCTACCTTGACCTATCCTCCAAT
AGGCTTACAGTTGTATCCAAGAGTGTCTCCTGAACTGGCCAGCCTACCAGAAATGCCGGCAGCCTGAC
TGTGGGGCTGAGATTCTCTCCAGCCTGGTGGTGGCCCTGCATGACAACCCCTGGGTATGTGACTGTCCG
CTAAGGGGGCTTGTCCAGTTTGTCAAGTCCATTACCCTCCAGTCATCCTGGTGAATTCCTACCTGATA
TGTGAGGGCCCTCTGTCCAAGGCAGGGCAGCTTTTTCATGAAACTGAGCTTAGTGCTTGCATGAAGCCA
CAGATCTCAACCCCCAGTGCCAATATCACCATCCGGGCAGGACAGAATGTGACCCCTGCGATGCTTGCA
CAGGCCAGCCCCTCACCATCCATTGCATGGACTTATCCCCTGAGTATGTGGAGAGAATTTGATGTGTTG
ACATCTTCTACTGGAGAAGACACTGCTCTGTGACAGCTGGCCATACCTGCTGCCACCTGGTAGACAGT
GGTAATTACACCTGCATGGCCTCCAACCTCCATTGGCAAGAGCAACCTTGTAACTCTCTCCATGTCCAG
CCTGCCAGGCCCTACATGCACCTGATTCTCTTCCATCCCCTCGGAGGGCAATGCCTACATTGACCTG
CGGTTGTCAAGCAGACAGTGCATGGGATTTTGTGGAGTGGCTTGCAGTGGCTGACACCTTAAGGAG
GAGTGGTTCACCCTCTACATTGCATCGGATGAAGCCTTCAGGAAGGAGGTGGTTCACATTGGCCCCGGA
ATCAATACTTATGCTGTGGATGACCTCCTTCTGGCACAAAATATGAGGCCTGCCTCAGCCTAGAGGGC
CAGCCTCCACACCCAGGGCAGTGTGTAGCTTTTGTAAACAGGCAGAGATGCTGGTGGGCTAGAGGCACGT
GAGCACCTCCTGCATGTCACAGTGGTCTGTGTGGTGTGCTTGCAGTGCCTGTGGGCGCCTATGCC
TGGGCAGCCCAGGGCCCCTGCAGCTGCAGCAAGTGGTCTGCGCGGCTGTCTTCATCGCAGGAAAGCC
CCCAGCTGCACCCCTGCAGCCCCGAGTCCAAGGATGGCTCCTTTAGAGAACATCCAGCTGTCTGTGAT
GACGGTGAAGGGCACATAGACACTGAGGGGGACAAGGAGAAAGGAGGAACGGAAGACAACAGCTGA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

**Restriction Sites:**

Sgfl-MluI

**Plasmid Map:**



<b>ACCN:</b>	NM_001017924
<b>Insert Size:</b>	1653 bp
<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>OTI Annotation:</b>	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"> <li>1. Centrifuge at 5,000xg for 5min.</li> <li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li> <li>3. Close the tube and incubate for 10 minutes at room temperature.</li> <li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li> <li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li> </ol>
<b>RefSeq:</b>	<u><a href="#">NM_001017924.3</a></u>
<b>RefSeq Size:</b>	3157 bp
<b>RefSeq ORF:</b>	1653 bp
<b>Locus ID:</b>	340745
<b>UniProt ID:</b>	<u><a href="#">A6NDA9</a></u>
<b>Cytogenetics:</b>	10q23.1
<b>Protein Families:</b>	Transmembrane
<b>MW:</b>	60.2 kDa