

Product datasheet for **SC302093**

MARVELD3 (NM_001017967) Human Untagged Clone

Product data:


Product Type:	Expression Plasmids
Product Name:	MARVELD3 (NM_001017967) Human Untagged Clone
Tag:	Tag Free
Symbol:	MARVELD3
Synonyms:	MARVD3; MRVLDC3
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC302093 representing NM_001017967. Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTACTATAGGGCCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGAAGATCCGTGCGGGGCTCGCGAGCCCGGGCCCGGCCGAGAGAGCGGGACCCGGGACGGCGCC
CACCCAGACCAAGGCCGACCCACGATCGACCGGGACCGACCCGGGGACCCGCGCAGGAAGCGAAGC
AGCGACGGGAACCGGCAAGGGACGGGACCGGGACCCGGAGAGAGACCAGGAGAGGGACGGGAACCGC
GACCGGAACCGGGACCGGGAGAGGGAGAGAGAGAGGGAAAGAGACCCGGACCGAGGCCCCCGCGGGAC
ACACACAGGGACCGGGCCCTCGCGCAGGTGAACACGGAGTTGGGAAAAACCGGCCAAAGCCGGACG
CGGGACGGAGCCCGGGACTGACCTGGGACGCAGCCGCGCCTCTGGGCCCGCGCCTGGGAAGCCCG
GAGCCCGCCGAGCCGAGAGGAAGGGAGACCCCGGGCCGCGCAGACCCGAAAGTGAACCCCTTCGGAG
AGATATCTGCCCTCGACCCCGAGCCTGGACGAGAGGAGGTGGAATATTACAGTCAGAGCCGGAAGGA
CTCCTGGAATGCCACAAATGCAAACTACTTGTGCACTGGGAGAGGTGTGGTGCAGATAGTGGAGGTGGTC
TTGAATGGGATGGTTCTCATATGCATCGTGGCCTCCTACTTTGTCTTGGCGGATTCAGTGCCAGCTTT
TCCAGCGCGGTGGCTTTGGGAACAATACTACTACCCTTCGAGGGCACTGAGCTAGAGCAGGTTCCG
CAGCTGGACCAGCAGTACACGATCCTCCGCTCGCCCTGATACGCGTGGCGTGGCTGTTTCTCTGGGT
CTGGGGTCTCACCATGGGTGTTTTACTCCAAGGAGCCAAGAGTCAACAATGTTGTCAGGGAAGTGG
CTCCTCACGGAGCCGCTTACGCTCCTAGCGCAGTGGCTACTGCACAGGCATTGGTGTTCACCTC
CAGTGGCTCTGCAGATCAATAGCACCGCACTTGCAAAACAAGAGAGAGGCTCTATGCCCGCAAGGGT
CTCACCTGGATGGACTGCCAGCTGGCAGGCACCGAGCAGCAGCCACCTTTGCTTGTCTTCTCGTG
ATCATGTACGGCGCCAGCGTGGTGTGGCCTGCGTAGCTACCGAGAACAAGCGCTACAAAGGCAGC
CGAGAACAGCCCGAAGTTACAGTGTGCACCGGAATATCTGTGGTCTGGAACCTTTGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTAAACGGCCGGC
```

Restriction Sites: Sgfl-MluI



[View online »](#)

Plasmid Map:


ACCN: NM_001017967

Insert Size: 1233 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001017967.3](#)

RefSeq Size: 2273 bp

RefSeq ORF: 1233 bp

Locus ID: 91862

UniProt ID: [Q96A59](#)

Cytogenetics: 16q22.2

Protein Families: Druggable Genome, Transmembrane

MW: 46 kDa

Gene Summary: As a component of tight junctions, plays a role in paracellular ion conductivity.
[UniProtKB/Swiss-Prot Function]
Transcript Variant: This variant (1) represents the longer transcript and encodes the longest isoform (1).