

## Product datasheet for **SC302158**

### **CORO1B (NM\_001018070) Human Untagged Clone**

#### **Product data:**

Product Type:	Expression Plasmids
Product Name:	CORO1B (NM_001018070) Human Untagged Clone
Tag:	Tag Free
Symbol:	CORO1B
Synonyms:	CORONIN-2
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

**Fully Sequenced ORF:** >SC302158 representing NM\_001018070.  
 Blue=Insert sequence Red=Cloning site Green=Tag(s)


```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGCATCGCC
ATGTCCTTCCGCAAAGTGGTCCGGCAGAGCAAATTCGGCATGTGTTCCGGCAGCCGGTCAAGAACGAC
CAGTGCTATGAGGACATTCGCGTGTCCCGTGTACCTGGGACAGCACCTTCTGCGCCGTCAACCCCAAG
TTCTTGGCGGTGATTGTGGAGGCCAGTGGAGGGGGTGCCTTTCTGGTGTCCCCCTAAGCAAGACGGGC
CGCATTGACAAGGCCTACCCGACGGTGTGTGGGCACACGGGACCTGTCTGGACATCGACTGGTGTCTT
CACAACGACGAAGTCATAGCCAGCGGCTCGGAGGACTGCACGGTCATGGTGTGGCAGATCCAGAGAAC
GGGCTGACCTCCCGCTGACAGAGCCGGTGGTGGTACTGGAGGGGCACACCAAGCGAGTGGGCATCATC
GCCTGGCACCCACGGCCGAAACGTGTGCTCAGTGCAGGCTGCGACAACGTGGTACTCATCTGGAAAT
GTGGGCACAGCGGAGGAGCTGTACCGCTGGACAGCCTGCACCCTGACCTCATCTACAATGTCAGCTGG
AACCACAATGGCAGCCTGTTTTGCTCAGCATGCAAGGACAAGAGCGTGGCATCATCGACCCCGTCCGG
GGCACCCCTGGTGGCAGAGCGGGAGAAGGCTCATGAGGGGGCCCGCCCATGCGGGCCATCTTCTGGCA
GATGGCAAGGTGTTACCCACAGGCTTACGCCAATGAGCGAGCGGAGCTGGCGCTCTGGGACCCAGAA
AACCTCGAGGAACCCATGGCCCTGCAGGAAGTGGACTCGAGCAACGGGGCCCTGTGCCCTTCTACGAC
CCCAGACACAGTGTGGTCTACGTCTGCGCAAGGGTGGACTCCAGCATCCGGTACTTTGAGATCACAGAG
GAGCCTCCCTACATCCACTTCTGAACAGTTCACCAGCAAGGAGCCGACGCGGGGTATGGGCAGCATG
CCCAAGCGGGGCTGGAGGTGAGCAAGTGCAGATCGCCCGTCTACAAACTGCATGAGCGCAAGTGT
GAGCCCATCGTCATGACTGTGCCAAGAAAGTCGGACCTTCCAGGATGATCTGTACCCGACACAGCC
GGGCCCCGAGGCAGCCCTGAGGCTGAGGAGTGGGTGAGCGGGCGGGATGCCGACCCGATCCTCATCTCA
CTGCGGGAGGCCCTACGTGCCAGCAAGCAGCGGGACCTGAAGATCAGCCGGCGCAACGTGTTGTCTGAC
AGCCGGCCCGCCATGGCCCGGGCTCTCCACCTAGGGGCCCGCCCTCCACCACACTGCTGCTGAT
GCCACCCCGAGCGCAGCCTGGCCAGAGCCGGGAGGCTGGGAAGCTGGAGGAGGTGATGCAAGGAGCTG
CGGGCCCTGAGGGCGCTGGTCAAGGAGCAGGGCGACCGCATCTGCCGCTGGAGGAGCAGCTGGCCCGC
ATGGAGAACGGGGATGCGTAG
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

**Restriction Sites:**

SgfI-MluI

**Plasmid Map:**



**ACCN:**

NM\_001018070

<b>Insert Size:</b>	1470 bp
<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>OTI Annotation:</b>	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_001018070.2</a>
<b>RefSeq Size:</b>	1949 bp
<b>RefSeq ORF:</b>	1470 bp
<b>Locus ID:</b>	57175
<b>UniProt ID:</b>	<a href="#">Q9BR76</a>
<b>Cytogenetics:</b>	11q13.2
<b>MW:</b>	54.2 kDa
<b>Gene Summary:</b>	<p>Members of the coronin family, such as CORO1B, are WD repeat-containing actin-binding proteins that regulate cell motility (Cai et al., 2005 [PubMed 16027158]).[supplied by OMIM, Mar 2008]</p> <p>Transcript Variant: This variant (2) represents the longer transcript. Both variants 1 and 2 encode the same protein.</p>