

## Product datasheet for **SC302182**

### MYZAP (NM\_001018100) Human Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	MYZAP (NM_001018100) Human Untagged Clone
Tag:	Tag Free
Symbol:	MYZAP
Synonyms:	GCOM1; Gup; MYOZAP
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)


**Fully Sequenced ORF:** >SC302182 representing NM\_001018100.  
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAATACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGCTGCGCTCCACGTCCACGGTCACCCCTGCTCTCGGGCGGCCCGCCAGGACGCCCGGGCGCCAGC
AGGAGGGCAAATGTTTGCAGACTACGGCTGACCTCCTGAGAGTCCAGTTCCTGAGCAATGTGAA
AAGAAGATTGAGAGAAAAGAGCAGCTTCTTGACCTGAGCAATGGAGAACCCTACCAGGAACTTCCTCAG
GGTGTGTTTATGGTGTGGTGCAGAGATCAGATCAAAATCAGCAGAAAAGAAATGGTGGTGTATGGGTGG
TCCACCAGTCAGCTGAAAGAAGAGATGAACTACATCAAAGATGTGAGAGCCACTTTGAAAAGGTGAGA
AAGCGAATGTATGGAGACTATGATGAGATGAGACAGAAGATTCGACAGCTCACCCAGGAACTATCAGTT
TCCCATGCTCAGCAGGAGTATCTGGAGAATCACATCCAAACCCAGTCGTCTGCCCTGGATCGTTTTAAT
GCCATGAACTCAGCCTTGGCATCAGATTCCATTGGCCTGCAGAAAACCCTCGTGGATGTGACTTTGGAA
AACAGCAACATTAAGGATCAAATCAGAAATCTGCAGCAGACGTATGAAGCATCCATGGACAAGCTGAGG
GAAAAGCAGAGGCAGTTGGAGGTAGCGCAAGTTGAAAACCAGCTGCTAAAATGAAGGTGGAATCGTCC
CAAGAAGCCAATGCTGAGGTGATGCGAGAGATGACCAAGAAGCTGTACAGCCAGTATGAGGAGAAGCTG
CAGGAAGAACAGAGGAAGCACAGTGTGAGAAGGAGGCTCTTTTGGAAAGAACCAATAGTTTTCTGAAA
GCGATTGAAGAAGCCAATAAAAAGATGCAAGCAGCAGAGATCAGCCTAGAGGAGAAAAGACCAGAGGATC
GGGAGCTGGACAGGCTGATTGAGCGCATGGAAAAGGAACGTCATCAACTGCAACTTCAACTCCTAGAA
CATGAAACAGAAATGTCTGGGGAGTTAACTGATTCTGACAAGGAAAGGTATCAGCAGTTGGAGGAGGCA
TCAGCCAGCCTCCGTGAGCGGATCAGACACCTAGATGACATGGTGCATTGCCAGCAGAAGAAAGTCAAG
CAGATGGTCGAGGAGATTGAATCATTAAAGAAAAAGTTGCAACAGAAAACAGCTCTTAATACTCAGCTT
TTAGAAAAGATACTTTTCTAGAAGGAGAGAATAATGAACTACAAGCAGGTTGGACTATTTAACAGAA
ACCCAGGCCAAGACCGAAGTGGAAACCAGAGAGATAGGAGTGGGCTGTGATCTTCTACCCAGCCAACA
GGCAGGACTCGTGAAATTGTGATGCCTTCTAGGAACTACACCCATACACAAGAGTCTGGAGTTAACC
ATGAAGAAAACCTGACTTAG
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

**Restriction Sites:** SgfI-MluI

**Plasmid Map:**



**ACCN:** NM\_001018100

**Insert Size:** 1401 bp

<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>OTI Annotation:</b>	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_001018100.4</a>
<b>RefSeq Size:</b>	2408 bp
<b>RefSeq ORF:</b>	1401 bp
<b>Locus ID:</b>	100820829
<b>UniProt ID:</b>	<a href="#">POCAP1</a>
<b>Cytogenetics:</b>	15q21.3
<b>MW:</b>	54.2 kDa
<b>Gene Summary:</b>	<p>This gene encodes a protein that is abundantly expressed in cardiac tissue. The encoded protein localizes to intercalated discs in cardiomyocytes and functions as an activator of Rho-dependent serum-response factor signaling. Alternative splicing results in multiple transcript variants. Readthrough transcription also exists between this gene and the neighboring downstream gene POLR2M (polymerase (RNA) II (DNA directed) polypeptide M) and is represented with GeneID: 145781. [provided by RefSeq, Mar 2014]</p> <p>Transcript Variant: This variant (1, also known as Gup1) represents the longer transcript and encodes the longer isoform (1).</p>