

Product datasheet for **SC302244**

TM2D2 (NM_001024381) Human Untagged Clone

Product data:


Product Type:	Expression Plasmids
Product Name:	TM2D2 (NM_001024381) Human Untagged Clone
Tag:	Tag Free
Symbol:	TM2D2
Synonyms:	BLP1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC302244 representing NM_001024381. Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGAAGAAATGCTTTTGCCCGTTTGGATTACGTGCATGCAAACAGCGATTTGCAAAGACCGTATGATG
ATGATCATGATCTTACTGGTGAATTACAGACCTGATGAATTTATAGAATGTGAAGACCCAGTGGATCAT
GTTGAAATGCAACTGCATCCCAGGAACCTGGTTATGGTTGTCTCAAGTTCGGCGGTGAGGCTACAGC
GACGTGGAACACACTTCAGTCCAGTGCCATGCCTTAGATGGAATTGAGTGTCCAGTCTTAGGACCTTT
CTACGAGAAAATAAACCTTGATAAAGTATACCGGACACTACTTCATAACCACTTTACTCTACTCCTTC
TTCCTGGGATGTTTGGTGTGGATCGATTCTGTTGGGACACACTGGCACTGCAGTAGGGAAGCTGTTG
ACGCTTGGAGGACTTGGGATTTGGTGGTTTGTGACCTTATTTTGCTAATTACTGGAGGGCTGATGCCA
AGTGATGGCAGCAACTGGTGCCTGTTTACTAA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGCGC
```

Restriction Sites: Sgfl-Mlul



[View online »](#)

Plasmid Map:


ACCN:	NM_001024381
Insert Size:	516 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001024381.1</u>
RefSeq Size:	3334 bp
RefSeq ORF:	516 bp
Locus ID:	83877
UniProt ID:	<u>Q9BX73</u>
Cytogenetics:	8p11.22
Protein Families:	Druggable Genome, Transmembrane
MW:	19 kDa
Gene Summary:	<p>The protein encoded by this gene contains a structural module related to that of the seven transmembrane domain G protein-coupled receptor superfamily. This protein has sequence and structural similarities to the beta-amyloid binding protein (BBP), but, unlike BBP, it does not regulate a response to beta-amyloid peptide. This protein may have regulatory roles in cell death or proliferation signal cascades. This gene has multiple alternatively spliced transcript variants which encode two different isoforms. [provided by RefSeq, Jul 2008]</p> <p>Transcript Variant: This variant (4) differs in the 5' UTR and coding region compared to variant 1. The resulting isoform (b) is shorter and has a distinct N-terminus compared to isoform a. Variants 2, 3, and 4 all encode isoform b.</p>