

Product datasheet for **SC302247**

RTL5 (NM_001024455) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	RTL5 (NM_001024455) Human Untagged Clone
Tag:	Tag Free
Symbol:	RTL5
Synonyms:	6430402L03Rik; MAR5; MART5; RGAG4; SIRH8
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC302247 representing NM_001024455.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGCACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGTCTGAGGCGTCAGGAAATCTCAACAGCCTCCGCATGGCGAATGTGGCCCTGCGCGAAGAATTAAT
GCCCTTCGCGGGGAGAATGCCAATTTGGCCCTCCAGCTCGGGAGAGCCCTGGCCGAGGTTAATTCCTTG
CGGGGCAATGTCTCCAGCTACATCCGCTGGCCAGTACCCATAGTGCCCGTCCTTGCGGAGGAGAATTTG
GAGTTCGCGCTCAGTGAGATCGAAGTCATTCCCGGGGAGAAGTCCCTTCTTGCCGGCCTCCCCCA
CGCGCGGAGCCGACTGTATCTCGGATGATCTTTTGATTAACGTGATCCAGGATCGCAGCACCCCGAC
GGGCCCCGCTGACCCCTCTGCTGCCATCCCGCCCCCGCGGCTGCCCTCGCCCGCTCAAAGGAA
CCGCCCCCGCAACCCCTCTGGCACCCTGGAGCGGCTGAGATAGAGCCCTTTTCAGGAGACCCAGTC
TACCTGGTGAATTCCTGATGCAGCTAGAGACCTTCATAGCCGACCATGAGGTTCAATTCCTCGGGGC
GCCGAGCGGGTGGCTTTCTGATCTCTTTTCTACTGGCGAAGCAAGGACTGGCCATCTCGGTACC
CAGGAAGGAAGCCCTTGCATGCCAATTCCTCGCTTCTGGATGAAATCCGTAAGGAATTCGTGGC
CCCATCCCCCACGTGTGGCTAAAAAGGCCATTGCAAGCTCAAGCAGGGCCACTGTACCTCGGCAGC
TATGCAGATGCTTTTCACTTCTGGCCCAATCTTGCTTGGGATGACTGCCGCCTCCAAAACCAATTC
CTCAAAGGCCTGTGGAATTCCTCCGAAAGAGCTTTATGGTCAACAGAAATGGCCGACCTGGATGAG
CTGATTCTCGAATGTGTGGAGATAGAAAGAAAAGTGCGCGTTCCTCAAGCCTATCCCGCTCCCTGGAGTT
CGAAATATCATCTTCCCTTTTGTCTCCGAGTCCCAATGAGGAAGAAAGTGAAGATGAAGAGTACTACAGT
GAAGATGAAGACCAAGAGGCACGCAGGCACAGGCTTCACTCAAGGACCAGAGGAAGCGCATGAGGGCT
TTTCAGCAAGAGATGAAGGAGAAGGAGGAGGAGGAGATGAAGAAGGAGGAAGAAATGAAGAAGAAGGAA
GAGAAGGAGGAGGAGGAGGAAGAGGAAATGAAACAGAAAGAGGAGGAGGAGGAGATAAGGAACAAGAA
GAGGAGGAAGGAGAGAGTAAGGATGAAGAAGATGAAGATGAGGATGGGGGCCAGAAACCAGAAGGAGAG
CCACAGCAGGATCCAGGACTGAGGAGACCTATGGTGAGGTGGAGGAGGAGGCACTGGATGAGGCCCAA
GATGATGACCTAGATGAGCTGATGGAGATGGAGCCGACCTTTGTTACGCGTCATCCAGACTTCTGGC
CCCACAAGTGTTACCACGCCGAAAACCTTCTGGGTGCATCACCTCCATAATACAGCCTAGCAGACGG
AGGAACCAGAATCGAGTCCCACTTCTGGAAGGCCTTCAGGCACCAATTCACCATCTACAGCTCCCCC
CAACTGATTCGCGCACAGGTCGTTTAGGACAACGCCAAGTTCGAAGACGCCCCCTGTGCTTTCCGC
CTCACTCCGAGACAGGGGGCCACCGAGCTGCTCGTGGCCGAATTCGAGTTTGA
AGCGGACCGACGCTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGAT
ATCCTGGATTACAAGGATGACGACGATAAGGTTTAA
  
```

Restriction Sites: Sgfl-RsrII

ACCN: NM_001024455

Insert Size: 1710 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001024455.3](#)

RefSeq Size: 4354 bp

RefSeq ORF: 1710 bp

Locus ID: 340526

UniProt ID: [Q5HYW3](#)

Cytogenetics: Xq13.1

MW: 64.7 kDa