

Product datasheet for **SC302257**

ZNF772 (NM_001024596) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ZNF772 (NM_001024596) Human Untagged Clone
Tag:	Tag Free
Symbol:	ZNF772
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)


Fully Sequenced ORF: >SC302257 representing NM_001024596.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGCGCGGCTGAGCCGATGGCCCGGCACAGGTTCCCATGAACCTCAGAAGTAATTGTGGACCCTATA
CAGGGCAGGTGAACCTTTGAGGACGTGTTCTGTACTTCTCCAGAGGAGTGGGTGCTCCTTGATGAG
GCTCAGAGGCTCCTGTACCGTGATGTGATGCTGGAGAACCTTGCACTTATGGCCTCTCTGGGACATACA
TCTTTCATGTCCACATAGTTGCCCTCACTGGTGATGGGTAGTGAGCCCTGGGTACCTGACTGGGTGGAC
ATGACCCTAGCTGTTGCCACAGAGACTCCAGGTGGGAGTGACCCTGGTTGTTGGCATGGAATGGAGGAT
GAAGAGATACCTTTTGGCAGAGCTTTTCTATAGGAATGTCACAGATCAGGATCCCAAGGGAGGTCCT
TCTACTCAGAAGGCTTACCCTGTGGGACATGTGGCCTGGTCTTAAAAGACATTTTGCACTTGGCTGAG
CACCAGAAACACACCCAGGGCAGAAACCATACATGTGTGTGCTGTGTGGAAACAGTTCTGCTTCAGT
GCAAACCTTACCAGCACCAGAAGCAGCACAGTGGAGAGAAACCTTTAGAAGTGATAAGAGCAGGCC
TTTCTTCTGAACAACTGTCTGTCAATCAATGGAGATGTCTTTGTGACAGGGGAGGCTTGTAAAGGAC
TTCTAGCCAGCTCAAGCATTTTCGAGCACCATGCCCTCACAATGAGTGGAAAGCCACACAGCAACACC
AAGTGTGAGGAGGCTCTCACTGTGAAAAAGGCATTACAAATGCAGTGAATGTGGGAAAACCTTCAGC
CGCAAAGACTCACTTGTTC AACACCAGAGAGTCCACACTGGAGAAAGGCCTTATGAGTGGCGTGAATGT
GGGAAAACCTTTAGCCGCAAACCCATACTTGTCTCAGCACCAGAGAATCCACACTGGAGAAATGCCTTAT
GAGTGTGGCATAATGTGGGAAAGTTTTTAAATCATAGCTCTAACCTTATTGTACATCAAAGAGTACACACT
GGAGCAAGGCCTTACAAGTGCAGTGAATGTGGGAAAGCCTATAGTCACAAATCTACACTTGTTCAGCAT
GAGAGTATCCATACTGGAGAAAGGCCATATGAGTGCAGCGAATGTGGAAAATACTTTGGTCACAAATAC
AGACTCATTAAACATTGGAGTGTTCATACTGGAGCAAGGCCTTATGAGTGCATCGCATGGGAAAGTTC
TTTAGCCAAAGCTCTGACCTTATTGCACATCAAAGAGTTCATAATGGTGAAGCCCTTATGTGTGCAGT
GAATGTGGAAAAGCCTTTAGCCACAAACATGACTTGTTCAGCATCATAGAATTCACACTGGAGAAAGG
CCATATAAGTGCAGTGAATGTGGGAAAGGCCTTCAGGCAGAGGGCTTCCCTCATCCGCCACTGGAAAATT
CACACTGGAGAAAGGCCTTAG
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGCGC
  
```

Restriction Sites: SgfI-MluI

Plasmid Map:



ACCN: NM_001024596

Insert Size: 1470 bp

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001024596.2
RefSeq Size:	5434 bp
RefSeq ORF:	1470 bp
Locus ID:	400720
UniProt ID:	Q68DY9
Cytogenetics:	19q13.43
MW:	55.4 kDa
Gene Summary:	May be involved in transcriptional regulation.[UniProtKB/Swiss-Prot Function] Transcript Variant: This variant (1) represents the longest transcript and it encodes the longest protein (isoform 1). Sequence Note: The RefSeq transcript and protein were derived from transcript and genomic sequence to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on alignments.