

## Product datasheet for **SC302266**

### **LRRC66 (NM\_001024611) Human Untagged Clone**

#### **Product data:**

|                           |                                            |
|---------------------------|--------------------------------------------|
| Product Type:             | Expression Plasmids                        |
| Product Name:             | LRRC66 (NM_001024611) Human Untagged Clone |
| Tag:                      | Tag Free                                   |
| Symbol:                   | LRRC66                                     |
| Mammalian Cell Selection: | Neomycin                                   |
| Vector:                   | pCMV6-Entry (PS100001)                     |
| E. coli Selection:        | Kanamycin (25 ug/mL)                       |



[View online »](#)

**Fully Sequenced ORF:** >SC302266 representing NM\_001024611.  
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGCATCGCC
ATGAAAAACCTCTATTTTCAGAGTCATTACCATAGTTATAGGTCTTTATTTTACTGGAATAATGACAAAT
GCATCAAGAAAAAGCAATATTTTATTCAATTCTGAATGCCAATGGAATGAATATATTCTGACAAATTTGT
TCTTTTACCGGAAAGTGTGATATACCTGTGGACATATCACAGACAGCAGCCACTGTGGATGTAAGTTTC
AATTTCTTTAGAGTTCTCTTACAGTCTCACACGAAAAAAGAGTGGAAAAATAAACATCTGGACCTC
AGTAACAATCTCATATCAAAAAAACCTTAAGCCCTTTTGCATATTTACATGCTTTGGAAGTGTTAAAC
CTCAGCAACAATGCCATCCACTCCCTCTCATTGGATCTACTCAGTCCTAAGTCTCATGGGTGAAACGC
CACAGAAGCAGCTTCAGAAACAGGTTTCCATTGCTGAAGGTGCTCATTCTTCAAAGAAATAAACTCAGT
GACACTCCCAAGGGACTGTGGAACTGAAGTCATTGCAGAGTTTGGATCTGTCATTCAATGGGATATTG
CAAATAGGGTGGTCTGATTTTCACAACTGCCTGCAACTGGAGAATCTCTGTTTAAAGAGCAACAAGATA
TTCAAATTTCCCCACAAGCCTTCAAGGACCTCAAAAAATTACAGGTCATAGACCTTAGCAACAATGCT
CTGATTACCATCTACCAATGATGATCATAGCTCTAGAATTTCCCATCTAGTGGTTGACTTGGCTGAT
AATAACTGGCAGTGTGATGATAGTGTGGCAGTCTTTCAAATTTTATTTCTGAATCCTGGAGGAAAAAG
TGGAATGTCAATTTGCAACAGGTCTATAGGGAGTGAGGAGGCCAACGGGGCACTCCCAAGAGCAGGATT
TCCAGGGAACCCGCCTTCTCCCATTCATCTGCATCGCATGAAAAGCCTCATAAGGAGCAAGCAGAG
AGGCCCCAGGGAGGAAGGCACACGGGCATTTCTACTCTGGGGAAGAAGGCAAGGCCGGCTCTGGTCTC
AGGAAGAAGCAGAGACGGCTGCCAAGGAGTGTTAGAAGCACCCGCGATGTGCAGGCTGCCGGCAAAAA
GAGGACGCTCCCAGGACCTGGCTCTGGCGGTGTGCTGTGAGTGTTCATCACATTCCTGTGCGCTTC
AGCCTGGGGCTTTTCAAGGCTTATGTTGACAGACTGTGGCAAAAAAGTGCCAGAGCAAAAGCCCT
GGCCTGGACAACCGGTATTCAAACGAGGGCTTCTACGATGACATGGAAGCTGCGGGGCACACACCACAC
CCAGAGACCCATCTGCGCCAAGTATTTCTCATCTAAGCCTCTACGAGAACCAGACCCCTTTCTGGGTG
ACACAGCCACACCCACGCCACCGTAATTCCTGATAGAACTCTGGGAAGGAGCAGAAAGGATCTGGC
AGTTTCGAGAGCCAGGACAGTGCGGGGACAACACCGGGCAGGAAGTGGAATGATGGTGCAGTCTAT
TCCATTCTCCAGAGACATCCACATGCCGTAAACGTGAACATATGTCAGCAGCGCAGGACCACATCCAT
AGGAATGATATTCTCGGAGAATGGACTTATGAACTGTGGCCAGGAAGAGCCTCTCAGTGCACATTCA
GTGGGCGTCTCTCTGTAGCTGGCAGTCTCACGCTGTCTCTGGCTCAAGCCGTTATGATTCCAATGAA
TTAGACCTTCCCTCTCCGAGAAATAACAGCTTCCCTCTGTAATGCTAACACATGCAGAAGCACAG
AGGACTGGAGATAGTAAGGAAAGAGGGGCACTGAACAGTCACTTTGGGACTCGCAGATGGAATTTCT
AAGGAAAGGCAAGTGAGTTCATCCATTGATTTGCTGAGCATACAGCAGCAAGGCTGTCCGGGGCAAGG
GCTGAGGAAGCGCTTTCAGCCCACTACAGCGAGGTTCCATACGGTGACCAAGAGACACAGGCCCATCA
GTCTTCTCTCAAGATGGGACAGTGGCCTGGATGTCACTCCTGCTAACAGGAACCAAGTGCAGAAATCC
ACTCTTCTGACACTTGCTGTGAGTTGGAGAGTGAAGTGAAGTCTGATGAGGGGTCTCTGTTCACTCTG
AGCTCCATAAGTTCAGAGAGTGAAGGAGCAAGACTGAAGAGGCAGTGCCTGATGAGGAGTCCCTGCAG
GACGAGAGCTCAGGGGCAAGCAAGGACAATGTGACGGCTGTAGACAGTCTTGAGGAAAATGTTACCTTC
CAAACAATTCAGGGAATGCAAGAATCAAGAAGATCCCTTTGAAAAACCTCTCATTCTGCTCCAGAC
TCTGGCATGTACAAGACTCATCTGGAAAAATGCCTCTGACACTGATAGATCTGAGGGCCTGTCAACCTGG
CCCAGGTCAACAGGAATAGTCCCTTAGGGGATGAGTTTCCGGGCATGTTCACTTATGATTATGACACA
GCTCTTCAATCAAGGCAGCAGAAATGGCATTGCTCACTTAGAGACTTAGAATTTTCAAATGTGGACGTT
TTACAGCAAAACACCACATGTTCTGCTGAAGTTCCCTCAGATCCTGATAAGGCTGCCTTCCATGAAAGA
GACTCAGACATTTTAAATAA
AGCGGACCGACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGAT
ATCCTGGATTACAAGGATGACGACGATAAGGTTTAA
  
```

**Restriction Sites:** SgfI-RsrII  
**ACCN:** NM\_001024611  
**Insert Size:** 2643 bp

|                               |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
|-------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>OTI Disclaimer:</b>        | Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).                                                                                                                                                                          |
| <b>OTI Annotation:</b>        | This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.                                                                                                                                                                                                                                                      |
| <b>Components:</b>            | The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).                                                                                                                                                                                                                                                                                                                                    |
| <b>Reconstitution Method:</b> | <ol style="list-style-type: none"> <li>1. Centrifuge at 5,000xg for 5min.</li> <li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li> <li>3. Close the tube and incubate for 10 minutes at room temperature.</li> <li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li> <li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li> </ol> |
| <b>RefSeq:</b>                | <u><a href="#">NM_001024611.2</a></u>                                                                                                                                                                                                                                                                                                                                                                                                                                                                                   |
| <b>RefSeq Size:</b>           | 3409 bp                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| <b>RefSeq ORF:</b>            | 2643 bp                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| <b>Locus ID:</b>              | 339977                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
| <b>UniProt ID:</b>            | <u><a href="#">Q68CR7</a></u>                                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |
| <b>Cytogenetics:</b>          | 4q12                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                    |
| <b>Protein Families:</b>      | Transmembrane                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |
| <b>MW:</b>                    | 97.8 kDa                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                |
| <b>Gene Summary:</b>          | This gene encodes a protein belonging to the leucine-rich repeat family of proteins, which are involved in diverse biological processes, including cell adhesion, cellular trafficking, and hormone-receptor interactions. [provided by RefSeq, Apr 2017]                                                                                                                                                                                                                                                               |