

Product datasheet for **SC302317**

SLC2A11 (NM_001024938) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	SLC2A11 (NM_001024938) Human Untagged Clone
Tag:	Tag Free
Symbol:	SLC2A11
Synonyms:	GLUT10; GLUT11
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>NCBI ORF sequence for NM_001024938, the custom clone sequence may differ by one or more nucleotides

```
ATGAGAGCGCTCCGAAGACTGATTCAGGGCAGGATCCTGCTCCTGACCATCTGCGCTGCCGGCATTGGTG
GGACTTTTCAGTTTGGCTATAACCTCTCTATCATCAATGCCCGACCTTGACATTCAGGAATTCACCAA
TGAGACATGGCAGGCGCTACTGGAGAGCCACTGCCCGATCACCTAGTCCTGCTTATGTGGTCCCTCATC
GTGTCTCTGTATCCCCTGGGAGGCTCTTTGGAGCACTGCTTGCAGGTCCTTGGCCATCACGCTGGGAA
GGAAGAAGTCCCTCCTGGTGAATAACATCTTTGTGGTGTGAGCAGCAATCCTGTTGGATTACGCCGAA
AGCAGGCTCCTTTGAGATGATCATGCTGGGAAGACTGCTCGTGGGAGTCAATGCAGGTGTGAGCATGAAC
ATCCAGCCATGTACCTGGGGGAGAGCGCCCTAAGGAGCTCCGAGGAGCTGTGGCCATGAGCTCAGCCA
TCTTTACGGCTCTGGGGATCGTGATGGGACAGGTGGTCGGACTCAGGGAGCTCCTAGGTGGCCCTCAGGC
CTGGCCCTGCTGCTGGCCAGCTGCCTGGTGGCCGGGGCGCTCCAGCTCGCCTCCCTGCCTCTGCTCCCT
GAAAGCCCGCGCTACCTCCTCATTGACTGTGGAGACACCGAGGCTGCCTGGCAGCACTACGGCGGCTCC
GGGGCTCCGGGGACTTGGCAGGGGAGCTGGAGGAGCTGGAGGAGGAGCGCGCTGCCTGCCAGGGCTGCCG
TGCCCGGCGCCATGGGAGCTGTCCAGCATCGGGCCCTGAGGAGACAGGTGACAAGCCTCGTGGTTCTG
GGCAGTGCCATGGAGCTCTGCGGGAATGACTCGGTGTACGCCCTACGCCTCCTCCGTGTTCCGGAAGGCAG
GAGTGCCGGAAGCGAAGATCCAGTACCGCATCATCGGGACTGGGAGCTGCGAGCTGCTCACGGCGGTTGT
TAGTTGTGTGGTAATCGAGAGGGTGGTGGCGCGTGTGCTCATCGGTGGGTACAGCCTGATGACCTGC
TGGGGGAGCATCTTCACTGTGGCCCTGTGCCTGCAGAGCTCCTTCCCCTGGACACTCTACCTGGCCATGG
CCTGCATCTTTGCCTTATCCTCAGCTTTGGCATTGGCCCTGCCGGAGTGACGGGGATCCTGGCCACAGA
GCTGTTTGACCAGATGGCCAGGCTGCTGCCTGCATGGTCTGCGGGGCGCTCATGTGGATCATGCTCATC
CTGGTCGGCCTGGGATTTCCCTTTATCATGGAGGCCTTGTCCCCTTCTCTATGTCCCTTTCTTGGTG
TCTGTGTCTGTGGGGCCATCTACACTGGCTGTTCTTCTGAGACCAAAGGCAAGACCTTCCAAGAGAT
CTCCAAGGAATTACACAGACTCAACTTCCCAGGGGGCCAGGGCCCCACGTGGAGGAGCTGGAGGTT
ATCCAGTCAACAGAACTTAG
```



[View online »](#)

Restriction Sites:	Please inquire
ACCN:	NM_001024938
OTI Disclaimer:	<p>Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.</p> <p>The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info</p>
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001024938.1 , NP_001020109.1
RefSeq Size:	1608 bp
RefSeq ORF:	1491 bp
Locus ID:	66035
UniProt ID:	Q9BYW1
Cytogenetics:	22q11.23
Protein Families:	Transmembrane

Gene Summary:

This gene belongs to a family of proteins that mediate the transport of sugars across the cell membrane. The encoded protein transports glucose and fructose. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Sep 2013]

Transcript Variant: This variant (3, also known as SLC2A11-a) differs in the 5' UTR, lacks a portion of the 5' coding region, and initiates translation at an alternate upstream start codon, compared to variant 1. The encoded isoform (c, also known as GLUT11-S) is shorter and has a distinct N-terminus, compared to isoform a. Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.