

Product datasheet for **SC302344**

ABCF1 (NM_001025091) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ABCF1 (NM_001025091) Human Untagged Clone
Tag:	Tag Free
Symbol:	ABCF1
Synonyms:	ABC27; ABC50
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)




[View online »](#)

Fully Sequenced ORF: >SC302344 representing NM_001025091.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGGATCGCC
ATGCCGAAGGCGCCCAAGCAGCAGCCGCCGAGCCGAGTGGATCGGGGACGGAGAGAGCAGCAGGCCCA
TCAGACAAGTGGTGAAGAAAGGGAAGAAGGACAAGAAGATCAAAAAACGTTCTTTGAAGAGCTGGCA
GTAGAAGATAAACAGGCTGGGGAAGAAGAGAAAGTGTCAAGGAGAAGGAGCAGCAGCAGCAACAG
CAACAGCAGCAAAAAAAGCGAGATACCCGAAAAGGCAGGCGGAAGAAGGATGTGGATGATGATGGA
GAAGAGAAAGAGCTCATGGAGCGTCTTAAGAAGCTCTCAGTGCCAACAGTGATGAGGAGGATGAAGTA
CCCGCCCCAAAACCCCGGAGGGAAGAAAACCAAGGGTGGTAATGTTTTGCAGCCCTGATTCAGGAT
CAGAGTGAGGAAGAGGAGGAGGAAGAAAACATCCTCCTAAGCCTGCCAAGCCGGAGAAGAATCGGATC
AATAAGGCCGTATCTGAGGAACAGCAGCCTGCACTCAAGGGCAAAAAGGGAAGGAGAAGTCAAAA
GGGAAGGCTAAGCCTCAAAATAAATTCGCTGCTCTGGACAATGAAGAGGAGGATAAAGAAGAAGAAAT
ATAAAGGAAAAGGAGCCTCCCAAACAAGGGAAGGAGAAGGCCAAGAAGGCAGAGCAGGGTTCAGAGGAA
GAAGGAGAAGGGGAAGAAGAGGAGGAGGAAGGAGGAGAGTCTAAGGCAGATGATCCCTATGCTCATCTT
AGCAAAAAGGAGAAGAAAAGCTGAAAAACAGATGGAGTATGAGCGCCAAGTGGCTTCATTAAGCA
GCCAATGCAGCTGAAAATGACTTCTCCGTGTCCAGGCGGAGATGCTCCTCCCGCCAAGCCATGTTAGAA
AATGCATCTGACATCAAGCTGGAGAAGTTCAGCATCTCCGCTCATGGCAAGGAGCTGTTGTCATGCA
GACCTGTACATTTAGCCGGCCGCCGCTACGGGCTGGTAGGCCAATGGCAAGGCAAGACCACACTC
CTCAAGCACATTCGCAACCGAGCCCTGAGCATCCCTCCCAACATTGATGTGTTGCTGTGTGAGCAGGAG
GTGGTAGCAGATGAGACACCAGCAGTCCAGGCTGTTCTTCGAGCTGACACCAAGCGATTGAAGCTGCTG
GAAGAGGAGCGGGCTTCAGGGACAGCTGGAACAAGGGGATGACACAGCTGCTGAGAGGCTAGAGAA
GTGTATGAGGAATTGCGGGCCACTGGGGCGGACAGCTGCAGAGGCCAAAGCACGGCGGATCCTGGCTGGC
CTGGGCTTTGACCCTGAAATGCAGAATCGACCCACACAGAAGTTCTCAGGGGGCTGGCGCATGCGTGTC
TCCCTGGCCAGGGCACTGTTTATGGAGCCCACTGCTGATGCTGGATGAGCCACCAACCACCTGGAC
CTCAACGCTGTACATCTGGCTTAATAACTACCTCCAGGGCTGGCGGAAGACCTTGCTGATCGTCTCCCAT
GACCAGGGCTTCTGGATGATGTCTGCACTGATATCATCCACCTCGATGCCAGCGGCTCCACTACTAT
AGGGGCAATTACATGACCTTCAAAAAGATGTACCAGCAGAAGCAGAAAAGAACTGCTGAAACAGTATGAG
AAGCAAGAGAAAAAGCTGAAGGAGCTGAAGGCAGGCGGGAAGTCCACCAAGCAGGCGAAAAACAAACG
AAGGAAGCCCTGACTCGAAGCAGCAGAAATGCCGACGGAACCAAGATGAGGAATCCAGGAGGCC
CCTGAGCTCCTGAAGCGCCCTAAGGAGTACACTGTGCGCTTCACTTTTCCAGACCCCCACCACTCAGC
CCTCCAGTGTGGTCTGCATGGTGTGACATTCGGCTACCAGGGACAGAAACCACTCTTTAAGAAGTGG
GATTTTGGCATCGACATGGATTCAAGGATTTGCATTGTGGGCCCTAATGGTGTGGGGAAGAGTACGCTA
CTCCTGTCTGACTGGCAAGCTGACACCGACCCATGGGGAAATGAGAAAAGAACCCGGCTGAAAATT
GGCTTCTTCAACCAGCAGTATGCAGAGCAGCTGCGCATGGAGGAGACGCCACTGAGTACCTGCAGCGG
GGCTTCAACCTGCCCTACCAGGATGCCCGCAAGTGCCTGGGCCGCTTCGGCTGGAGAGTACGCCCCAC
ACCATCCAGATCTGCAAACTCTCTGGTGGTCAAGAGCGCGAGTTGTGTTGCTGAGCTGGCCTGTCCG
GAACCTGATGTCCTCATCTTGGACGAGCCAACCAATAACCTGGACATAGAGTCTATTGATGCTCAGGG
GAGGCCATCAATGAATACAAGGGTGTGTGATCGTTGTCAGCCATGATGCCCGACTCATCACAGAAACC
AATTGCCAGTGTGGTGGTGGAGGAGCAGAGTGTAGCCAAATCGATGGTACTTTGAAGACTACAAG
CGGGAGGTGTTGGAGGCCCTGGGTGAAGTCATGGTCAAGCCGCCCCGAGAGTGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: SgfI-MluI

Plasmid Map:


ACCN: NM_001025091

Insert Size: 2538 bp

OTI Disclaimer: Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq:	<u>NM_001025091.1</u>
RefSeq Size:	3474 bp
RefSeq ORF:	2538 bp
Locus ID:	23
UniProt ID:	<u>Q8NE71</u>
Cytogenetics:	6p21.33
Protein Families:	Druggable Genome
MW:	95.9 kDa
Gene Summary:	<p>The protein encoded by this gene is a member of the superfamily of ATP-binding cassette (ABC) transporters. ABC proteins transport various molecules across extra- and intra-cellular membranes. ABC genes are divided into seven distinct subfamilies (ABC1, MDR/TAP, MRP, ALD, OABP, GCN20, White). This protein is a member of the GCN20 subfamily. Unlike other members of the superfamily, this protein lacks the transmembrane domains which are characteristic of most ABC transporters. This protein may be regulated by tumor necrosis factor-alpha and play a role in enhancement of protein synthesis and the inflammation process. [provided by RefSeq, Jul 2008]</p> <p>Transcript Variant: This variant (1) represents the longer transcript and encodes the longer isoform (a).</p>