

Product datasheet for **SC302398**

VEGFA (NM_001025370) Human Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: VEGFA (NM_001025370) Human Untagged Clone
Tag: Tag Free
Symbol: VEGFA
Synonyms: MVCD1; VEGF; VPF
Mammalian Cell Selection: None
Vector: pCMV6-XL4
E. coli Selection: Ampicillin (100 ug/mL)

Fully Sequenced ORF: >OriGene ORF sequence for NM_001025370 edited
 GTCGTGGCGCTGGGGCTAGCACCAGCGCTCTGTCGGGAGGCGCAGCGGTAGGTGGACC
 GGTACGGGACTCACCGCCAGGGCGCTCGGTGCTGGAATTTGATATTCATTGATCCGGG
 TTTTATCCCTCTTCTTTTTCTTAAACATTTTTTTTTTAAACTGTATTGTTCTCGTTTT
 AATTTATTTTTGCTTGCCATTCCCCACTTGAATCGGGCCGACGGCTTGGGAGATTGCTC
 TACTTCCCAAATCACTGTGGATTTTGGAAACCAGCAGAAAGAGGAAAGAGGTAGCAAGA
 GCTCCAGAGAGAAGTCGAGGAAGAGAGACGGGGTCAGAGAGAGCGCGGGCGTGC
 GCAGCGAAAGGGACAGGGCAAAGTGAGTGACCTGCTTTTGGGGTGACCGCGGAGCGC
 GGCGTGAGCCCTCCCCCTGGGATCCCGCAGCTGACCAGTCGCGCTGACGGACAGACAGA
 CAGACACCGCCCCAGCCCCAGTACCACCTCCTCCCCGGCCGGCGGACAGTGGACG
 CGCGGGCAGCCGCGGGCAGGGCCGGAGCCCGCCGGAGGCGGGGTGGAGGGGGTCCG
 GGCTCGCGGCGTCGCACTGAAACTTTTCGTCCAACCTTCTGGGCTGTTCTCGCTTCGAG
 GAGCCGTGGTCCGCGCGGGGAAGCCGAGCCGAGCGGAGCCGCGAGAAGTCTAGCTCGG
 GCCGGGAGGAGCCGAGCCGGAGGAGGGGAGGAGGAAGAAGAGAAGGAAGAGGAGAGGG
 GGCCGAGTGGGACTCGGCGCTCGGAAGCCGGGCTCATGGACGGGTGAGGCGCGGTGT
 GCGCAGACAGTGTCCAGCCGCGCGCTCCCCAGGCCCTGGCCCGGGCTCGGGCCGGG
 GAGGAAGAGTAGCTCGCCGAGGCGCCGAGGAGAGCGGGCCGCCACAGCCCGAGCCGGA
 GAGGGAGCGGAGCCGCGCCGGCCCGGTCGGGCTCCGAAACATGAACCTTCTGCTGT
 CTTGGGTGCATTGGACCTTGCCTTGTGCTCTACCTCCACCATGCCAAGTGGTCCCAGG
 CTGCACCCATGGCAGAAGGAGGAGGGCAGAATCATCACGAAGTGGTGAAGTTCATGGATG
 TCTATCAGCGCAGCTACTGCCATCCAATCGAGACCCTGGTGGACATCTTCCAGGAGTACC
 CTGATGAGATCGAGTACATCTTCAAGCCATCCTGTGTGCCCTGATGCGATGCGGGGGCT
 GCTGCAATGACGAGGGCCTGGAGTGTGTGCCACTGAGGAGTCCAACATCACCATGCAGA
 TTATGCGGATCAAACCTCAACAAGGCCAGCACATAGGAGAGATGAGCTTCTACAGCACA
 ACAATGTGAATGCAGACCAAAGAAAGATAGAGCAAGACAAGAAAAATGTGACAAGCCGA
 GCGGTGA



[View online »](#)

3' Read Nucleotide Sequence:	>OriGene 3' read for NM_001025370 unedited NNNGGGCCATGGGGATGGCAACTTCCAGGGCCAGNANGAGCACTGGGGNAGGGTCACAGG GATGCCACCCGGGATCTGTTTCAGGAAAAGCTATGACCGCGCCGCAATCTAGCTACAGGA ATCCCAGAAATAAACTCTCTAATCTTCCGGGCTCGGTGATTTAGCAGCAAGAAAAATAA AATGGCGAATCCAATCCAAGAGGGACCGTGTGGGTACCCGCCGGGAATGCTTCCGC CGGAGTCTCGCCCTCCGGACCCAAAGTGTCTGCGCAGAGTCTCCTCTTCCTTCATTTC GGTTTCTGGATTAAGGACTGTTCTGTGATGGTGTGGTGGTGGCGGCAGCGTGGTTT CTGTATCGATCGTTCGTATCAGTCTTCTGGTGAGAGATCTGGTTCCCGAAACCCTGA GGGAGGCTCCTTCTCTGCCCCGCTCACCGCCTCGGCTTGTCACATTTTCTTGTCTTG CTCTATCTTTCTTGGTCTGCATTACATTTGTTGTGCTGTAGGAAGCTCATCTCTCCTA TGTGCTGGCCTTGGTGGGTTTGTATCCGCATAATCTGCATGGTGTGTTGGACTCCTCAG TGGGCACACACTCCAGGCCCTCGTATTGCAGCAGCCCCGCATCGCATCAGGGGCACAC AGGATGGCTTGAAGATGTACTCGATCTCATCAGGGTACTCTGGGAGATGTCCACCAAGG TCTCGATTGGATGGCAGTAGCTGCGCTGGTAGACATCCATGGACCTTCACCACTTCGTGA AGGATCTGCCCCCTCCTTCTGCCATGGGGTGCAGCCTGGGGCCACTTGGCATGGGGGG AGGTAGAACAACCAGGAAGGGGTCCAAGCCACCAAGAA
Restriction Sites:	Please inquire
ACCN:	NM_001025370
Insert Size:	2100 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	The ORF of this clone is found to be a perfect match to NM_001025370.1.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001025370.1 , NP_001020541.2
RefSeq Size:	3410 bp
RefSeq ORF:	984 bp
Locus ID:	7422
UniProt ID:	P15692
Cytogenetics:	6p21.1
Protein Families:	Druggable Genome, Secreted Protein

Protein Pathways:	Bladder cancer, Cytokine-cytokine receptor interaction, Focal adhesion, mTOR signaling pathway, Pancreatic cancer, Pathways in cancer, Renal cell carcinoma, VEGF signaling pathway
Gene Summary:	<p>This gene is a member of the PDGF/VEGF growth factor family. It encodes a heparin-binding protein, which exists as a disulfide-linked homodimer. This growth factor induces proliferation and migration of vascular endothelial cells, and is essential for both physiological and pathological angiogenesis. Disruption of this gene in mice resulted in abnormal embryonic blood vessel formation. This gene is upregulated in many known tumors and its expression is correlated with tumor stage and progression. Elevated levels of this protein are found in patients with POEMS syndrome, also known as Crow-Fukase syndrome. Allelic variants of this gene have been associated with microvascular complications of diabetes 1 (MVCD1) and atherosclerosis. Alternatively spliced transcript variants encoding different isoforms have been described. There is also evidence for alternative translation initiation from upstream non-AUG (CUG) codons resulting in additional isoforms. A recent study showed that a C-terminally extended isoform is produced by use of an alternative in-frame translation termination codon via a stop codon readthrough mechanism, and that this isoform is antiangiogenic. Expression of some isoforms derived from the AUG start codon is regulated by a small upstream open reading frame, which is located within an internal ribosome entry site. The levels of VEGF are increased during infection with severe acute respiratory syndrome coronavirus 2 (SARS-CoV-2), thus promoting inflammation by facilitating recruitment of inflammatory cells, and by increasing the level of angiotensin II (Ang II), one of two products of the SARS-CoV-2 binding target, angiotensin-converting enzyme 2 (ACE2). In turn, Ang II facilitates the elevation of VEGF, thus forming a vicious cycle in the release of inflammatory cytokines. [provided by RefSeq, Jun 2020]</p> <p>Transcript Variant: This variant (6) lacks two alternate in-frame exons in the 3' coding region, compared to variant 1. This variant can initiate translation from four non-AUG (CUG) sites, and also from a downstream, in-frame AUG. The isoform (f) represented in this RefSeq is derived from the 5'-most CUG start codon, and is shorter than isoform a. CCDS Note: This CCDS representation starts at an upstream non-AUG (CUG) translation start codon that is supported by conservation data and experimental evidence, including that in PMIDs:11352659, 11563986 and 11731620. Most of the publicly available long mRNAs at this locus do not extend far enough 5' to include the CUG start codon, including BC065522.1 which supports this variant. However, BC058855.1, several unspliced partial transcripts and the aforementioned publications all support transcription initiation further upstream. Hence, this CCDS representation includes the extended 5' CDS. Alternative translation initiation from a downstream AUG codon (Met-181), which results in a shorter isoform with a predicted signal peptide, has also been shown for this gene. The shorter isoform encoded by this variant is represented by CCDS 55014.1.</p>