

Product datasheet for **SC302422**

ARFIP1 (NM_001025593) Human Untagged Clone

Product data:


Product Type:	Expression Plasmids
Product Name:	ARFIP1 (NM_001025593) Human Untagged Clone
Tag:	Tag Free
Symbol:	ARFIP1
Synonyms:	HSU52521
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC302422 representing NM_001025593. Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTGTAGTAACCGTCAGAATTTTGTAAACGACTACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGCTCAAGAATCTCCCAAAAATTCAGCAGCAGAAAATCCAGTGACTAGTAATGGAGAAGTTGATGAC
TCTCGTGAACATAGCTTTAATAGGGATTGAAGCATTACCATCTGGACTTGGTCTCTCAGAAAACC
CAAATTACATCTCATGGCTTTGACAATACCAAAGAGGGTGTATTGAAGCAGGAGCATTTCAGGTGGC
CAGAGAACACAGACAAAAAGTGGACCAGTTATTCTAGCAGATGAAATAAAAATCCTGCAATGGAAAAG
TTAGAACTTGTTAGAAAATGGAGTCTAAACACCTATAAGTGTACTCGACAGATTATCTGAGAAGCTA
GGCCGTGGCTCAAGAAGTGGACCTTGAAGTGAAGCTCAGATTGATATATTAAGGATAACAAGAAA
AAATATGAAAATATTTTAAAAGTGGCTCAAACATTGTCGACCCAGCTTTCCAGATGGTACATACCCAA
AGGCAACTTGAGATGCATTTGCTGACCTGAGTTTGAAGTCACTAGAAGTTCATGAAGAATTTGGCTAT
AATGCCGATACCCAGAACTGCTGGCTAAAAATGGAGAGACTCTTCTTGGGGCCATTAATTTTTTTCATT
GCTAGTGTGAACACTTTGGTGAATAAAACCATTGAAGATACATTAATGACTGTGAAACAGTATGAAAGT
GCCAGGATTGAATATGATGCATATCGCACTGATTTGGAAGAACTGAATCTTGGACCACGTGACGCAAAC
ACTCTGCCAAAGATTGAGCAGTCACAGCATCTCTCAAGCACATAAGGAAAAATATGATAAAATGCCG
AATGATGTTTCTGTCAAATGAAATTTCTAGAAGAAAAAAGGTTAAAGTATTGCACAATCAGCTGGTC
CTTTTCCACAATGCCATTGCCGTTACTTTGCTGGGAATCAGAAGCAGCTTGAACAGACACTTAAACAG
TTCCATATCAAATGAAAACCCCTGGAGTGGATGCCCATCTTGGCTTGAAGAACAGTAA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites: Sgfl-Mlul



[View online »](#)

Plasmid Map:


ACCN: NM_001025593

Insert Size: 1026 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001025593.2](#)

RefSeq Size: 3031 bp

RefSeq ORF: 1026 bp

Locus ID: 27236

UniProt ID: [P53367](#)

Cytogenetics: 4q31.3

MW: 38.6 kDa

Gene Summary: Putative target protein of ADP-ribosylation factor.[UniProtKB/Swiss-Prot Function]
Transcript Variant: This variant (2) lacks an in-frame exon in the central coding region, compared to variant 1. The encoded isoform (2) is shorter, compared to isoform 1. Variants 2, 3, and 6 encode the same isoform.