

Product datasheet for **SC302428**

Rho GTPase activating protein 24 (ARHGAP24) (NM_001025616) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Rho GTPase activating protein 24 (ARHGAP24) (NM_001025616) Human Untagged Clone
Tag:	Tag Free
Symbol:	ARHGAP24
Synonyms:	FILGAP; p73; p73RhoGAP; RC-GAP72; RCGAP72
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)




[View online »](#)

Fully Sequenced ORF: >SC302428 representing NM_001025616.
Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGCATCGCC
ATGGAGGAGAACAATGACTCCACGGAGAACCCCAACAAGGCCAAGGGCGGCAGAATGCCATCAAGTGT
GGGTGGGTGAGGAAGCAAGGAGGCTTTGTCAAGACTTGGCATACTCGCTGGTTTTGTCTCAAGGGGGAT
CAGCTCTATTATTTCAAAGATGAAGATGAAACCAAGCCCTTGGTACTATTTTTCTGCCTGGAAATAAA
GTTTCTGAGCATCCCTGCAATGAAGAGAACCAGGGAAGTTCTTTTTGAAGTAGTTCAGGAGGCGAT
CGAGATCGGATGACAGCAATCATGAAAGCTACCTCCTCATGGCAAGCACCCAGAATGATATGGAAGAC
TGGGTGAAGTCAATCCGCCGAGTCATATGGGGACCTTTCCGGAGGAGGCATTTTTGGACAGAACTGGAG
GATACTGTTCTGTTATGAGAAGAGATATGGGAACCGTCTGGCTCCGATGTTGGTGGAGCAGTGCCTGGAC
TTTATCCGACAAAGGGGGCTGAAAGAAGAGGGTCTTTTCGACTGCCAGGCCAGGCTAATCTTGTTAAG
GAGCTCCAAGATGCCTTTGACTGTGGGGAGAAGCCATCATTTGACAGCAACACAGATGTACACACGGTG
GCATCACTTCTTAAGCTGTACCTCCGAGAACTTCCAGAACCAGTTATTCCTTATGCGAAGTATGAAGAT
TTTTTGTGATGTGCCAACTGCTCAGCAAGGAAGAGGAAGCAGGTGTTAAGGAATTAGCAAAGCAGGTG
AAGAGTTTGCCAGTGGTAAATTACAACCTCCTCAAGTATATTTGCAGATTCTTGATGAAGTACAGTCC
TACTCGGGAGTTAACAAAATGAGTGTGCAGAACTTGGCAACGGTCTTTGGTCTTAATATCCTGCGCCCC
AAAGTGGAAAGTCTTTGACTATCATGGAGGGCACTGTGGTGGTCCAGCAGTTGATGTCAGTGTGATT
AGCAAACATGATTGCCTTTCCCAAAGATGCAGAACTACAAAGCAAGCCCCAAGATGGAGTGAGCAAC
AACAAATGAAATTCAGAAGAAAGCCACCATGGGGCAGTTACAGAACAAGGAGAACAATAACACCAAGGAC
AGCCCTAGTAGGCAGTGTCTCTGGGACAAGTCTGAGTCACCCAGAGAAGCAGCATGAACAATGGATCC
CCCACAGCTCTATCAGGCAGCAAAACCAACAGCCCAAAGAACAGTGTTCACAAGCTAGATGTGCTAGA
AGCCCCCTCTCATGGTCAAAAAGAACCAGCCTTTAATAGGGTAGTGGGATAGTTACCAATGGGTCC
TTCAGCAGCAGTAAATGCAGAAGGTCTTGAGAAAACCCAAACCACCCCAATGGGAGCCTACAGGCCAGA
AGGAGCTCTTCACTGAAGGTATCTGGTACCAAAATGGGCACGCACAGTGTACAGAATGGAACGGTGGCG
ATGGGCATTTTGAACAGCGACACACTCGGGAACCCACAATGTTGAAACATGAGCTGGCTGCCAAT
GGCTATGTGACCCTGAGGGATAACAAGCAGAAAGAACAAGCTGGAGAGTTAGGCCAGCACAAACAGACTG
TCCACCTATGATAATGTCCATCAACAGTCTCCATGATGAACCTTGATGACAAGCAGAGCATTGACAGT
GCTACCTGGTCCACTTCTCTGTGAAATCTCCCTCCCTGAGAACTCCAACCTCTGTGCTCTTCTACC
ACCACCTGCCAGAGCAAGACTTTTTTGGGGGAACTTTGAGGACCCTGTTTTGGATGGGCCCCCGCAG
GACGACCTTTCCACCCAGGGACTATGAAAGCAAAAGTGACCACAGGAGTGTGGGAGGTGGAAGTAGT
CGTGCCACCAGTAGCAGTGACAACAGTGAGACATTTGTGGGCAACAGCAGCAGCAACCACAGTGCCTG
CACAGTTTAGTTTCCAGCCTGAAACAGGAAATGACCAAACAGAAGATAGAGTATGAGTCCAGGATAAAG
AGCTTAGAACAGCGAAACTTGACTTTGGAACAGAAATGATGAGCCTCCATGATGAACTGGATCAGGAG
AGGAAAAAGTTTACAATGATAGAAATAAAATGCGAAATGCCGAGCGAGCAAAAGAAGATGCCGAGAAA
AGAAATGACATGCTACAGAAAGAAATGGAGCAGTTTTTTTTCCACGTTTGGAGAAGTACAGTGGAAACC
AGGAGAACCAGAGAGGAAACACAATATGGATTCACTAGTA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGCGC
```

Restriction Sites: SgfI-MluI

Plasmid Map:



ACCN: NM_001025616

Insert Size: 2247 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001025616.2](#)

RefSeq Size: 4644 bp

RefSeq ORF: 2247 bp

Locus ID: 83478

UniProt ID: [Q8N264](#)

Cytogenetics: 4q21.23-q21.3

MW: 84.3 kDa

Gene Summary: This gene encodes a Rho-GTPase activating protein, which is specific for the small GTPase family member Rac. Binding of the encoded protein by filamin A targets it to sites of membrane protrusion, where it antagonizes Rac. This results in suppression of lamellae formation and promotion of retraction to regulate cell polarity. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Sep 2016]
Transcript Variant: This variant (1) represents the longest transcript and encodes the longest isoform (1).