

## Product datasheet for SC302541

### RBFOX2 (NM\_001031695) Human Untagged Clone

#### Product data:


Product Type:	Expression Plasmids
Product Name:	RBFOX2 (NM_001031695) Human Untagged Clone
Tag:	Tag Free
Symbol:	RBFOX2
Synonyms:	dj106l20.3; Fox-2; FOX2; fxb; HNRBP2; HRNBP2; RBM9; RTA
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC302541 representing NM_001031695. Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTGTAGTGAACCGTCAGAATTTGTAAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGAGAAAAAGAAAATGGTAACTCAGGGTAACCGAGGAGCCGACAACAACCTCTGACGCAATGGTTCAG
CCTTTTACTACCATCCCATTTCACCACCTCCGCAGAATGGAATCCCACAGAGTATGGGGTGCCACAC
ACTCAAGACTATGCCGGCCAGACCGGTGAGCATAACCTGACACTCTACGGAAGTACGCAAGCCCACGGG
GAGCAGAGCAGCAACTCACCCAGCACACAAAATGGATCTTTCAGACAGAAGGTGGAGCACAGACAGAC
GGCCAGCAGTACAGACACAAAGTAGTAAAATTCAGAGAGTAAATCTACCCGAAACGGCTGCATGTC
TCTAATATTCCTTCCGCTTCCGGGACCCTGACCTCCGGCAGATGTTTGGGCAGTTTGGCAAATCCTA
GATGTAGAAATAATCTTTAATGAACGTGGCTCTAAGGGATTCCGGTTCGTAACCTTCGAGAATAGTGCT
GATGCAGACAGGGCCAGGGAGAAATACACGGCACCGTGGTAGAGGGCCGTAAAATCGAGGTGAATAAT
GCTACAGCACGTGTAATGACCAATAAGAAGATGGTACACCATATGCAATGGTTGGAAATTAAGCCCA
GTAGTTGGAGCTGTATATGGTCCGGAGTTATATGCAGCATCCAGCTTTCAGCAGATGTGTCCTTAGGC
AATGATGCAGCAGTGCCCTATCAGGAAGAGGGGTATCAACACTTACATTCCTTAATCAGTCTCCCT
TTAGTTCTGGCTTCCCTTACCCTACTGCAGCCACCACGGCAGCCGCTTTCAGAGGAGCCATTTGAGG
GGCAGAGGGCGGACAGTATATGGTGCAGTCCGAGCGGTACCTCCAACAGCCATCCCCGCCTATCCAGT
GTGGTTTACCAGGACGGATTTTACGGTGTGCTGACCTCTATGGTGGATATGCAGCTACAGATATGCACAG
CCTGCTACTGCAACCGCAGCCAGCTGCTGCAGCCGCTGCAGCCGCTTACAGTACAGGTTATGGCAGG
GTGTACACAGCCGACCCCTACCATGCCCTTGCCCTGCCCTAGCTATGGAGTTGGCGCTGTGGCGAGT
TTATACCGAGGTGGCTACAGCCGATTTGCCCTACTGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites: SgfI-MluI



[View online »](#)

**Plasmid Map:**


**ACCN:** NM\_001031695

**Insert Size:** 1143 bp

**OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

**OTI Annotation:** This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_001031695.2](#)

**RefSeq Size:** 6997 bp

**RefSeq ORF:** 1143 bp

**Locus ID:** 23543

**UniProt ID:** [O43251](#)

**Cytogenetics:** 22q12.3  
**Protein Families:** Druggable Genome, Transcription Factors  
**MW:** 40.4 kDa

**Gene Summary:** This gene is one of several human genes similar to the *C. elegans* gene Fox-1. This gene encodes an RNA binding protein that is thought to be a key regulator of alternative exon splicing in the nervous system and other cell types. The protein binds to a conserved UGCAUG element found downstream of many alternatively spliced exons and promotes inclusion of the alternative exon in mature transcripts. The protein also interacts with the estrogen receptor 1 transcription factor and regulates estrogen receptor 1 transcriptional activity. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jul 2008]

Transcript Variant: This variant (1) encodes isoform 1. Sequence Note: This RefSeq record was created from transcript and genomic sequence data because no single transcript was available for the full length of the gene. The extent of this transcript is supported by transcript alignments.