

Product datasheet for **SC302566**

ZNF613 (NM_001031721) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ZNF613 (NM_001031721) Human Untagged Clone
Tag:	Tag Free
Symbol:	ZNF613
Vector:	<u>pCMV6 series</u>
Fully Sequenced ORF:	>NCBI ORF sequence for NM_001031721, the custom clone sequence may differ by one or more nucleotides

```

ATGATCAAGTCCCAGGAATCACTGACCCTGGAGGATGTGGCTGTGGAGTTCACTTGGGAG
GAGTGGCAGCTCCTCGGCCCTGCTCAGAAGGACCTGTACCGAGACGTGATGTTGGAGAAC
TATAGCAACCTCGTGTCACTGAGGATCAAGCCAGCAAACAGATGCACTCTTCAAGTTG
GAACAAGGAGAGCCATGGACAGTAGAAAATGAAATCCACAGCCAAATCTGTCCAGAAATC
AAGAAAGTTGACAATCATCTACAGATGCACTCACAAAAGCAAAGATGTCTGAAGAGAGTG
GAACAATGCCATAAACATAATGCATTTGGAAACATCATTATCAGAGGAAAAGTGATTTT
CCTTTAAGGCAAAATCATGATACATTTGACTTACATGGGAAAATACTGAAATCAAATTTA
AGTTTAGTCAACCAGAACAAAAGGTATGAAATCAAGAATTCTGTGGGGGTTAATGGAGAT
GGGAAATCCTTCTTCATGCCAAGCATGAACAATTTTATAATGAAATGAACCTTCCCCGAA
GGTGGAAATTCGTGAATACAAATTCACAATTCATTAAGCATCAGCGAACTCAAAACATA
GATAAACCCCATGTATGCACTGAGTGTGGGAAGGCTTCTCAAGAAAGTCTCGCCTCATC
TATCATCAGAGAGTTCACACTGGGGAGAAACCTCATGGATGCAGTATATGTGGGAAAGCC
TTCTCCAGAAAGTCCGGGCTCACTGAACACCAGAGAAACACACAGGAGAGAAACCTAT
GAATGCACTGAATGTGACAAAGCATTCCGCTGGAAATCACAGCTCAATGCACACCAGAAA
ATTCATACAGGAGAGAAGTCATATATATGCAGTGATTGTGGAAAAGGCTTCATCAAGAAG
TCTCGGCTCATTAAATCATCAGAGAGTTCATACAGGAGAGAAACCACATGGATGCAGCCTG
TGTGGGAAGGCCTTCTCCAAAAGGTCCAGGCTCACTGAACACCAGAGAACTCATACAGGA
GAGAAACCTATGAATGCACTGAATGTGACAAAGCATTCCGCTGGAAATCACAGCTCAAT
GCACATCAGAAAGCTCACACAGGAGAGAAAGTCATATATATGCCGTGATTGTGGAAAAGGC
TTCATTCAGAAGGGAATCTCATTGTACATCAGCGAATTCATACTGGAGAAAAACCTAT
ATATGCAATGAATGTGGAAAAGGCTTCATCCAAAAGGGCAACCTCCTTATTCATCGACGT
ACTCACACTGGAGAGAAACCTATGTATGCAATGAATGTGGAAAAGGCTTCAGCCAGAAG
ACATGTTTAATATCCCATCAGAGATTTACACAGGAAAGACACCTTTGTATGTACTGAG
TGTGGAAAATCCTGCTCACACAAGTCAGGTCTCATTAAACCACAGAGAAATTCACACAGGA
GAGAAACCTATACATGCAGTGAATGTGGAAAAGCTTTCAGAGATAAATCATGTCTCAAC
AGACATCGGAGAACTCATACAGGGGAGAGACCGTATGGATGCTCTGATTGTGGAAAAGCT
TTCTCCCACTTGTGATGCCTTGTATATCATAAGGGAATGCTGCATGCAAGAGAGAAATGT
GTAGGTTCACTCAAATTGGAAAATCCTTGTCTCAGAGAGTCATAGCTTATCACATACACGT
GATCTCATACAGGATAAAGACTCTGTTAACATGGTGACTCTGCAGATGCCTTCTGTGGCA
GCTCAGACCTCATTAACTAACAGTGCCTTCCAAGCAGAGAGCAAAGTAGCCATTGTGAGC
CAGCCTGTTGCCAGAAGTTCAGTCTCAGCAGATAGTAGAATTTGCACAGAATAA

```



[View online »](#)

Restriction Sites:	Please inquire
ACCN:	NM_001031721
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001031721.1</u> , <u>NP_001026891.1</u>
RefSeq Size:	2308 bp
RefSeq ORF:	1854 bp
Locus ID:	79898
UniProt ID:	<u>Q6PF04</u>
Cytogenetics:	19q13.41
Protein Families:	Transcription Factors
Gene Summary:	May be involved in transcriptional regulation.[UniProtKB/Swiss-Prot Function] Transcript Variant: This variant (1) represents the longer transcript and encodes the longer isoform (1).