

## Product datasheet for **SC302603**

### CSNK1G3 (NM\_001031812) Human Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	CSNK1G3 (NM_001031812) Human Untagged Clone
Tag:	Tag Free
Symbol:	CSNK1G3
Synonyms:	CKI-gamma 3; CSNK1G3L
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC302603 representing NM_001031812. Blue=Insert sequence Red=Cloning site Green=Tag(s)


```
GCTCGTTTGTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGAAAATAAAAAGAAAGACAAGGACAAATCAGATGATAGAATGGCAGCAGCTAGTGGTCGATCGGGA
CACAACTCGAGGAACTGGGTCTTCATCGTCTGGAGTTTTAATGGTTGGACCTAACTTTAGAGTTGGA
AAAAAATTTGGATGTGGCAATTTGGAGAATTACGATTAGGGAAAAATTTATACACAAATGAATATGTG
GCAATTAAGTTGGAGCCCATGAAATCAAGAGCACCACAGCTACATTTGGAATACAGATTCTATAAGCAG
TTAGGATCTGGAGATGGTATACCTCAAGTTTACTATTTTCGGCCCTTGTGGTAAATACAATGCTATGGTG
CTGGAAGTCTGGGACCTAGTTTGAAGACTTGTGGTACTTGTGTGACAGAACATTTTCTCTAAAACA
GTTCTCATGATAGCTATACAAGTATTTCTCGCATGGAATATGTCCATCAAGAAGTATGATATACAGA
GATGTAACCTGAGAACTTCTAATAGGACGACCAGGAAACAAAACCCAGCAAGTTATTCACATTATA
GATTTTGGTTTGGCAAAGGAATATATTGATCCGGAGACAAAGAAACACATACCATACAGAGAACAAG
AGCCTTACAGGAACAGCTAGATATATGAGCATAAACACACATTTAGGAAAAGAACAAAGTAGAAGAGAC
GATTTAGAAGCTTTAGGTCATATGTTTCTGAGAGGCAGTCTTCTTGGCAAGGCTTAAAG
GCTGACACATTAAGGAGAGGTATCAGAAAATGGAGATACAAAACGGGCTACACCAATAGAAGTGTTA
TGTGAAAATTTCCAGAAATGGCAACATATCTTCGTTATGTAAGAAGGCTAGATTTTTTTGAAAACCA
GACTATGACTACTTAAGAAAGCTTTTTACTGACTTGTGGTATCGAAAAGGATATGTTTGTATTGAA
TATGACTGGATTGGTAAACAGTTGCCTACTCCAGTGGGTGAGTTCAGCAAGATCCCTGCTGTGCATCA
AACAGAGAAGCACATCAACACAGAGATAAGATGCAACAATCCAAAACCCAGGTTGTAAGTTCTACAAT
GGAGAGTTAAACACAGATGACCCACCAGGACGTTCAAATGCACCCATCACAGCCCCTACTGAAGTA
GAAGTGATGGATGAAACCAACTGCCAGAAAGTGTGAACATGTGGTGCTGCTGTTTTTCAAACGAAG
AAAAGGAAAACCATACAGCGCCACAAA TGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```



[View online »](#)

Restriction Sites: SgfI-MluI

Plasmid Map:



ACCN: NM\_001031812

Insert Size: 1272 bp

**OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

**OTI Annotation:** This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM\\_001031812.3](#)

RefSeq Size: 4651 bp

RefSeq ORF: 1272 bp

Locus ID: 1456

<b>UniProt ID:</b>	<u>Q9Y6M4</u>
<b>Cytogenetics:</b>	5q23.2
<b>Protein Families:</b>	Druggable Genome, Protein Kinase
<b>Protein Pathways:</b>	Hedgehog signaling pathway
<b>MW:</b>	48.9 kDa
<b>Gene Summary:</b>	<p>This gene encodes a member of a family of serine/threonine protein kinases that phosphorylate caseins and other acidic proteins. A related protein in the African clawed frog participates in the transmission of Wnt/beta-catenin signaling. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene. [provided by RefSeq, Jul 2012]</p> <p>Transcript Variant: This variant (2) differs in the 5' UTR and has multiple differences in the coding region, but maintains the reading frame, compared to variant 1. The encoded isoform (2) is shorter than isoform 1. Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.</p>