

Product datasheet for **SC302638**

GGT1 (NM_001032365) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	GGT1 (NM_001032365) Human Untagged Clone
Tag:	Tag Free
Symbol:	GGT1
Synonyms:	CD224; D22S672; D22S732; GGT; GGT 1; GTG
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)




[View online »](#)

Fully Sequenced ORF: >SC302638 representing NM_001032365.
Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGAAGAAGAAGTTAGTGGTGCTGGGCTGCTGGCCGTGGTCCTGGTGCTGGTCATTGTCCGGCTCTGT
CTCTGGCTGCCCTCAGCCTCAAGGAACCTGACAACCATGTGTACACCAGGGCTGCCGTGGCCCGGAT
GCCAAGCAGTGCTCGAAGATTGGGAGGGATGCACTGCGGGACGGTGGCTCTGCGGTGGATGCAGCCATT
GCAGCCCTGTTGTGTGTGGGGCTCATGAATGCCACAGCATGGGCATCGGGGGTGGCCTCTTCTCACC
ATCTACAACAGCACACACGAAAAGCTGAGGTCATCAACGCCCGCAGGTGGCCCCAGGCTGGCCTTT
GCCACCATGTTCAACAGCTCGGAGCAGTCCCAGAAGGGGGGGTGTGGTGGCGGTGCCTGGGGAGATC
CGAGGCTATGAGCTGGCACACCAGCGCATGGGCGGTGCCCTGGGCTCGCCTCTCCAGCCAGCATC
CAGCTGGCCCGCAGGGCTTCCCGTGGCAAGGGCTTGGCGGCAGCCCTGAAAAACAAGCGGACCGTC
ATCGAGCAGCAGCTGTCTTGTGTGAGGTGTTCTGCCGGATAGAAAGGTGCTTCGGGAGGGGGAGAGA
CTGACCCTGCCGAGCTGGCTGACACCTACGAGACGCTGGCCATCGAGGGTGCCAGGCCTTCTACAAC
GGCAGCCTCACGCCCAGATTGTGAAGGACATCCAGGCGCCGGGGGATTGTGACAGCTGAGGACCTG
AACTACTACCGTCTGAGCTGATCGAGCACCCGCTGAACATCAGCCTGGGAGACGTGGTGTGTACATG
CCCAGTGCGCCGCTCAGCGGGCCCGTGCTGGCCCTCATCTCAACATCCTCAAAGGTACAACCTTCTCC
CGGGAGAGCGTGGAGAGCCCCGAGCAGAAGGGCCTGACGTACCACCGCATCGTAGAGGCTTTCGGTTT
GCCTACGCCAAGAGGACCCTGCTTGGGGACCCCAAGTTTGTGGATGTGACTGAGGTGGTCCGCAACATG
ACCTCCGAGTTCTTCGCTGCCAGCTCCGGGCCAGATCTCTGACGACACCACTACCCGATCTCCTAC
TACAAGCCCGAGTTCTACACGCCGATGACGGGGCAGTGTACCTGTCTGTGTCGTCAGAGGACGGC
AGTGTGTGTCCGCCACCAGCACCATCAACCTCTACTTTGGCTCCAAGGTCCGCTCCCGGTGAGCGGG
ATCCTGTTCAATAATGAAATGGACGACTTCAGCTCTCCAGCATACCAACGAGTTTGGGTACCCCC
TCACTGCCAATTTATCCAGCCAGGAAGCAGCCGCTCTCGTCCATGTGCCCGACGATCATGGTGGGC
CAGGACGGCCAGGTCCGGATGGTGGTGGGAGCTGCTGGGGGCACACAGATCACCACGGCCACTGCACTG
GCCATCATCTACAACCTCTGGTTCGGCTATGACGTGAAGCGGGCCGTGGAGGAGCCCCGGCTGCACAAC
CAGCTTCTGCCAACGTACGACAGTGGAGAGAAAATTGACCAGGCAGTGACTGCAGCCCTGGAGACC
CGGCACCATCACACCAGATCGCGTCCACCTTCATCGCTGTGGTGAAGCCATCGTCCGCACGGCTGGT
GGCTGGGAGCTGCTCGGACTCCAGAAAGCGGGGAGCCTGCCGGCTACTGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites: SgfI-MluI

Plasmid Map:


ACCN: NM_001032365

Insert Size: 1710 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001032365.2](#)

RefSeq Size: 2393 bp

RefSeq ORF: 1710 bp

Locus ID: 2678

Cytogenetics: 22q11.23

Protein Families: Protease, Transmembrane

Protein Pathways:	Arachidonic acid metabolism, Cyanoamino acid metabolism, Glutathione metabolism, Metabolic pathways, Selenoamino acid metabolism, Taurine and hypotaurine metabolism
MW:	61.4 kDa
Gene Summary:	<p>The enzyme encoded by this gene is a type I gamma-glutamyltransferase that catalyzes the transfer of the glutamyl moiety of glutathione to a variety of amino acids and dipeptide acceptors. The enzyme is composed of a heavy chain and a light chain, which are derived from a single precursor protein. It is expressed in tissues involved in absorption and secretion and may contribute to the etiology of diabetes and other metabolic disorders. Multiple alternatively spliced variants have been identified. There are a number of related genes present on chromosomes 20 and 22, and putative pseudogenes for this gene on chromosomes 2, 13, and 22. [provided by RefSeq, Jan 2014]</p> <p>Transcript Variant: This variant (5) differs in the 5' UTR compared to variant 1. Variants 1, 3, 4, and 5 all encode the same protein.</p>