

Product datasheet for **SC302640**

ZNF226 (NM_001032372) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ZNF226 (NM_001032372) Human Untagged Clone
Tag:	Tag Free
Symbol:	ZNF226
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)




[View online »](#)

Fully Sequenced ORF: >SC302640 representing NM_001032372.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCACGATCGCC
ATGAATATGTTCAAGGAAGCAGTGACCTTCAAGGACGTGGCTGTGGCCTTCACGGAGGAGGAATTGGGG
GTGCTGGGCCCTGCCAGAGGAAGCTGTACCGAGATGTGATGGTGGAGAAGCTTTAGAACCTGCTGTCA
GTGGGGCATCCACCCTTCAAACAAGATGTATCACCTATAGAAAAGAAATGAGCAGCTTTGGATAATGACG
ACAGCAACCCGAAGACAGGGAAATTTAGGAGAGAAAAATCAAAGTAAGTTAATTACTGTTCAAGACAGA
GAATCAGAAGAAGAGCTTTCTTGTGGCAAATCTGGCAACAAATGCAAATGACTTAACCAGGTGTCAA
GACTCCATGATCAATAATTCTCAGTGTCAACAAAGGTGATTTCCCTTACCAGGTAGGGACAGAAGCTG
TCTATTCAAATTTCTGAAGATGAGAAGTATATAGTAAATAAAGCAGATGGTCCCAATAATACTGGGAAT
CCAGAGTTTCTATCTTGAGAACCCAGGATTCTTGAGGAAAACATTCTGACTGAGTACAGAGATTG
AACAGAGATCAGCAAATTTCCATAAAAAATAAATTATGTCAATGAAGAAGGGTGTGATCCCATCGGT
TGGATTTACATCATGATGGTCAAGAGTACACAAAAGTAAAAATCTTATAGACCAATGATTATGAA
AAAGACAACATGAAGATTTTGACATTTGATCACAATAGCATGATTCACACAGGACAGAAATCGTACCAG
TGTAATGAGTGTAAAAACCCCTTCAGTGATCTCTCCAGCTTTGATCTTCATCAGCAGTTACAATCAGGA
GAGAAGTCTTTACATGTGTTGAGCGTGGAAAAGGCTTCTGTTACAGCCAGTTCTTCTGTTTCATCAG
AAAGTACATGTGGGAGAAAAACTTAAGTGTGATGAGTGTGGTAAGGAATTCAGTCAGGGCGCTCATCTA
CAGACCCATCAGAAAGTCCACGTGATAGAGAAACCATACAAATGTAAGCAATGTGGGAAAGGTTTCAGT
CGTAGATCAGCACTTAATGTTTCATTGCAAGGTCCACACGGCAGAGAAACCTTATAATTGTGAGGAGTGT
GGGAGGGCCTTCAGTCAGGCCTCTCATCTTCAGGACCATCAGAGACTCCACACTGGGAGAGGCCATTC
AAATGTGATGCATGTGGTAAGAGCTTCAGTCGGAATTCACATCTCAATCCCATCAAAGATTCATACA
GGAGAGAAACCATACAAATGTGAGGAGTGTGGTAAGGGCTTCATTTGAGCTCAAATCTTTACATTCAT
CAGAGAGTCCACACAGGAGAAAAACCTATAAATGTGAGGAATGTGGTAAAGGCTTTAGTCGGCCTTCA
AGTCTTCAGGCCATCAGGGAGTTCACACTGGAGAGAAGTACATACATATGTAAGTGTGTTGAGGAAAGGC
TTTACTCTGAGTTCAAATCTTCAAGCCATCAGAGAGTCCACACTGGAGAGAAGCCATACAAATGCAAT
GAGTGTGGGAAGAGCTTCAGGAGGAATCCCATTAAGTTCATCTAGTGGTCCACACAGGAGAGAGAAA
CCCTATAAATGTGAGATATGTGGGAAGGGCTTCAGTCAAAGTTCGTATCTTCAAATCCATCAGAAGGCC
CACAGTATAGAGAAACCTTTAAGTGTGAGGAGTGTGGCAGGGTTTCAATCAGAGCTCAGACTTCAG
ATTCACCAGCTGATCCATACGGGTGAGAAACCATACAAATGTGAAGAGTGTGGCAAGGGATTTAGTCGT
AGAGCAGATCTTAAATCACTGTAGGATCCACACAGGAGAGAAACCATATAATTGTGAGGAGTGTGGG
AAGGTCTTCAGGCAGGCCTCAAATCTTTTGGCCATCAGAGAGTCCACAGTGGAGAAAAACCATTCAA
TGTGAAGAATGTGGGAAGAGTTTCGGTCCGAGTGCACATCTTCAAGCCATCAAAGTCCACACTGGGA
GATAAGCCATACAAATGTGATGAGTGTGGGAAGGGCTTCAAGTGGAGCTTGAACCTTGACATGCATCAG
AGGGTGCACACAGGAGAGAAAAACCATATAAATGTGGGAGTGTGGTAAGTACTTCAGTCAGGCCTCAAGT
CTTCAACTTCATCAGAGTGTCCACACAGGAGAGAAACCATACAAATGTGATGTGTGGTAAAGTCTTC
AGTCGGTCTTCAACTACAGTCTCATCAGCGAGTTTACACTGGGAGAGAAACCTTATAAATGTGAGATA
TGTGGTAAGAGCTTCAGTTGGCGATCAAATCTTACAGTTCATCACAGAATCCATGTTGGTGATAAATCC
TATAAAAGTAATAGGGGTGGTAAGAACATCAGAGAATCCACACAGGAAAAAATCTATAAAATGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: SgfI-MluI

Plasmid Map:


ACCN: NM_001032372

Insert Size: 2412 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001032372.1](#)

RefSeq Size: 2623 bp

RefSeq ORF: 2412 bp

Locus ID: 7769

UniProt ID: [Q9NNT6](#)

Cytogenetics: 19q13.31

Protein Families: Transcription Factors

MW: 91.9 kDa

Gene Summary: May be involved in transcriptional regulation.[UniProtKB/Swiss-Prot Function]
Transcript Variant: This variant (1) encodes the longer isoform (a). Variants 1, 2, 7, 8, 9, and 10 all encode the same isoform (a).