

Product datasheet for **SC302674**

FBX07 (NM_001033024) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	FBX07 (NM_001033024) Human Untagged Clone
Tag:	Tag Free
Symbol:	FBX07
Synonyms:	FBX; FBX07; FBX7; PARK15; PKPS
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC302674 representing NM_001033024. Blue=Insert sequence Red=Cloning site Green=Tag(s)


```
GCTCGTTTGTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGCCCGGCCTCCCGGGGCTCTGGTCCCTCCTCGATTAGAGCATTCTTCACTCCAGAATAATGAG
CAACCCTCTTTGGCCACCAGCTCCAATCAGACTAGCATGCAGGATGAACAACCAAGTGATTCATTCCAA
GGACAGGCAGCCAGTCTGGTGTGGAAATGACGACAGTATGTTAGGGCCTAGTCAAAATTTGAAGCT
GAGTCAATTCAGATAATGCGCATATGGCAGAGGGCACAGTTTCTATCCCTCAGAACCCATGCTCTGT
AGTGAATCGGTGGAAGGGCAAGTCCACATTCATTAGAGACCTTGTATCAATCAGCTGACTGTTCTGAT
GCCAATGATGCCCTTGATAGTGTGATACATCTTCTCATGTTGGAGTCAGGTTACATACCTCAGGGCACC
GAAGCCAAAGCACTGTCCATGCCGAGAAGTGAAGTTGAGCGGGGTGATAAGCTGCAGTACATGCAT
CCTCTCTGCGAGGGCAGCTCCGCTACTCTCACCTGTGTGCCTTTGGGAAACCTGATTGTTGAAATGCT
ACACTAAAAATCAACAATGAGATTAGAAGTGTGAAAAGATTGCAGCTGCTACCAGAATCTTTATTTGC
AAAGAGAAACTAGGGGAAAATGTAGCCAACATATACAAAGATCTTCAGAAACTCTCTCGCCTCTTTAAA
GACCAGCTGGTGTATCCTCTTCTGGCTTTTACCCGACAAGCACTGAACCTACCAGATGATTTGGGTTG
GTCGTCTCCATTGGAAGTAACTACGGATCTCCGACTTCTGGATGTTCCGTTCCGCTTGTCTTTG
TCTGCGGTTTGTGCGTACCTCTTACTGCTTCAAATGACCACTCCTGTGGAGGTTTTATATCTGCGT
GATTTTCGAGACAATACTGTGAGATTCAAGACACAGATTGGAAGAAGTGTACAGGAAGAGGCACATA
CAAAGAAAAGAATCCCCGAAAGGGCGGTTTGTGATGCTCCTGCCATCGTCAACTCACACCATTCATTTC
TATCCCAACCCCTTGACCCCTAGGCCATTTCTAGCTCCCGCCTTCTCCAGGAATTATCGGGGGTAA
TATGACCAAAGACCAACACTTCCCTATGTTGGAGACCAATCAGTTCACTCATTCTGGTCTGGGGAG
ACGCCCAGCCAGTTTCTCCACTGAGACCAGCTTTGATCCAGTTGGCCCACTTCCAGGACCTAACCCC
ATCTTGCCAGGGCGAGGGCGCCCAATGACAGATTTCCCTTTAGACCCAGCAGGGGTGCGCCAAGTAT
GGCCGGCTGTATTGATGA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTAAACGGCCGCGC
```



[View online >](#)

Restriction Sites: SgfI-MluI

Plasmid Map:



ACCN: NM_001033024

Insert Size: 1332 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001033024.1](#)

RefSeq Size: 1758 bp

RefSeq ORF: 1332 bp

Locus ID: 25793

UniProt ID: [Q9Y311](#)

Cytogenetics: 22q12.3

Protein Families: Druggable Genome

MW: 49.4 kDa

Gene Summary: This gene encodes a member of the F-box protein family which is characterized by an approximately 40 amino acid motif, the F-box. The F-box proteins constitute one of the four subunits of the ubiquitin protein ligase complex called SCFs (SKP1-cullin-F-box), which function in phosphorylation-dependent ubiquitination. The F-box proteins are divided into 3 classes: Fbws containing WD-40 domains, Fbls containing leucine-rich repeats, and Fbxs containing either different protein-protein interaction modules or no recognizable motifs. The protein encoded by this gene belongs to the Fbxs class and it may play a role in regulation of hematopoiesis. Alternatively spliced transcript variants of this gene have been identified with the full-length natures of only some variants being determined. [provided by RefSeq, Jul 2008]
Transcript Variant: This variant (2) contains a different 5' UTR and 5' coding region, compared to variant 1. This results in a shorter protein (isoform 2) with a distinct N-terminus, compared to isoform 1.