

Product datasheet for **SC302718**

cleavage stimulation factor (CSTF1) (NM_001033521) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	cleavage stimulation factor (CSTF1) (NM_001033521) Human Untagged Clone
Tag:	Tag Free
Symbol:	CSTF1
Synonyms:	CstF-50; CstFp50
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC302718 representing NM_001033521. Blue=Insert sequence Red=Cloning site Green=Tag(s)


```
GCTCGTTTGTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGTACAGAACCAAGTGGGCTTGAAGGACCGCCAGCAGCTCTACAAGCTGATCATTAGCCAGCTGCTA
TATGACGGCTACATCAGCATCGCCAATGGCCTCATCAATGAAATCAAGCCTCAGTCTGTGTGTGCACCC
TCGGAGCAGTCTCTGCATCTCATCAAACCGGAATGGAAAACGATGACACCGCAGTTCAGTATGCAATT
GGTCGTTTCAAGTACTGTTGCCCTGGCACAGGGATTGACCTGGAATTTGATGCAGATGTTTCAAGTATG
TCCCCAGAGGCTTCTGAGTACGAAACATGCTATGTCACATCACATAAAGGACCATGCCGTGTAGCTACC
TATAGTAGAGATGGACAGTTAATAGCTACTGGTCTGCTGATGCTTCGATAAAGATACTTGACACAGAG
AGGATGTTGGCCAAAAGTGCCATGCCAATAGAGGTCATGATGAATGAGACCGCACAAACAAATATGGAA
AACCACCCAGTGATTCGAACTCTTTATGACCATGTGGATGAAGTACGTCGCTTGTCTTCCACCCAAACA
GAACAGATCCTGGCTTCTGGTCAAGGGATTACTCTTAAATATTTGATTATTTCAAACCATCAGCA
AAAAGAGCCTTCAAATACATTCAAGGAGCTGAAATGTTACGTTCCATCTCTTTTCATCCTTCTGGAGAC
TTTATACTTGTGCGAACTCAGCATCCTACTCTTCGCCTTTATGATATCAACACCTTTCAATGTTTGTGTC
TCTTGAATCCTCAAGATCAACACCCGATGCTATATGTTCCGTTAATTACAATTCTAGTGCCAATATG
TACGTAACCTGGAAGCAAGGACGGCTGCATCAAATATGGGATGGTGTTCAAATCGATGCATCACAAT
TTTGAGAAAGCAGATGACGGTGTGAAGTTTGTCTGCCATTTTTTCCAAAAATCTAAATACATTCTC
TCAAGTGGAAAAGACTCTGTAGCTAAACTTTGGGAAATATCAACGGGACGAACACTGGTCAGATACACG
GGCGCGGGTTTAAAGTGGACGCCAGGTGCACCGGACACAGGCTGTGTTTAAACCACCCGAGGACTATGTG
TTGCTGCCGACGAGAGGACGATCAGTCTTTGCTGCTGGGACTCGAGGACAGCCGAGCGGAGAAACCTG
CTGTCGTTGGGGCACAACAATATTGTACGCTGCATAGTGCCTCCCCACCAACCCCGGTTTCATGACG
TGCAGCGATGACTTCAGAGCGCGGTTTTGGTACCGGAGATCGACCACTGACTGA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```



[View online »](#)

Restriction Sites: SgfI-MluI

Plasmid Map:



ACCN: NM_001033521

Insert Size: 1296 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001033521.1](#)

RefSeq Size: 2295 bp

RefSeq ORF: 1296 bp

Locus ID: 1477

UniProt ID: [Q05048](#)

Cytogenetics: 20q13.2-q13.31

MW: 48.4 kDa

Gene Summary: This gene encodes one of three subunits which combine to form cleavage stimulation factor (CSTF). CSTF is involved in the polyadenylation and 3'end cleavage of pre-mRNAs. Similar to mammalian G protein beta subunits, this protein contains transducin-like repeats. Several transcript variants with different 5' UTR, but encoding the same protein, have been found for this gene. [provided by RefSeq, Jul 2008]
Transcript Variant: This variant (1) has a 5' UTR that is further upstream than the other variants. Variants 1, 2 and 3 encode the same protein.