

Product datasheet for **SC302722**

TOM1L2 (NM_001033551) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	TOM1L2 (NM_001033551) Human Untagged Clone
Tag:	Tag Free
Symbol:	TOM1L2
Vector:	<u>pCMV6 series</u>
Fully Sequenced ORF:	>NCBI ORF sequence for NM_001033551, the custom clone sequence may differ by one or more nucleotides ATGGAGTTCCTCCTGGGGAACCCGTTTCAGCACACCAGTGGGGCAGTGCCTCGAAAAGGCA ACAGATGGCTCCCTGCAAAGTGAGGATTGGACGTTGAATATGGAGATCTGTGACATCATC AATGAGACGGAGGAAGGGCCAAAGGATGCCATTCGAGCCCTGAAGAAGCGGCTCAACGGG AACCGAACTACAGAGAGGTGATGCTGGCATTAAACAGCATGGGCTGATGCCTTTGCAAGC AGTCTGATCTCACCGCGTGTGTGCACATATATGAGGAGCTGAAGAGGAAAGGGTTGAA TTTCCCATGGCAGACTTGGACGCTCTGTCTCCCATACACACACCACAGCGGAGTGTCCCT GAAGTGGATCCAGCTGCGACCATGCCAGTCCCAATCACAGCAGAGGACAAGTGCCTGGT TCCTATTCTCGCCGCTCCTGCTCCCTACTCCGCACCGCAGGCCAGCTCTGAGTGTG ACTGGCCCATCACAGCAATTCAGAACAGATTGCCAGGCTGCGGAGTGAAGTGGACGTC GTTGAGGAAACACAAAAGTCATGTCTGAGATGTTAACAGAAATGGTCCCTGGACAGGAG GATTCATCTGATCTGGAGTTGCTGCAGGAGCTCAACAGGACCTGTCGGGCCATGCAGCAG CGCATCGTGGAGCTCATCTCCCGCTGTCCAATGAGGAGGTCACCGAGGAGCTGCTGCAT GTGAACGATGACCTCAACAACGTTCTTCTCGATACGAGAGGTTGCAACGATACAGGCTC GGCCGATCCGTTCAAAATGCCAGTAATGGAGTACTGAATGAAGTAACCGAAGACAACCTTA ATAGACCTGGGGCCAGGGTCTCCAGCCGTGGTGAAGCCCAATGGTGGGGAACACAGCGCCC CCATCTTCCCTCCTCCAGCTTGCAGGCTTAGACTTGGGACAGAGAGCGTCAGTGGC ACCCTCAGTTCCTCAGCAATGTAATCCCGTGACGGCTTTGACATGTTTGCCAGACG AGAGGAACTCCTTGGCTGAGCAGCGCAAGACGGTAACCTATGAGGATCCTCAGGCTGTC GGAGGACTTGCTTCTGCACTAGACAATCGAAAACAGAGTTCAGAAAGGGATCCCCGTTGCG CAGCCATCTGTGATGGACGACATTGAGGTGTGGCTCAGGACCGACCTGAAGGGTGTGAT CTGGAGGAGGGTGTCACAAGTGAAGAGTTTGATAAATTCCTTGAAGAAAGAGCCAAAGCT GCTGAAATGGTCCCGACCTCCCTCGCCCCCATGGAGGCTCCTGCCCCAGCCTCAAAC CCTTCTGGCCGGAAGAAGCCAGAGCGGTCAGAGGATGCCCTCTTCGCCCTGTGA
Restriction Sites:	Please inquire
ACCN:	NM_001033551



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001033551.1</u> , <u>NP_001028723.1</u>
RefSeq Size:	2296 bp
RefSeq ORF:	1374 bp
Locus ID:	146691
UniProt ID:	<u>Q6ZVM7</u>
Cytogenetics:	17p11.2
Gene Summary:	<p>This gene belongs to a small gene family whose members have an N-terminal VHS domain followed by a GAT domain; domains which typically participate in vesicular trafficking. The canonical protein encoded by this gene also has a C-terminal clathrin binding motif. This protein has been shown to interact with Tollip, clathrin and ubiquitin and is thought to play a role in endosomal sorting. This gene resides in the 3.7 Mb deletion of chromosome region 17p11.2 that is associated with Smith-Magenis syndrome. Alternative splicing results in multiple transcript variants encoding distinct proteins. [provided by RefSeq, Apr 2017]</p> <p>Transcript Variant: This variant (1) lacks an alternate in-frame exon in the 5' coding region, compared to variant 3. It encodes isoform 1, which is shorter than isoform 3.</p>