

Product datasheet for SC302726

PPM1B (NM_001033556) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	PPM1B (NM_001033556) Human Untagged Clone
Tag:	Tag Free
Symbol:	PPM1B
Synonyms:	PP2C-beta-X; PP2CB; PP2CBETA; PPC2BETAX
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC302726 representing NM_001033556. Blue=Insert sequence Red=Cloning site Green=Tag(s)


```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGGTGCATTTTTGGATAAACCCAAAACCTGAAAAACATAATGCTCATGGTGCTGGGAATGGTTTACGT
TATGGCCTGAGCAGCATGCAAGGATGGAGAGTGGAAATGGAAGATGCACACACAGCTGTTGTAGGTATT
CCTCACGGCTTGAAGACTGGTCATTTTTTGCAGTTTATGATGGTCATGCTGGATCCCGAGTGGCAAAT
TACTGCTCAACACATTTATTAGAACACATCACTACTAACGAAGACTTTAGGGCAGCTGGAAAATCAGGA
TCTGCTCTTGAGCTTTCAGTGGAAAATGTTAAGAATGGTATCAGAAGTGGATTTTTGAAAATTGATGAA
TACATGCGTAACTTTTAGACCTCAGAAACGGGATGGACAGGAGTGGTTCACTGCAGTGGGAGTTATG
ATTTACCTAAGCATATCTACTTTATCAACTGTGGTGATTCACGTGCTGTTCTGTATAGGAATGGACAA
GTCTGCTTTTCTACCCAGGATCACAACCTTGCAATCCAAGGGAAAAGGAGCGAATCCAAAATGCAGGA
GGCAGCGTGATGATACAACGTGTTAATGGTTCATTAGCAGTATCTCGTGCTCTGGGGACTATGATTAC
AAGTGTGTTGATGGCAAGGGCCCAACAGAACTTGTCTCCAGAGCCTGAGGTTTATGAAAATTTA
AGAGCAGAAGAGGATGAATTTATCATCTTGCTTGTGATGGGATCTGGGATGTTATGAGTAAATGAGGAG
CTCTGTGAATATGTTAAATCTAGGCTTGAGGTATCTGATGACCTGGAAAATGTGTGCAATTGGGTAGTG
GACACTTGTTTACACAAGGGAAGTCGAGATAACATGAGTATTGTACTAGTTTGCCTTTTCAAATGCCTCC
AAGGTCTCAGATGAAGCGGTGAAAAAAGATTGAGAGTTGGATAAGCACTTGAATCACGGGTTGAAGGT
AAGACAAATGCTTTTAA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites: SgfI-MluI



[View online »](#)

Plasmid Map:



ACCN: NM_001033556

Insert Size: 984 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001033556.1](#)

RefSeq Size: 1497 bp

RefSeq ORF: 984 bp

Locus ID: 5495

Cytogenetics:	2p21
Protein Families:	Druggable Genome, Phosphatase, Stem cell - Pluripotency
Protein Pathways:	MAPK signaling pathway
MW:	36.3 kDa
Gene Summary:	<p>The protein encoded by this gene is a member of the PP2C family of Ser/Thr protein phosphatases. PP2C family members are known to be negative regulators of cell stress response pathways. This phosphatase has been shown to dephosphorylate cyclin-dependent kinases (CDKs), and thus may be involved in cell cycle control. Overexpression of this phosphatase is reported to cause cell-growth arrest or cell death. Alternative splicing results in multiple transcript variants encoding different isoforms. Additional transcript variants have been described, but currently do not represent full-length sequences. [provided by RefSeq, Jul 2008]</p> <p>Transcript Variant: This variant (4) uses an alternate in-frame splice site in the 5' coding region, compared to variant 1. The resulting protein (isoform 4) has a distinct and shorter C-terminus, compared to isoform 1.</p>