

## Product datasheet for **SC302734**

### **MIRO1 (RHOT1) (NM\_001033566) Human Untagged Clone**

#### **Product data:**

Product Type:	Expression Plasmids
Product Name:	MIRO1 (RHOT1) (NM_001033566) Human Untagged Clone
Tag:	Tag Free
Symbol:	RHOT1
Synonyms:	ARHT1; MIRO-1; MIRO1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)




[View online »](#)

Fully Sequenced ORF: >SC302734 representing NM\_001033566.  
Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGAAGAAAGACGTGCGGATCCTGCTGGTGGGAGAACCCTAGAGTTGGGAAGACATCACTGATTATGTCT
CTGGTCAGTGAAGAATTCAGAAAGAGGTTCTCCCGGGCAGAAGAAATCACCATTCCAGCTGATGTC
ACCCAGAGAGAGTTCCAACACACATTGTAGATTACTCAGAAGCAGAACAGAGTGATGAACAACTTCAT
CAAGAAATATCTCAGGCTAATGTCATCTGTATAGTGTATGCCGTTAAACAACAAGCATTCTATTGATAAG
GTAACAAGTCGATGGATTCTCTCATAAATGAAAGAACAGACAAAGACAGCAGGCTGCCTTTAATATTG
GTTGGGAACAAATCTGATCTGGTGAATATAGTAGTATGGAGACCATCCTTCCTATTATGAACCAGTAT
ACAGAAATAGAAACCTGTGTGGAGTGTTGAGCGAAAAACCTGAAGAACATATCAGAGCTCTTTATTAC
GCACAGAAAGCTGTTCTTCATCTACAGGGCCCTGTACTGCCAGAGGAGAAGGAGATGAAACCAGCT
TGTATAAAAGCCCTTACTCGTATATTTAAAATATCTGATCAAGATAATGATGGTACTCTCAATGATGCT
GAACTCAACTTCTTCAGAGGATTTGTTTCAACTCCATTAGCTCCTCAAGCTCTGGAGGATGTCAAG
AATGTAGTCAGAAAACATATAAGTGATGGTGTGGCTGACAGTGGGTTGACCCTGAAAGGTTTTCTCTTT
TTACACACACTTTTTATCCAGAGAGGGAGACACGAAACTACTTGGACTGTGCTTCGACGATTTGGTTAT
GATGATGACCTGGATTTGACACCTGAATATTTGTTCCCCCTGCTGAAAATACCTCCTGATTGCACTACT
GAATTAATCATCATGCATATTTATTTCTCAAAGCACCTTTGACAAGCATGATTTGGATAGAGACTGT
GCTTTGTCACCTGATGAGCTTAAAGATTTTAAAGTTTTCCCTTACATACCTTGGGGGCCAGATGTG
AATAACACAGTTTGTACCAATGAAAGAGGCTGGATAACCTACCAGGGATTCCTTTCCAGTGGACGCTC
ACGACTTATTTAGATGTACAGCGGTGCCGGAATATTTGGGCTATCTAGGCTATTCAATATTGACTGAG
CAAGAGTCTCAAGCTTCAGCTGTTACAGTGACAAGAGATAAAAAGATAGACCTGCAGAAAAACAACT
CAAAGAAATGTGTTCCAGATGTAATGTAATTGGAGTGAAAAACTGTGGGAAAAGTGGAGTTCTTCAGGCT
CTTCTTGAAGAAACTTAATGAGGCAGAAGAAAATTCGTGAAGATCATAAATCCTACTATGCGATTAAC
ACTGTTTATGTATATGGACAAGAGAAAATACTTGTGTTGCATGATATCTCAGAATCGGAATTTCTAACT
GAAGCTGAAATCATTGTTGATGTTGATGCCTGGTATATGATGTCAGCAATCCCAAATCCTTTGAATAC
TGTGCCAGGATTTTAAAGCAACTTTATGGACAGCAGAATACCTTGCTTAATCGTAGCTGCAAAGTCA
GACCTGCATGAAGTTAAACAAGAATACAGTATTTACCTACTGATTTCTGCAGGAAACACAAAATGCCT
CCACCACAAGCCTTCACTTGCAATACTGCTGATGCCCCAGTAAGGATATCTTTGTTAAATTGACAACA
ATGGCCATGTATCCCATGCCCGTTACGCTGTATGTGCACCTGCAACAGGTGTACATTTTGCATCTGT
CAGAACTTCTCAACTCAGACTTGTGCAATCTGTAAGAACAAAATCTTCACTGCAGTTCTTAACAGG
CACGTGACACAAGCTGACCTCAAGAGCTCCACGTTTTGGCTTCGAGCAAGTTTTGGTGCTACTGTTTTT
GCAGTTTTGGGCTTTGCTATGTACAAAGCATTATTGAAACAGCGATGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites: SgfI-MluI

**Plasmid Map:**


**ACCN:** NM\_001033566

**Insert Size:** 1980 bp

**OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

**OTI Annotation:** This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_001033566.2](#)

**RefSeq Size:** 3263 bp

**RefSeq ORF:** 1980 bp

**Locus ID:** 55288

**UniProt ID:** [Q8IXI2](#)

**Cytogenetics:** 17q11.2

<b>Protein Families:</b>	Transmembrane
<b>MW:</b>	75.5 kDa
<b>Gene Summary:</b>	<p>Mitochondrial GTPase involved in mitochondrial trafficking. Probably involved in control of anterograde transport of mitochondria and their subcellular distribution.[UniProtKB/Swiss-Prot Function]</p> <p>Transcript Variant: This variant (2) lacks an alternate in-frame exon in the 3' coding region, compared to variant 1. It encodes isoform 2, which is shorter than isoform 1.</p>