

## **Product datasheet for SC302748**

## PKC zeta (PRKCZ) (NM 001033581) Human Untagged Clone

**Product data:** 

**Product Type:** Expression Plasmids

Product Name: PKC zeta (PRKCZ) (NM\_001033581) Human Untagged Clone

Tag: Tag Free
Symbol: PRKCZ

**Synonyms:** PKC-ZETA; PKC2

Mammalian Cell

Selection:

Neomycin

Vector: pCMV6-Entry (PS100001)

E. coli Selection: Kanamycin (25 ug/mL)

Fully Sequenced ORF: >SC302748 representing NM\_001033581.

Blue=Insert sequence Red=Cloning site Green=Tag(s)

GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC

GAGGAGACAGATTGCTTACATTTCCTCATCCCGGAAGCATGACAGCATTAAAGACGACTCGGAG GACCTTAAGCCAGTTATCGATGGGATGGATGGAATCAAAATCTCTCAGGGGCTTGGGCTGCAGGACTTT GACCTAATCAGAGTCATCGGGCGCGGGAGCTACGCCAAGGTTCTCCTGGTGCGGTTGAAGAAGAATGAC CAAATTTACGCCATGAAAGTGGTGAAGAAGAGCTGGTGCATGACGAGGATATTGACTGGGTACAG ACAGAGAAGCACGTGTTTGAGCAGGCATCCAGCAACCCCTTCCTGGTCGGATTACACTCCTGCTTCCAG ACGACAAGTCGGTTGTTCCTGGTCATTGAGTACGTCAACGGCGGGGACCTGATGTTCCACATGCAGAGG CAGAGGAAGCTCCCTGAGGAGCACGCCAGGTTCTACGCGGCCGAGATCTGCATCGCCCTCAACTTCCTG CACGAGAGGGGGATCATCTACAGGGACCTGAAGCTGGACAACGTCCTCCTGGATGCGGACGGGCACATC AAGCTCACAGACTACGGCATGTGCAAGGAAGGCCTGGGCCCTGGTGACACAACGAGCACTTTCTGCGGA ACCCCGAATTACATCGCCCCCGAAATCCTGCGGGGAGAGGAGTACGGGTTCAGCGTGGACTGGTGGGCG CTGGGAGTCCTCATGTTTGAGATGATGGCCGGGCGCTCCCCGTTCGACATCATCACCGACAACCCGGAC ATGAACACAGAGGACTACCTTTTCCAAGTGATCCTGGAGAAGCCCATCCGGATCCCCGGTTCCTGTCC GTCAAAGCCTCCCATGTTTTAAAAGGATTTTTAAATAAGGACCCCAAAGAGAGGCTCGGCTGCCGGCCA CAGACTGGATTTTCTGACATCAAGTCCCACGCGTTCTTCCGCAGCATAGACTGGGACTTGCTGGAGAAG AAGCAGGCGCTCCCTCCATTCCAGCCACAGATCACAGACGACTACGGTCTGGACAACTTTGACACACAG TTCACCAGCGAGCCCGTGCAGCTGACCCCAGACGATGAGGATGCCATAAAGAGGATCGACCAGTCAGAG

AGCGGACCGACGCGTACGCGGCCCCCCGAGCAGAAACTCATCTCAGAAGAGAGATCTGGCAGCAAATGAT

TTCGAAGGCTTTGAGTATATCAACCCATTATTGCTGTCCACCGAGGAGTCGGTGTGA

ATCCTGGATTACAAGGATGACGACGATAAGGTTTAA

**Restriction Sites:** Sgfl-Rsrll




**OriGene Technologies, Inc.** 9620 Medical Center Drive, Ste 200

CN: techsupport@origene.cn

Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com



## Plasmid Map:



ACCN: NM 001033581

**Insert Size:** 1230 bp

**OTI Disclaimer:** Due to the inherent nature of this plasmid, standard methods to replicate additional amounts

of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at <a href="mailto:customercom">customercom</a> or by

calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing

variants is recommended prior to use. <u>More info</u>

**OTI Annotation:** This TrueClone is provided through our Custom Cloning Process that includes sub-cloning

into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube

containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:** 1. Centrifuge at 5,000xg for 5min.

2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.

3. Close the tube and incubate for 10 minutes at room temperature.

4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid

at the bottom.

5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of

shipping when stored at -20°C.



**RefSeq:** <u>NM 001033581.1</u>

 RefSeq Size:
 2147 bp

 RefSeq ORF:
 1230 bp

 Locus ID:
 5590

 UniProt ID:
 Q05513

 Cytogenetics:
 1p36.33

**Protein Families:** Druggable Genome, Protein Kinase

**Protein Pathways:** Chemokine signaling pathway, Endocytosis, Insulin signaling pathway, Tight junction, Type II

diabetes mellitus

**MW:** 46.6 kDa

**Gene Summary:** Protein kinase C (PKC) zeta is a member of the PKC family of serine/threonine kinases which

are involved in a variety of cellular processes such as proliferation, differentiation and secretion. Unlike the classical PKC isoenzymes which are calcium-dependent, PKC zeta exhibits a kinase activity which is independent of calcium and diacylglycerol but not of phosphatidylserine. Furthermore, it is insensitive to typical PKC inhibitors and cannot be activated by phorbol ester. Unlike the classical PKC isoenzymes, it has only a single zinc finger module. These structural and biochemical properties indicate that the zeta subspecies is related to, but distinct from other isoenzymes of PKC. Alternative splicing results in multiple

transcript variants encoding different isoforms. [provided by RefSeq, Jul 2008]

Transcript Variant: This variant (2) differs in the 5' UTR and has multiple coding region differences, compared to variant 1. These differences cause translation initiation from a downstream ATG and an isoform (2) with a shorter N-terminus, compared to isoform 1.

Variants 2 and 3 encode the same isoform (2).