

## Product datasheet for **SC302764**

### TRPT1 (NM\_001033678) Human Untagged Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** TRPT1 (NM\_001033678) Human Untagged Clone  
**Tag:** Tag Free  
**Symbol:** TRPT1  
**Vector:** pCMV6 series  
**Fully Sequenced ORF:** >NCBI ORF sequence for NM\_001033678, the custom clone sequence may differ by one or more nucleotides

```
ATGAACTTCTCTGGAGGAGGGAGGCAGGAAGCAGCAGGGTCCAGGGGTAGAAGGGCTCCC
AGACCCCGAGAACAGGACCGAGACGTGCAGCTGTCCAAGGCTCTGTCTATGCCCTGCGC
CATGGGGCCTTGAAGCTGGGGCTTCCCATGGGAGCTGATGGCTTCGTGCCCTGGGCACC
CTCCTGCAGTTGCCCCAGTTCGCGGGCTTCTCTGCTGAAGATGTGCAGCGCGTGGTGGAC
ACCAATAGGAAGCAGCGGTTCCGCCCTGCAGCTGGGGGATCCCAGCACTGGCCTTCTCATC
CGGGCCAAACCAGGGCCATTCCCTGCAGGTACCTAAGTTGGAGCTGATGCCCTGGAGACA
CCGCAGGCCCTGCCCCGATGCTAGTCCATGGTACATTCTGGAAGCACTGGCCATCCATC
CTACTCAAAGGCCTGTCTCGCCAGGGAAGGACGCACATTACCTGGCCCCAGGACTGCCT
GGAGACCCCGGTATCATCAGTGGCATGCGGTCCCATTGTGAAATAGCTGTGTTTCATCGAT
GGACCCCTGGCTCTGGCAGATGGAATACCCTTCTTCCGCTCTGCCAATGGGGTGATTCTG
ACTCCAGGGAATACTGATGGCTTCTCCTTCCCAAGTACTTCAAGGAGGCCCTGCAGCTA
CGCCCTACCCGAAAGCCCTTTCCTTGGCTGGTGATGAAGAGACAGAGTGTGAGAGTAGC
CCCAAGCACAGCTCCAGAGAAAGGAGGAGGATCCAACAATAA
```

**Restriction Sites:** Please inquire  
**ACCN:** NM\_001033678  
**OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).  
**OTI Annotation:** This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.  
**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).



[View online »](#)

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_001033678.1](#), [NP\\_001028850.1](#)

**RefSeq Size:** 963 bp

**RefSeq ORF:** 762 bp

**Locus ID:** 83707

**UniProt ID:** [Q86TN4](#)

**Cytogenetics:** 11q13.1

**Gene Summary:** Catalyzes the last step of tRNA splicing, the transfer of the splice junction 2'-phosphate from ligated tRNA to NAD to produce ADP-ribose 1''-2'' cyclic phosphate.[UniProtKB/Swiss-Prot Function]