

## Product datasheet for **SC302765**

### **NOP2 (NM\_001033714) Human Untagged Clone**

#### **Product data:**

Product Type:	Expression Plasmids
Product Name:	NOP2 (NM_001033714) Human Untagged Clone
Tag:	Tag Free
Symbol:	NOP2
Synonyms:	NOL1; NOP120; NSUN1; p120
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

**Fully Sequenced ORF:** >SC302765 representing NM\_001033714.  
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGCATCGCC
ATGGGGCGCAAGTTGGACCTACGAAGGAGAAGCGGGGCCAGGCCGAAAGGCCCGAAGCAGAAGGGT
GCCGAGACAGAACTCGTCAGATTCTTGCCTGCAGTAAGTGACGAAAATTCCAAGAGGCTGTCTAGTCGT
GCTCGAAAAGAGGGCAGCCAAGAGGAGATTGGGCTCTGTTGAAGCCCTAAGACAAAATAAGTCTCCTGAG
GCCAAACCATTGCCTGGAAGCTACAAAAGGAGCTGTCCAGACAGCTGGTAAGAAGGGACCCAGTCC
CTATTTAATGCTCCTCGAGGCAAGAAGCGCCAGCACCTGGCAGTGATGAGGAAGAGGAGGGAAGAC
TCTGAAGAAGATGGTATGGTGAACCACGGGACCTCTGGGGCTCCGAGGACGATGCTGATACGGTAGAT
GACTATGGAGCTGACTCCAACCTGAGGATGAGGAGGAAGGTGAAGCGTTGCTGCCATTGAAAGAGCT
GCTCGAAGCAGAAGGCCGGGAAGCTGCTGCTGGGATCCAGTGGAGTGAAGAGGAGACCGAGGACGAG
GAGGAAGAGAAAAGTGAACCCCTGAGTCAGGCCCCCAAGGTGAAGAGGCAGATGGGGCTGCAG
ATCAATGTGGATGAGGAACATTTGTGCTGCCCTGCTGGGAGATGGAGCAGGATGCCAGGCTCCA
GACCTGCAACGAGTTCACAAGCGGATCCAGGATATTGTGGGAATTCTGCGTGATTTTGGGGCTCAGCGG
GAGGAAGGGCGGTCTCGTTCTGAATACCTGAACCGGCTCAAGAAGGATCTGGCCATTTACTACTCTAT
GGAGACTTCTGCTTGGCAAGCTCATGGACCTCTCCCTCTGCTGAGCTGGTGGAGTCTTAGAAGCT
AATGAGGTGCCTCGGCCCGTCACCTCCGGACCAATACCTTGAACCCGACGCCGAGACCTTGCACAG
GCTCTAATCAATCGTGGGGTTAACCTGGATCCCCTGGCAAGTGGTCAAAGACTGGACTAGTGGTGTAT
GATTCTTCTGTGCCCATTTGGTGTACCCCGAGTACCTGGCTGGGCACTACATGCTGCAGGGAGCCTCC
AGCATGTTGCCGTCATGGCCTTGGCACCCAGGAACATGAGCGGATCCTGGACATGTGTGTGCCCT
GGAGAAAAGACCAGCTACATGGCCAGCTGATGAAGAACACGGGTGATCCTTGCCAATGACGCCAAT
GCTGAGCGGCTCAAGAGTGTGGGGCAACTTGCATCGGCTGGGAGTACCAACACCAATTATCAGCCAC
TATGATGGGCGCCAGTTCCTCAAGGTGGTGGGGGCTTTGACCGAGTACTGCTGGATGCTCCCTGCAGT
GGCACTGGGGTCACTCCAAGGATCCAGCCGTGAAGACTAACAAGGATGAGAAGGACATCCTGCGCTGT
GCTCACCTCCAGAAGGAGTTGCTCCTGAGTGTATTGACTCTGTCAATGCGACCTCCAAGACAGGAGGC
TACCTGGTTTACTGCACCTGTTCTATCACAGTAGAAGAGAATGAGTGGTGGTAGACTATGCTCTGAAA
AAGAGGAATGTGCGACTGGTGCCACGGGCTAGACTTTGGCCAGGAAGGTTTTACCCGCTTTCGAGAA
AGGCGCTTCCACCCAGTCTGCGTCTACCCGACGTTCTACCCTCATACCCACAATATGGATGGGTTC
TTCATTGCCAAGTCAAGAAAATTTCCAATTCTATCCCTCAGTCCAGACAGGAAATTTGAAACAGCC
ACACCTACAAATGTAGACTTGCCTCAGGTCATCCCAAGTCTGAGAACAGCAGCCAGCCAGCCAAGAAA
GCCAAGGGGGCTGCAAAGACAAAGCAGCAGCTGCAGAAACAGCAACATCCAAGAAGGCCTCCTCCAG
AAGCTGAATGGCATCTCCAAGGGGCGAGCTCAGAATTGTCCACTGTACCTTCTGTACAAAAGACCCAA
GCTTCTCCAGCTTCCAGGATAGCAGTCAGCCAGCTGGAAGGCCGAAAGGATCAGGGAGCCAAAGGTG
ACTGGGAAGCTAAAGCAACGATCACCTAAATTACAGTCTCCAAGAAAGTTGCTTTCTCAGGCAGAAAT
GCCCTCCCAAGGGCACAGACACAAAACACCGGCTGTGTATCCCCATCAAGACTCAGGCCACCTG
AAACCTAAGGACCATCATCAGCCCTTGAAGGGCCAAGGGGTTGAGAAGCAGCAGTTGCCAGAGCAG
CCTTTTGAGAAAGCTGCCTTCCAGAAACAGAATGATACCCCAAGGGGCTCAGCCTCCCCTGTGTCT
CCCATCCGTTCCAGCGCCCCCACCAGCAAAGAGGAAGAAATCTCAGTCCAGGGCAACAGCCAGCTG
CTGCTATCTAG
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTAAACGGCCGGC
  
```

**Restriction Sites:** Sgfl-Mlul  
**ACCN:** NM\_001033714  
**Insert Size:** 2427 bp

<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>OTI Annotation:</b>	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"> <li>1. Centrifuge at 5,000xg for 5min.</li> <li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li> <li>3. Close the tube and incubate for 10 minutes at room temperature.</li> <li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li> <li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li> </ol>
<b>RefSeq:</b>	<u><a href="#">NM_001033714.2</a></u>
<b>RefSeq Size:</b>	2772 bp
<b>RefSeq ORF:</b>	2427 bp
<b>Locus ID:</b>	4839
<b>UniProt ID:</b>	<u><a href="#">P46087</a></u>
<b>Cytogenetics:</b>	12p13.31
<b>MW:</b>	89 kDa
<b>Gene Summary:</b>	<p>Involved in ribosomal large subunit assembly (PubMed:24120868). S-adenosyl-L-methionine-dependent methyltransferase that specifically methylates the C(5) position of cytosine 4447 in 28S rRNA (Probable). May play a role in the regulation of the cell cycle and the increased nucleolar activity that is associated with the cell proliferation (Probable).[UniProtKB/Swiss-Prot Function]</p> <p>Transcript Variant: This variant (2) differs in the 5' UTR compared to variant 1. Variants 1 and 2 encode the same isoform (1).</p>