

Product datasheet for **SC302766**

ZNF404 (NM_001033719) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ZNF404 (NM_001033719) Human Untagged Clone
Tag:	Tag Free
Symbol:	ZNF404
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC302766 representing NM_001033719.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGCATCGCC
ATGGCCCGGGTGCCATTGACATTCAGCGATGTTGCCATAGACTTCTCTCAGGAGGAGTGGGAATATTTA
AACTCGGATCAGAGGGATTTGTACAGAGATGTGATGTTGGAGAATTATACTAACTTGGTCTATTGGAC
TTTAATTTCACAACTGAAAGCAACAAGTTATCTTCAGAAAAAAGAAATTATGAAGTAAATGCGTACCAT
CAGGAGACATGGAAAAAGAAATAAAACCTTCAACCTTATGAGGTTTATTTTCAGAAGTACCCACAGTAC
ACAATTGAATTTGGGAGACAACAGAGACCTAAAGTGGGATGTTTTAGTCAAATGATATTCAAAAACAT
AAATCCCTTCTCTACATAAGAGAAATAACACAAGAGAGAAATCATATGAGTGAAGGAATAAAGAAG
GGCTTTAGAAAATATTTGCACCTTACTGAACATCTGAGAGACCATACTGGTGTGATACCCTATGAATGT
AATGAATGTGGAAAAGCATTGTAGTTTTCCAGCATTATTTAGACATCGAAAAATCCACACTGATTTG
AAACCCTATGAATGCAATGGATGTGAGAAGGCCTTTAGGTTTTATTCACAGCTTATTCAGCATCAGATA
ATTCATACTGGTATGAAACCTATGAATGTAAGCAATGCGGGAAGGCTTTTAGACGTCATTCTCACCTT
ACAGAACATCAGAAAATTCATGTTGGCTTGAAACCTTTGAATGTAAGGAATGTGGGAAACGTTTAGA
TTATATCGACATATGTGCTGCATCAGAAAATTCATCATGGTGTGAAACCTACAATGTAAGAATGT
GGAAAAGCCTTTGGTTCATCGTTCAAGTCTTTACCAACATAAGAAAATTCATTCTGGTGAGAAACCATAT
AAATGTGAACAATGTGAAAAGGCCTTTGTTGCGAGCTATCTACTTGTGAAACATCAAAGAAGTCATACT
GGTGAGAAACCTCATGAATGCATGGAATGTGAAAAGGCCTTTTGGTAAGGGCTCAAGCCTTCTAAACAT
AAGAGAATTCATAGTAGTGAGAACTCTATGATTGTAAGGATTGTGAAAAGGCCTTTTGTAGAGGCTCT
CAACTACACAGCATCAGAGAATTCATACTGGTGAAGCCACATGAATGTAAGAATGTGGGAAGACT
TTTAAGCTTCATTCATATCTTATTCAACATCAGATAATTCATACTGATTTGAAGCCATATGAATGTAAG
CAATGTGGGAAAAGCCTTCAGTCGTGTTGGAGACCTTAAGACACATCAATCAATTCATGCTGGGGAGAAA
CCCTATGAATGTAAGGAATGTGAAAACCTTTAGACTTAATTCTCAACTAATTTATCATCAGACAATT
CATACTGGTTTGAACCCATGTATGTAAGAATGTAAGAAGGCCTTTCTGTTCTATCTCAGGCTTTTCT
CAACATAAGAGAATTCATACTGGTGAAAAACCTATGAATGTAAGAATGTGATAAGGCCTTTAATCGC
AGTGATCGACTTACTCAACATGAGACAATTCATACTGGTGTGAAACCACAGAAATGCAAAGAATGTGGT
AAGGCCTTAGTCATTGCTATCAACTTAGTCAACATCAAAGATTTACCATGGTGAGAGACTCTTAATG
TAA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTAAACGGCCGCGC
  
```

Restriction Sites: Sgfl-Mlul

ACCN: NM_001033719

Insert Size: 1659 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001033719.2](#)

RefSeq Size: 1868 bp

RefSeq ORF: 1659 bp

Locus ID: 342908

UniProt ID: [Q494X3](#)

Cytogenetics: 19q13.31

MW: 65.4 kDa

Gene Summary: May be involved in transcriptional regulation.[UniProtKB/Swiss-Prot Function]