

Product datasheet for **SC302802**

ITPRIPL2 (NM_001034841) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ITPRIPL2 (NM_001034841) Human Untagged Clone
Tag:	Tag Free
Symbol:	ITPRIPL2
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)


Fully Sequenced ORF: >SC302802 representing NM_001034841.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGTCGGTGCACTACACCCTCAATCTACGCGTCTTCTGGCCCTGGTGACCGCCTGTGCACCCGCTG
GTGTGCCCTACCATGTCTGCGGGGAAGCGGGGCCCGGGCCGAGCCCGCCACGGCGTGGATGGC
GGCTTCCCGTTGCTTAAGGTGGCCGTCTGCTCCTCCTCAGCTATGTCTCCTGCGCTGTGCCACGCT
GTCCGGCAGCGCTTCTGCCCGGTCTCCCGTCTGGAGGTCACGCCGCTTCTCCTCGAGACTTC
CGAGAGCCGGCCTCAGCATCTGCTGGAGAGTTACTACGAGCATGAGGTGCGCCTGTCTCCGCACGTG
TTGGGCCACAGCAAGGCGCACGTGAGCCGATCGTGGCGAGCTGGTGGGGTGGCCGCGCCGGGGG
TCCCCCGTCTATTCTGGGGAGCGCTGGCCTTGGCCTTCCGCGGAGACTTATCCAGGTGGGCAGC
GCCTACGAGCAACATAAAATCCGCCGGCCGACAGCTTCGACGTGCTGGTCCACTGCGCCTCCCGCCG
CTTGTGGCGCTGGAGCCACGGAGCTGGCGGAGGCCAGCGCTGGCCCCGCCCTCCGCGCTGCTTC
TTGTGCGCCTCAAGGCACCACCTCACCATCGGGGGCTCGGGGGCCACTGGCTTCGGGACTGCAAA
CCTTTGCTGATGCCTTCTGCGTGGATGTGCGCGGGCGGCGTCACCTCTCTGCTACTCTGGTGTGCGC
TGGTTCAGTCGCATCTGCAGCGCTCCTTGGCCACTGTGCGTTACAGCCTGGAGGGGCGCTGCGGGTC
ACCTTGACCCAGGTGGCCTGGAACAGCCCCACCTTACACATCTTGCCTGCCGACTGACTACGGC
TGCTGCCGCTTCTATGGCTGTGCGTCTATCCCCGTGTCCATCTGGGAGATGGGGTCTTCTTGTG
GCGCCACCACCGCCACCTTGGCCAGCGCGCCCTGTTGGAGCTCCCTGAGGGCTGCGTGCGGAGGCA
CTGTGGGTGTGAACACAGCACGCCAGGAGCAGAAGCTGCTGAGTTGGCTGCAGGAACGGCAGCTCCA
GGTGCCTGCTACCTCAAGTGCTGCAGTTGCTTAAGCTCTGCGGATCTGGGGCCCGTGGCTGGAC
TCAGCGGCCGCCACCCAGTGGGACGCATCCTATCCTCATATGTGCTCAAGACAGTGTGCTGGCAGTG
CTGCTGCGCAAGGGGGCCCTGGGCAAGGCTGGGACGAGGACACCTGGGAAGGTGTTGGAGGATTG
GTGCAATTCTTAGGGACTGCCTGCTGCGACGCCATACGCTTCTCCACTGCGTCTGGGCCCTGGTGGG
GCGGCTGCCGAGGTGGTCCCCTGCCAAGGCACTGAGGGAAGCCGCCAGTTGACCTCTGGCCGCT
TTCGACGGGCACGCCGGAACTTGACGACGCGGTTGCTGTCCAGTGGCAAAGGCTGCCACGCTT
CTCCGGGCTACGGGGTCCCCTACCTTGCAGGTGCCCCACCCCGGAGTCAGCGCACCCAGGGC
TTCTTGAAGTGAACCGTAA
AGCGGACCGACGCTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGAT
ATCCTGGATTACAAGGATGACGACGATAAGGTTAA
  
```

Restriction Sites: SgfI-RsrII

Plasmid Map:



ACCN:	NM_001034841
Insert Size:	1608 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001034841.3</u>
RefSeq Size:	7722 bp
RefSeq ORF:	1608 bp
Locus ID:	162073
UniProt ID:	<u>Q3MIP1</u>
Cytogenetics:	16p12.3
Protein Families:	Transmembrane
MW:	58.4 kDa