

Product datasheet for **SC302825**

RGL3 (NM_001035223) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	RGL3 (NM_001035223) Human Untagged Clone
Tag:	Tag Free
Symbol:	RGL3
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC302825 representing NM_001035223.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGCATCGC
ATGGAGCGCACAGCAGGCAAAGAGCTGGCCCTGGCACCCTGCAGGACTGGGGTGAAGAGACCGAGGAC
GGCGCGGTGTACAGTGTCTCCCTGCGGCGGCAGCGCAGTCAGCGCAGGAGCCCGCGGAGGGCCCCGGG
GGCAGCCAGGCTCCAGCCCCATTGCCAATACCTTCTCCACTATCGAACCAGCAAGGTGAGGGTGCTG
AGGGCAGCGCGCTGGAGCGGCTGGTGGGAGAGTTGGTGTGGAGACCGTGAGCAGGACCCAGCTTC
ATGCCCGCTTCTGGCCACCTACCGGACCTTTGTACCCACTGCCTGCCTGCTGGGCTTTCTGCTGCCA
CCAATGCCACCGCCCCACCTCCCGGGGTAGAGATCAAGAAGACAGCGGTACAAGATCTGAGCTTCAAC
AAGAACCTGAGGGCTGTGGTGTCACTGCTGGGCTCCTGGCTGCAGGACACCCCTCAGGATTTCCGAGAC
CACCTGCCCATTCGGACCTGGGACGTGCCGAACCTTTCTGGGTGGGCGGCCAGGGAGTGCTGAG
GCTCAAAAGCAGAGAAGCTTCTGGAAGATTTTGGAGGAGGCTGAGCGAGAGCAGGAAGAGGAGCCG
CCTCAGGTGTGGACAGGACCTCCAGAGTTGCCAAACTTCTGACCCAGACTCTCAGAGGCTGCGCG
GAGGAAGAGGAAGGGCTCATGCCTCAAGGTCCCAGCTCCTGGACTTCAGCGTGACGAGGTGGCCGAG
CAGCTGACCTCATAGACTTGGAGCTCTCTCCAAGGTGAGGCTCTACGAGTGCTTGGGCTCCGTGTGG
TCGAGAGGGACCGCGGGGGCTGCAGGCGCCTCCCCACTGTGCGCGCCACCGTGGCCAGTTCAAC
ACCGTGACCGGCTGTGTGCTGGGTTCCGTGCTCGGAGCACCGGCTTGGCCGCCCGCAGAGGGCGCAG
CGGCTGGAGAAGTGATCCGCATCGCCAGCGCTGCCGAGAACTGCGGAACCTTCTCTCTTGCAGGCC
ATCCTGTCCGCCCTGCAATCTAACCCATCTACCGGCTCAAGCGCAGCTGGGGGGCAGTGAGCCGGGAA
CCGCTATCTACTTTCAGGAACTTTCGAGATTTTCTCCGATGAGAACAACACCTCAGCAGCAGAGAG
ATTCTTTTCCAGGAGGAGGCCACTGAGGGATCCCAAGAAGAGGACAACACCCAGGCAGCTGCCCTCA
AAACACCCCGAGGCCCTGTCCCTACCTTGGCACCTTCTTACGGACCTGGTTATGCTGGACACAGCC
CTGCCGGATATGTTGGAGGGGATCTCATTAACTTTGAGAAGAGGAGGAAGGAGTGGGAGATCCTGGCC
CGCATCCAGCAGCTGCAGAGGCGCTGTCAAGCTACACCTGAGCCCCACCCGCCCATCCTGGTGCC
CTGCATGCCCAGAACAGCTCACCGAGGAGCAGAGCTACCGGCTCTCCCGGGTATTGAGCCACAGCT
GCCTCTGCCCGAGTCCCCACGCATCCGACGGCGGATCAGCCTACCAAGCGTCTCAGTGCAAGCTT
GCCCAGAGAGAAAGCTCATCACCTAGTGGGAGTCCCGGGGACCCCTCATCCCCACCTCCAGTGTGTCC
CCAGGGTACCCCCCTCAAGTCTAGAAGCAGAGATGCTCCTGCTGGCAGTCCCCGGCCTCTCCAGGG
CCCCAGGGCCCCAGCACCAAGCTGCCCTGAGCCTGGACCTGCCAGCCCCCGGCCCTTCGCTTTGCCT
CTGGGCAGCCCTCGAATCCCCCTCCCGGCGCAGCAGAGCTCGGAGGCCGTGTATCCGCGTCAGCATC
GACAATGACCACGGGAACCTGTATCGAAGCATCTTGTGACCAAGTCAGGACAAAGCCCCAGCGTGGTC
CGGCGAGCCTTGCAAGACACAATGTGCCCCAGCCCTGGGCTGTGACTATCAGCTCTTTCAAGTCTT
CCTGGGACCGGGTGCTCCTGATTCTGACAATGCCAACGTCTTCTATGCCATGAGTCCAGTCGCCCC
AGAGACTTCATGCTGCGGCGAAAGAGGGGACCGGAACACTCTGTCTGTCTCCCAAGCTGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: SgfI-MluI

ACCN: NM_001035223

Insert Size: 2133 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001035223.3</u>
RefSeq Size:	2553 bp
RefSeq ORF:	2133 bp
Locus ID:	57139
UniProt ID:	<u>Q3MIN7</u>
Cytogenetics:	19p13.2
MW:	78.1 kDa
Gene Summary:	Guanine nucleotide exchange factor (GEF) for Ral-A. Potential effector of GTPase HRas and Ras-related protein M-Ras. Negatively regulates Elk-1-dependent gene induction downstream of HRas and MEKK1 (By similarity).[UniProtKB/Swiss-Prot Function] Transcript Variant: This variant (2) uses an alternate in-frame splice site in the 3' coding region compared to variant 1. This results in a shorter protein (isoform 2) compared to isoform 1.