

## Product datasheet for **SC302871**

### GGPS1 (NM\_001037277) Human Untagged Clone

#### Product data:


Product Type:	Expression Plasmids
Product Name:	GGPS1 (NM_001037277) Human Untagged Clone
Tag:	Tag Free
Symbol:	GGPS1
Synonyms:	GGPPS; GGPPS1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC302871 representing NM_001037277. Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTT TAGTGAACCGTCAGAATTTTGT AATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCC GCGATCGCC
ATGGAGAAGACTCAAGAAACAGTCAAAGAATTCTTCTAGAACCCATAAATACTTACTTCAGTTACCA
GGTAAACAAGTGAGAACCAAACCTTTCACAGGCATTTAATCATTGGCTGAAAGTTCAGAGGACAAGCTA
CAGATTATTATTGAAGTGACAGAAATGTTGCATAATGCCAGTTTACTCATCGATGATTTGAAGACAAC
TCAAACCTCCGACGTGGCTTCCAGTGGCCACAGCATCTATGGAATCCCATCTGTCATCAATTCTGCC
AATTACGTGATTTCTTGGCTTGGAGAAAGTCTTAACCTTGATCACCCAGATGCAGTGAAGCTTTTT
ACCCGCCAGCTTTTGGAACTCCATCAGGGACAAGGCCTAGATATTTACTGGAGGGATAATTACACTTGT
CCCACTGAAGAAGATATAAAGCTATGGTCTGCAGAAAACAGGTGGACTGTTTGGATTAGCAGTAGGT
CTCATGCAGTTGTTCTCTGATTACAAAGAAGATTTAAAACCGCTACTTAATACACTTGGGCTCTTTTTC
CAAATTAGGGATGATTATGCTAATCTACACTCAAAGAATATAGTGAAAACAAAAGTTTTTGTGAAGAT
CTGACAGAGGGAAAGTTCATTTCCTACTATTCATGCTATTTGGTCAAGGCCTGAAAGCACCCAGGTG
CAGAATATCTTGCGCCAGAGAACAGAAAACATAGATATAAAAAATACTGTGTACATTATCTTGAGGAT
GTAGGTTCTTTGAATACACTCGTAATACCTTAAAGAGCTTGAAGCTAAAGCCTATAAACAGATTGAT
GCACGTGGTGGGAACCTGAGCTAGTAGCCTTAGTAAAACACTTAAGTAAGATGTTCAAAGAAGAAAAA
GAATAA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites: SgfI-MluI



[View online »](#)

**Plasmid Map:**


**ACCN:** NM\_001037277

**Insert Size:** 903 bp

**OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

**OTI Annotation:** This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_001037277.1](#)

**RefSeq Size:** 2757 bp

**RefSeq ORF:** 903 bp

<b>Locus ID:</b>	9453
<b>UniProt ID:</b>	<a href="#">O95749</a>
<b>Cytogenetics:</b>	1q42.3
<b>Protein Pathways:</b>	Metabolic pathways, Terpenoid backbone biosynthesis
<b>MW:</b>	34.9 kDa
<b>Gene Summary:</b>	<p>This gene is a member of the prenyltransferase family and encodes a protein with geranylgeranyl diphosphate (GGPP) synthase activity. The enzyme catalyzes the synthesis of GGPP from farnesyl diphosphate and isopentenyl diphosphate. GGPP is an important molecule responsible for the C20-prenylation of proteins and for the regulation of a nuclear hormone receptor. Alternate transcriptional splice variants, both protein-coding and non-protein-coding, have been found for this gene. [provided by RefSeq, Sep 2010]</p> <p>Transcript Variant: This variant (2) encodes the functional protein.</p>