

## Product datasheet for **SC302927**

### DEFB110 (NM\_001037728) Human Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	DEFB110 (NM_001037728) Human Untagged Clone
Tag:	Tag Free
Symbol:	DEFB110
Synonyms:	DEFB-10; DEFB-11; DEFB111
Vector:	<u>pCMV6 series</u>
Fully Sequenced ORF:	>NCBI ORF sequence for NM_001037728, the custom clone sequence may differ by one or more nucleotides ATGAAGATTCAACTTTTTTTCTTTATTCTGCACTTTTGGGTCACAATTTTACCAGCCAGA AGCAATTTTGAACCAAAGTATAGATTTGAAAGATGCGAAAAAGTGAGAGGAATATGTAAA ACGTTTTGTGATGATGTTGAGTATGATTATGGATACTGCATTAATGGAGAAGTCAGTGC TGCGTATAA
Restriction Sites:	Please inquire
ACCN:	NM_001037728
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>



[View online »](#)

RefSeq:	<u>NM_001037728.1, NP_001032817.1</u>
RefSeq Size:	189 bp
RefSeq ORF:	189 bp
Locus ID:	245913
UniProt ID:	<u>Q30KQ9</u>
Cytogenetics:	6p12.3
Protein Families:	Secreted Protein
Gene Summary:	<p>Defensins form a family of antimicrobial and cytotoxic peptides made by neutrophils. Defensins are short, processed peptide molecules that are classified by structure into three groups: alpha-defensins, beta-defensins and theta-defensins. All beta-defensin genes are densely clustered in four to five syntenic chromosomal regions. Two transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Oct 2014]</p> <p>Transcript Variant: This variant (2) lacks a 3' UTR and differs in the 3' coding region, compared to variant 1, resulting in an isoform (b) that has a distinct C-terminus and is shorter than isoform a.</p>