

Product datasheet for **SC302939**

eIF3s8 (EIF3C) (NM_001037808) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	eIF3s8 (EIF3C) (NM_001037808) Human Untagged Clone
Tag:	Tag Free
Symbol:	EIF3C
Synonyms:	eIF3-p110; EIF3CL; EIF3S8
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)




[View online »](#)

Fully Sequenced ORF: >SC302939 representing NM_001037808.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGTCGCGGTTTTCCACCACCGTTCCGACAGCGAGTCCGAGTCGTCCTTGTCCGGGGAGGAGCTCGTC
ACCAAACCTGTCGGAGGCAACTATGGCAAACAGCCATTGTTGCTGAGCGAGGATGAAGAAGATACCAAG
AGAGTTGTCGCGAGTGCCAAGGACAAGAGTTTGAGGAGCTGACCAACCTTATCCGGACCATCCGTAAT
GCCATGAAGATTCGTGATGTACCAAGTGCCTGGAAGAGTTTGAGCTCCTGGGAAAAGCATATGGGAAG
GCCAAAAGCATTGTGGACAAAGAAGGTGTCCCCGGTTCTATATCCGCATCCTGGCTGACCTAGAGGAC
TATCTTAATGAGCTTTGGGAAGATAAGGAAGGGAAGAAGAAGTGAACAAGAACAATGCCAAGGCTCTG
AGCACCTTGCCTCAGAAGATCCGAAAATACAACCGTGATTTGAGTCCCATATCACAAGCTACAAGCAG
AACCCCGAGCAGTCTGCGGATGAAGATGCTGAGAAAAATGAGGAGGATTCAGAAGCTCTTCAGATGAG
GATGAGGATGAGGACGGAGTCAGTGTGCAACTTTCTTGAAGAAGAAATCAGAAGCTCCTCTGGGGAG
AGTCGCAAGTTCCCTAAAAAGATGGATGATGAAGATGAGGACTCAGAAGATTCGGAAGATGATGAAGAC
TGGGACACAGGTTCCACATCTTCCGATTCGACTCAGAGGAGGAAGAAGGAAAACAAACCCGCTGGCC
TCAAGATTTCTTAAAAAGCACCCACCACAGATGAGGACAAGAAGGCAGCCGAGAAGAAACGGGAGGAC
AAAGCTAAGAAGAAGCACGACAGGAAATCCAAGCGCTGGATGAGGAGGAGGAGGACAATGAAGGCGGG
GAGTGGGAAAGGGTCCGGGGCGGAGTCCGTTGGTTAAGGAGAAGCCAAAAATGTTTGCCAAGGGAAT
GAGATCACCCATGCTGTTGTTATCAAGAAACTGAATGAGATCCTACAGGCACGAGGCAAGAAGGGAAT
GATCGTCTGCCAGATTGAGCTGCTGCAACTGCTGGTTCAGATTGCAGCGGAAAACAACCTGGGAGAG
GGCGTCATTGTCAAGATCAAGTTCAATATCATCGCCTCTCTATGACTACAACCCCAACCTGGCAACC
TACATGAAGCCAGAGATGTGGGGAAGTGCCTGGACTGCATCAATGAGCTGATGGATATCCTGTTGCA
AATCCCAACATTTTTGTTGGAGAGAATATTCTGGAAGAGAGTGAGAACCTGCACAACGCTACCAGCCA
CTGCGTGTCCGTGGCTGCATCCTAACTCTGGTGAACGAATGGATGAAGAATTTACCAAAAATAATGCAA
AATACTGACCCTCACTCCCAAGAGTACGTGGAGCACTTGAAGGATGAGGCCAGGTGTGTGCCATCATC
GAGCGTGTGCAGCGCTACCTGGAGGAGAAGGGCACTACCGAGGAGGTCTGCCGCATCTACCTGCTGCGC
ATCCTGCACACCTACTACAAGTTTGATTACAAGGCCCATCAGCGACAGCTGACCCCGCTGAGGGCTCC
TCAAAGTCTGAGCAAGACCAGGCAGAAAATGAGGGCGAGGACTCGGCTGTGTTGATGGAGAGACTGTGC
AAGTACATCTACCCAAGGACCCACAGACCGGATCCGCACATGTGCCATCCTCTGCCACATCTACCAC
CATGCTCTGCACTCGCGTGGTACCAGGCCCGCAGCTCATGCTCATGAGCCACTTGAGGACAACATT
CAGCATGCAGACCCGCCAGTGCAGATCCTTTACAACCGCACCATGGTGCAGCTGGGCATCTGTGCCCTC
CGCCAAGGCTGACCAAGGACGCACACAACGCCCTGCTGGACATCCAGTCGAGTGGCCGAGCCAAGGAG
CTTCTGGGCCAGGGCTGCTGCTGCGCAGCCTGCAGGAGCGCAACCAGGAGCAGGAGAAGGTGGAGCGG
CGCCGTGAGTCCCTTCCACCTGCACATCAACCTGGAGCTGCTGGAGTGTGTCTACCTGGTGTCTGCC
ATGCTCCTGGAGATCCCTACATGGCCGCCATGAGAGCGATGCCCGCCGACGCATGATCAGCAAGCAG
TTCCACCACCAGCTGCGCGTGGCGGAGCAGACAGCCCTGCTGGTCCCTGAGTCCATGCGGGAACAT
GTGGTCTGCTGCCCAAGGCCATGAAGATGGGTGACTGGAAGACCTGTACAGTTTTATCATCAATGAG
AAGATGAATGGAAAGTGTGGACCTTTTCCCGAGGCTGACAAAGTCCGACCATGCTGGTTAGGAAG
ATCCAGGAAGAGTCACTGAGGACCTACCTCTTACCTACAGCAGTGTCTATGACTCCATCAGCATGGAG
ACGCTGTGAGACATGTTTGGAGCTGGATCTGCCACTGTGCACTCCATCATCAGCAAAAATGATCATTAA
TGGAGCTGATGGCTCCCTGGACCAGCCAACACAGACAGTGGTGTGATGCACCGCACTGAGCCCACTGCC
CAGCAGAACCTGGCTCTGCACTGGCCGAGAAGCTGGGCAGCCTGGTGGAGAACAACGAACGGGTGTTT
GACCACAAGCAGGGCACCTACGGGGGCTACTTCCGAGACCAGAAGGACGGCTACCGCAAAAACGAGGGC
TACATGCGCCGCGGTGGTACCGCCAGCAGCAGTCTCAGACGGCTACTGA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTAAACGGCCGGC
  
```

Restriction Sites: SgfI-MluI

Plasmid Map:


ACCN: NM_001037808

Insert Size: 2742 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001037808.2](#)

RefSeq Size: 3145 bp

RefSeq ORF: 2742 bp

Locus ID: 8663

UniProt ID: [Q99613](#)

Cytogenetics: 16p11.2

Protein Families: Druggable Genome

MW: 105.3 kDa

Gene Summary: Component of the eukaryotic translation initiation factor 3 (eIF-3) complex, which is required for several steps in the initiation of protein synthesis (PubMed:17581632, PubMed:25849773, PubMed:27462815). The eIF-3 complex associates with the 40S ribosome and facilitates the recruitment of eIF-1, eIF-1A, eIF-2:GTP:methionyl-tRNAⁱ and eIF-5 to form the 43S pre-initiation complex (43S PIC). The eIF-3 complex stimulates mRNA recruitment to the 43S PIC and scanning of the mRNA for AUG recognition. The eIF-3 complex is also required for disassembly and recycling of post-termination ribosomal complexes and subsequently prevents premature joining of the 40S and 60S ribosomal subunits prior to initiation (PubMed:17581632). The eIF-3 complex specifically targets and initiates translation of a subset of mRNAs involved in cell proliferation, including cell cycling, differentiation and apoptosis, and uses different modes of RNA stem-loop binding to exert either translational activation or repression (PubMed:25849773).[UniProtKB/Swiss-Prot Function]
Transcript Variant: This variant (2) differs in the 5' UTR compared to variant 1. Variants 1-4 encode the same isoform (a).