

Product datasheet for **SC302941**

ZNF284 (NM_001037813) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ZNF284 (NM_001037813) Human Untagged Clone
Tag:	Tag Free
Symbol:	ZNF284
Synonyms:	ZNF284L
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC302941 representing NM_001037813.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGCATCGCC
ATGACCATGTTCAAGGAGGCAGTGACCTTCAAGGATGTGGCTGTGGTCTTCACCGAGGAGGAGCTGGGG
CTGCTGGACGTTTCCAGAGGAAGCTGTATCGAGATGTCATGCTGGAGAACTTCAGAAACCTGCTATCA
GTGGGGCATCAACTTTCCACCGAGATACTTTTCACTTCCAAAGAGAAGAAAAGTTTGGATCATGGAG
ACAGCAACCCAAAGAGAAGGAATTGAGGAGGCAAGATCCAAACTGAGTTGGAGTCTGTTCCAGAAACA
GGACCACATGAAGAGTGGTCTTGCCAGCAAACTGGAACAACTGCAAGTGAGTTAACTAGACCTCAA
GACTCCATAAGTAGCTCTCAGTTCTCCACACAAGGTGATGTCCCTCCAGGTTGACGCAGGACTATCT
ATAATTCACATAGGAGAGACACCTTCTGAGCATGGGAAGTGTAAAAATCTTTCAGTGATGTCTCCATC
CTTGATCTTCATCAACAATTACACTCAGGAAAGATATCCCATACATGTAATGAGTACAGGAAGAGATTTC
TGTTATAGCTCAGCTCTTTGTCTTCATCAGAAAGTTCACATGGGAGAGAAACGCTATAAGTGTGATGTG
TGTAGTAAGGCATTTAGTCAGAACTCACAACCTGCAAACTCATCAGAGAATCCACACTGGAGAGAAACCA
TTCAAATGTGAGCAGTGTGGGAAAAGTTTCAGCCGTAGATCAGGAATGTATGTTCAATTGCAAATTACAC
ACAGGAGAAAAACCTCATATTTGTGAGGAATGTGGGAAGGCCTTCATTCACAATCCAGCTTCGGGAA
CATCAAAGAATCCATACTGGGGAAGCCATTCAAATGTTATATATGTGGTAAGAGCTTCCATAGTAGA
TCAAATCTTAATAGGCATTCCATGGTCCACATGCAAGAGAAAATCATTTAGATGTGATACCTGTAGTAAT
AGCTTTGGTCAGAGATCAGCACTTAATAGTCATTGCATGGACCACACAAAAGAGAACTATACAAATGT
GAAGAATGTGGAAGGAGCTTCACTTGTAGGCAAGATCTTTGTAAGCATCAGATGGACCATACAGGAGAC
AAACCATATAATTGTAATGTATGTGGGAAGGGCTTCAGGTGGTCCTCATGTCTTTCAAGACATCAGCGG
GTCCACAATGGAGAAACAACATTCAAGTGCAGCGATGTGGGAAGAGATTTATATGAATTCACAGGGC
CATTACATCAGAGAGCCTATAGAGAAGAAGAACTGTATAAATGTCAGAAAGTGTGGGAAGGGCTACATT
AGTAAGTTTAATCTTGACTTGCACAGAGGGTCCACACGGGAGAGAGACCTTATAATTGTAAGGAATGT
GGAAAGAGCTTCAGGTGGGCCTCAGGTATTTTGAGACATAAGAGACTCCATACTGGAGAAAAACCATTC
AAATGTGAAGAGTGTGGGAAGAGGTTTACTGAGAATTCAAAACCTCGTTTCCATCAAAGAATTCACACT
GGAGAAAAAGCCTTACAAATGTGAGGAGTGTGGAAGGGATTGAGATGGGCCTCAACTCATCTAACCCAT
CAAAGACTCCACAGCAGAGAAAACTATTCCAATGTGAGGATTGTGGGAAGAGCAGTGAGCACAGTTCA
TGCCTTCAAGACCAACAAAGCGACCACAGTGGAGAAAAACATCCAATGTGAGGACTGTGGGAAGCGC
TACGAGAGGCGCTGAATCTAGATATGATTTTATCATTATTTTAAATGATATATAA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: Sgfl-MluI

ACCN: NM_001037813

Insert Size: 1782 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001037813.3</u>
RefSeq Size:	4368 bp
RefSeq ORF:	1782 bp
Locus ID:	342909
UniProt ID:	<u>Q2VY69</u>
Cytogenetics:	19q13.31
MW:	69 kDa
Gene Summary:	May be involved in transcriptional regulation.[UniProtKB/Swiss-Prot Function]