

Product datasheet for **SC302944**

DIXDC1 (NM_001037954) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	DIXDC1 (NM_001037954) Human Untagged Clone
Tag:	Tag Free
Symbol:	DIXDC1
Synonyms:	CCD1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)




[View online »](#)

Fully Sequenced ORF: >SC302944 representing NM_001037954.
Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGCATCGCC
ATGCTAGCCTGCCTGACCCGAGGGAACCTACTGGACGTCCTCGAGGAGGGCTTCAATGAGCAACAGCTG
CAGGCCTATGTGGCCTGGGTGAATGCACAGCTGAAGAAGAGGCCAGCAGTGAAGCCTGTGCAGGACCTG
CGACAAGATCTCCGGGATGGGGTATCCTGGCATACTCATCGAGATTGTTGCAGGAGAAAAGCTGAGT
GGGTACAGCTGAGTCCCGTAACCAACAGGAGATGAAGAATAATGTGGAGAAAAGTCTACAGTTTGTG
GCCTCTAAAAAGATTTCGTATGCACCAGACTTCGGCTAAAGATATTGTGGATGAAAACCTGAAGTCTATC
ATGAGGCTGGTCTTGCCTTGGCAGCTCATTTCAAACCTGGCTCCAGCAGGACGGTGAACCAAGGACGG
GACTCCAGAGCCCTCTGCAAAGTACCAGCACACTGTGCCACTGCTGTTGCCAGGGAGCAGCTGCT
GCTCTGGCCGATGTGTGCATGACATGTCCGATCAGGACGGGATGCTTTTCGATATAGACAGAGGAAC
AGCAGCATGGATGAGGAAATTGAGAATCCATACTGGAGTGTGCGGGCCCTAGTGCAGCAGTACGAAGGG
CAACAAAGTCCCGTCTGAATCCAGCTGCTCCAGCCTGACTTCAACCAGTCCAATCCACAGTGAAGG
AGCGAGTCCATTATAACCGAGTCAAGAGAAGGCAGATTTTGTGATTATCCCGCTGAAGGAATAGAG
AACAGAACAGAAGGGACAGATTCTCCATTATCTCGAGACTGGCGGCCAGGGAGCCCTGGAACCTATCTG
GAGACCTCATGGGAAGAACAGCTGTTGGAACAACAAGATATTTAGAAAAAGAAATGGAGGAAGCAAAG
AAAATGATATCAGGACTACAGGCCTTACTGCTCAATGGATCCTTACCTGAAGATGAACAGGAGAGGCC
TTGGCCCTCTGTGAACAGGTGTCAATCCCGAGGAACAACCTGATTATAATCCAAAGTCGTCTGGATCAG
AGTATGGAGGAGAATCAGGACTTAAAGAAGGAAGTGTGAAATGTAAACAAGAAGCCAGAAAACCTACAG
GGGATAAAGGATGCCTTGCAGCAGAGATTGACTCAGCAGGACACATCTGTTCTTCAGCTCAAACAAGAG
CTACTGAGGGCAAATATGGACAAAGATGAGCTGCACAACCAGAATGTGGATCTGCAGAGGAAAGTAGAT
GAGAGGAACCGGCTCTTGGGAGAATATAAAAAAGAGCTGGGGCAGAAGGATCGCCTTCTTCAGCAGCAC
CAGGCCAAGTTAGAAGAAGCACTCCGAAAACCTCTGATGTGAGTTACCACCAGGTGGATCTAGAGCGA
GAGCTAGAACACAAGATGTCTCTTGGCTCACTGTATGAAAAGAGAGGCAGATGAGGCGACCAACTAC
AACAGTCACAACTCTCAAAGCAATGGTTTTCTCTTCCAACGGCAGGAAAAGGAGCTACTTCAGTCAGC
AACAGAGGGACCAGCAGCTGCAGCTTGTTCGAGATGCTCTCCGAGCCTGCGCAACAGCTTCAGTGGC
CACGATCCTCAGCACCACTATTGACAGCTTGGAGCAGGGCATTCTAGCCTCATGGAGCGCCTGCAT
GTTATGGAGAGCAGAGAAGAAACAAGAAAGAAAGGTTCCGGTCAAGTCAACCAGAAGTCAAGTAGGTAGT
GAATACCGGGAGTCTGGCCCCCTAACTCAAAGTTGCCTCACTCACAGAGCTCTCAACTGTCAGCAGC
ACCTGTACTAAAGTGTCTATTTCACTGACCGGCTCACTTACGCCCTTCATGGTCAATATACCAAAGAGG
TTGGAGGAGGTGACGTTAAAGGATTTTAAAGCAGCTATTGATCGGGAAGGAAATCACCGGTATCACTTC
AAAGCACTGGATCCTGAGTTTGGCACTGTCAAAGAGGAGATTTTCCATGATGATGATGCCATCCCTGGA
TGGGAAGGAAAAATTGTAGCTTGGGTGGAAGAAGACCATGGAGAGAATTAA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites: SgfI-MluI

Plasmid Map:


ACCN: NM_001037954

Insert Size: 2052 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001037954.3](#)

RefSeq Size: 5963 bp

RefSeq ORF: 2052 bp

Locus ID: 85458

UniProt ID: [Q155Q3](#)

Cytogenetics: 11q23.1

MW: 77.5 kDa

Gene Summary: The protein encoded by this gene is a positive regulator of the Wnt signaling pathway. The encoded protein is found associated with gamma tubulin at the centrosome. Alternative splicing results in multiple transcript variants encoding different isoforms. [provided by RefSeq, Jun 2013]
Transcript Variant: This variant (1) represents the longest transcript and encodes the longest isoform (a, also known as I-DIXDC1).